

# SilverBoom 300



Shown: Interchangeable appearance panel;  
comes with 3 interchangeable appearance panels:



The SilverBoom 300 easily fits into most vehicles and is a perfect complement to a Jazzy® Power Chair or Pride® Scooter. Engineered to provide smooth, efficient operation, maximum durability, and sporty style, the SilverBoom 300 makes it easy and affordable to take your mobility product with you wherever you go.

- Lifting capacity 300 lbs.
- Low headroom required - 30"
- Telescoping base legs and extending boom arm length for maximum installation adjustability
- Boom arm end cap to protect vehicle headliners
- End swivel hook to prevent strap twisting
- Stop warning on strap to identify lifting limit
- Lift is easily disassembled for stowing
- Integrated boom arm securement when not stowed
- Remote hand control operation
- Integrated dual protective devices to prevent lift or vehicle electrical system damage
- Interchangeable appearance panel to complement your vehicle's interior
- **COMING SOON!** Optional power rotating base
- **COMING SOON!** Wireless remote operation

## SilverStar™ Lifts

A Division of **Pride** Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • (800) 800-8586 (U.S.) • (888) 570-1113 (Canada)  
[www.silverstarmobility.com](http://www.silverstarmobility.com) • [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

# SilverBoom 300



Lift ready to load the product.



Loading the product.



Product loaded into vehicle.

***Can easily lift the following Jazzy® Power Chairs<sup>1</sup>***

- Jazzy® 600
- Jazzy® 610
- Jazzy® 1103 Ultra
- Jazzy® 1107
- Jazzy® 1113 ATS
- Jazzy® 1121
- Jazzy® 1122
- Jazzy® Select
- Jazzy® Select 14
- Jazzy® Select 14 XL

***Can easily lift the following Quantum® Power Chairs<sup>1</sup>***

- Dynamo® 1107
- Quantum® 600
- Quantum® 610
- Quantum® 1107
- Quantum® 1121
- Quantum® 1402
- Quantum® 1420
- Vibe®

***Can easily lift the following Pride® Scooters***

- Celebrity® X
- Go-Go Elite Traveller™
- Go-Go® Ultra X
- Hurricane®
- Legend®
- Legend® XL
- Maxima
- Rally®
- Revo®
- Sonic®
- Sundancer®
- Victory®

***Can easily lift the following Jet™ Power Chairs<sup>1</sup>***

- Jet™ 3 Ultra
- Jet™ 7
- Jet™ 10 Ultra

<sup>1</sup> Seating system dependent

## SilverStar™ Lifts

A Division of **Pride®** Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • (800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

[www.silverstarmobility.com](http://www.silverstarmobility.com) • [www.pridemobility.com](http://www.pridemobility.com)