

NEW!

The Super Low "FLOOR HUGGER"

*Lower than a
standard wheelchair
for easy transfers*



Standard 2 or 3 Function Package

- ❖ Head Adjustable
- ❖ Foot Adjustable
- ❖ Hi-Low Adjustable (2 Function model will not have this operation)
- ❖ 6 Button Hand Control



Standard 4 Function Package

- ❖ Head Adjustable
- ❖ Foot Adjustable
- ❖ Hi-Low Adjustable
- ❖ Reverse Trendelenburg
- ❖ 8 Button Sleek Hand Control



Deluxe 3 Function Package

- ❖ Head Adjustable
- ❖ Foot Adjustable
- ❖ Hi-Low Adjustable
- ❖ 6 Button Sleek Hand Control
- ❖ Deluxe Illuminated Hand Control
- ❖ Single Variable Massage (operated only by Deluxe Control)
- ❖ Soft Illuminated Floor Light Feature



Quick Facts

- ❖ Super Low 9.5 inches from deck to the floor
- ❖ Height range - 12 inches
- ❖ 300 pound weight capacity
- ❖ Robust, modular control box design
- ❖ Tool-free, easy setup split bed
- ❖ Single or dual hand controls
- ❖ ECU compatible**
- ❖ Compatible with most Patient Lift systems**

Available Sizes**

Twin	38"x80"
Twin SL	38"x84"
Full	53"x80"
Full SL	53"x84"
Queen	60"x80"
Also available in Dual King	