

SHOPRIDER



K0823

STREAMER SPORT [888WA]

Built only with the best technology, the Streamer Sport includes excellent features such as puncture proof 10" tires, off board charger and a dramatically sloped chassis cover. The Streamer Sport is also equipped with self diagnosing Dynamic controllers. The Streamer Sport combines power, comfort and reliability.

Toll-Free: 800.743.0772
www.shoprider.com




Pan Seat Option

REAR WHEEL DRIVE

- Adjustable angle foot plate
- Easy maneuverability
- 10" wheel
- 300 lb. capacity
- Height and angle adjustment arm rests
- Optional Swing away controller arm bracket

STREAMER SPORT



Red Blue

| | | |
|--|--------------|----------------|
| Overall Dimensions (L x W x H) | in | 36 x 24 x 46 |
| Number, Size of Tires (front; rear) | Front (rear) | 2, 6"; 2, 10" |
| Maximum User Weight (on a level road) | lb | 300 |
| Total Weight (w/26 Ah batteries, seat) | lb | 154 |
| Each Battery Module Weight | lb | 21 |
| Battery Module Capacity | _V_Ah x Pcs | 12V26Ah x 2 |
| Charger | | Off Board (5A) |
| Output Power of Motor | hp | 0.60 x2 |
| Maximum Speed | mph | 5 |
| Max. Climbing Angle/Safe Climbing Angle(based on user weight of 90kg/200lb.) | % (°) | 10 (6) |
| Range (per charge) (After the battery and mechanical moving parts fully break in) | mile | 10 |
| Turning Radius | in | 17.7-25.5 |

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical moving parts fully break in

h) Etc.