

We've Got Your Back



Seating and Positioning **JAY® J3™**

JAY®
SEATING



When buying a back, most people expect it to work without complication: the shell fits just right, the hardware mounts easily, contouring matches exactly, and moving the back on and off the chair is effortless. For most products, this is the ideal not the reality....until the J3™ Back. Never before have the creative minds of JAY® spent more time studying the subtleties of fitting the back to the chair and the user to the back. Their result: The J3 Back. The most comprehensive combination of adjustable hardware, individual shell sizes and positioning capabilities ever assembled.

Now more than ever, we've got your back.

Seating and Positioning JAY® J3™

Multiple mounting locations without compromising back position



- **JAYMount™ hardware designed to fit:** Engineered to mount on crowded back posts ranging from 1 1/8" to 3/4", JAYMount hardware virtually eliminates incompatibility. All three dimensions of height, angle and depth can be adjusted while the user remains in the chair.



- **A multitude of shapes and sizes:** "We've got your back". Providing the correct size back for each individual is critical to preventing problems down the road. For this reason, the J3 back offers virtually every combination of depth, height and width that an adult user may require.

- **Spine-Align Posture Conforming System:** The Spine Align System is simple to use and creates the individual shape necessary to best support the user. Result: Optimal postural support, pressure distribution and comfort.

For more information, please visit J3.SunriseMedical.com

© 2007 Sunrise Medical®, JAY® and J3™ are registered trademarks of Sunrise Medical.