

R82[®] 



swanTM



Bath and toilet chair
cleanness




The Swan provides helpers with the best working conditions and allows users to preserve their dignity to the maximum. The seating comfort is outstanding and the highly adjustable frame adds flexibility.

A wide range of accessories guarantee the correct support for every single child.

Upholstery
Tray
Push handles
Pan
Foot support
*Vest and belts
*Head support



The Swan comes in seven different sizes - the sizes 4, 5, 6 and 7 with a perfect anatomical shape. The Swan fits over most toilets and can also be equipped with a pan.


swan™ by **R82**®

Measurements

	Size 1 cm (inch)	Size 2 cm (inch)	Size 3 cm (inch)	Size 4 cm (inch)	Size 5 cm (inch)	Size 6 cm (inch)	Size 7 cm (inch)
Height, with wheels (H)	68-84 (26½-32¼")	68-84 (26½-32¼")	79-95 (30¾-37")	71-87 (27¾-34")	78-94 (30½-36¾")	86-102 (33½-39¾")	86-102 (33½-39¾")
Width, with wheels (W)	49½ (19¾")	49½ (19¾")	53½ (20¾")	53½ (20¾")	53½ (20¾")	57½ (22½")	65½ (25½")
Length, with wheels (L)	46 (18")	46 (18")	57 (22")	57 (22")	57 (22")	62 (24¾")	62 (24¾")
Seat height, with wheels (SH)	37-53 (14½-20¾")	37-53 (14½-20¾")	48-64 (18½-25")	48-64 (18½-25")	48-64 (18½-25")	48-64 (18½-25")	48-64 (18½-25")
Back width, shoulders (BWS)	22½ (8¾")	22½ (8¾")	25 (9¾")	31 (12")	39 (15¼")	46 (18")	49 (19")
Back width, hips (BWH)	24 (9¾")	24 (9¾")	25½ (10")	28 (11")	31 (12")	35½ (13¾")	43 (16¾")
Back height (BH)	38 (14¾")	38 (14¾")	42 (16½")	40 (15½")	48 (18½")	57 (22")	57 (22")
Seat width, back (SWB)	20 (7¾")	20 (7¾")	23 (9")	25 (9¾")	30 (11¾")	34 (13¼")	41 (16")
Seat width, front (SWF)	30 (11¾")	30 (11¾")	34 (13¼")	36 (14")	38 (14¾")	43 (16¾")	48 (18½")
Seat depth (SD)	20 (7¾")	27 (10½")	32½ (12¾")	32½ (12¾")	38 (14¾")	46 (18")	46 (18")
Angle (A)	18°	18°	18°	27°	22°	20°	20°
Height, without wheels	59-75 (23-29¾")	59-75 (23-29¾")	70-86 (27¼-33½")	62-78 (24¼-30¾")	69-85 (27-33¾")	77-93 (30-36¾")	77-93 (30-36¾")
Width, without wheels	46 (18")	46 (18")	50 (19½")	50 (19½")	50 (19½")	54 (21")	62 (24¾")
Length, without wheels	39 (15")	39 (15")	50 (19½")	50 (19½")	50 (19½")	55 (21½")	55 (21½")
Seat height, without wheels	28-44 (11-17")	28-44 (11-17")	39-55 (15¼-21½")	39-55 (15¼-21½")	39-55 (15¼-21½")	39-55 (15¼-21½")	39-55 (15¼-21½")
Weight	7,5 kg (16,5 lb)	8 kg (17,5 lb)	9,5 kg (21 lb)	8,5 kg (18,5 lb)	9,5 kg (21 lb)	10,5 kg (23 lb)	11,5 kg (25,5 lb)
Max load	30 kg (66 lb)	30 kg (66 lb)	60 kg (132 lb)	70 kg (154 lb)	80 kg (176 lb)	90 kg (198 lb)	110 kg (242 lb)

R82 A/S

Parallevej 3
8751 Gedved
Denmark
Tel +45 79 68 58 88
Fax +45 75 66 51 92
E-mail: r82@r82.com
Web: www.r82.com

R82 UK Ltd.

Unit D4A
Coombswood Business Park East
Coombswood Way, Halesowen
West Midlands B62 8BH
Great Britain
Tel 0121 561 2222
Fax 0121 559 5437
E-mail: info@r82.com
Web: www.r82.com

Distributor: