

Invacare® Top End® Excelerator™ Stock Handcycle



PROVIDER / CUSTOMER INFORMATION:			
Provider Information:		End User Information:	
Date:	Name:	Height:	
IVC Acct #:	Ship to Address:	Weight (lbs.):	
Contact:	Address Line 2:	Leg Length (in.):	
Phone #:	City:	Hip Width (in.):	
Fax/Email:	State:	Zip:	Phone #:
<input type="radio"/> Ship to Provider	<input type="radio"/> Ship to End User		Email:

ORDER REQUEST TYPE:		
<input type="radio"/> Quote	<input type="radio"/> Order	PO # <input type="text"/>

BASE MODEL:		
<input checked="" type="radio"/> 10006	Excelerator Stock Handcycle	\$3,961

Assembly: <i>(see definitions)</i>		
<input checked="" type="radio"/> 10048	Assembled	\$0

SEAT WIDTH:		
<input type="radio"/> 15017	16" Seat Width - Simple Red Frame	\$0
<input type="radio"/> 15021	18" Seat Width - Royal Blue Frame	\$0
<input type="radio"/> 15025	20" Seat Width - Simple Red Frame	\$0

CAMBER ANGLE:		
<input checked="" type="radio"/> 45002	3 Degree Camber	\$0

WHEELBASE: <i>(Front to Rear Axle)</i>		
<input checked="" type="radio"/> 75085	40" Wheelbase	\$0

WHEEL SETS:		
<input checked="" type="radio"/> 50121	24" Front Wheels & 24" Rear Wheels with Quick Release Axles	\$0

TIRES:		
<input checked="" type="radio"/> 50213	24" X 1.75" Treaded Tires	\$0

FOOTREST:		
<input checked="" type="radio"/> 60002	Dual, Height/Angle Adjustable	\$0

DRIVE CONTROL:		
CRANKS:		
<input checked="" type="radio"/> 65054	Crank Arms: 165mm Long by 12" Wide	\$0
HAND PEDALS:		
<input checked="" type="radio"/> 65127	Medium 1-3/4" Ovalized	\$0
SHIFTING:		
<input checked="" type="radio"/> 65060	Twist Activated Shifter with Internal Gearing	\$0
BRAKING:		
<input checked="" type="radio"/> 65080	Reversing Activated Drum Brake	\$0
<input checked="" type="radio"/> 65090	Locking Parking Brake	\$0
CHAINRING:		
<input checked="" type="radio"/> 65100	Single 38 Tooth Chainring	\$0

Invacare® Top End® Excelerator™ Stock Handcycle



UPHOLSTERY:		
MATERIAL & COLOR:		
<input checked="" type="radio"/> 30001	Nylon, Black Upholstery	\$0
SEAT UPHOLSTERY:		
<input checked="" type="radio"/> 30014	Adjustable Tension Seat Upholstery	\$0
<input checked="" type="radio"/> 25060	Hook and Loop Positioning Strap	\$0
BACK UPHOLSTERY:		
<input checked="" type="radio"/> 30029	Fold - Over Back Upholstery	\$0
SEAT CUSHION:		
<input checked="" type="radio"/> 25063	Seat Cushion	\$0

STANDARD ACCESSORIES:		
<input checked="" type="radio"/> 25061	Safety Flag	\$0
<input checked="" type="radio"/> 25062	Chain Guard	\$0

UPGRADABLE OPTIONS AND ACCESSORIES:		
<i>(Accessories will not be installed and may be shipped separately - See Definitions)</i>		
HAND PEDAL UPGRADE OPTIONS: <i>(See definitions)</i>		
<input type="radio"/> 65128	Horizontal Hand Pedals	\$127
<input type="radio"/> 65061	Quad Twist Shifter Adaptor	\$181
<input type="radio"/> 65134	Tri-Pin Hand Pedals	\$361
<input type="radio"/> 65135	Quad Cuff Hand Pedals	\$481
<input type="radio"/> 65133	Quad-Grips Hand Pedals with Wrist Wraps & 3" Crank Width Extensions	\$962
CRANK KIT UPGRADE OPTIONS: <i>(See definitions)</i>		
<input type="radio"/> 65141	V - Crank Upgrade	\$360
<input type="radio"/> 65140	Mountain Drive	\$1,204
SAFETY:		
<input type="radio"/> 25015	Light, Front <i>(White)</i>	\$18
<input type="radio"/> 25016	Light, Rear <i>(Red)</i>	\$18
<input type="radio"/> 25097	Mirror, Rear View	\$37
HYDRATION:		
<input type="radio"/> 25021	Water Bottle and Cage	\$60
POSITIONING:		
<input type="radio"/> 25001	Click Strap, 8" Small	\$90
<input type="radio"/> 25002	Click Strap, 10" Medium	\$97
<input type="radio"/> 25003	Click Strap, 12" Large	\$109
PERFORMANCE ENHANCEMENTS:		
<input type="radio"/> 25023	Crank Width Extensions <i>(+1.5" of crank width per set)</i>	\$24
<input type="radio"/> 25024	Cordless Computer <i>(speed, distance, etc.)</i>	\$120
<input type="radio"/> 25022	Wheel Alignment Gauge	\$259
TRANSPORTATION:		
<input type="radio"/> 25006	Crutch Holder	\$181
<input type="radio"/> 25025	Tow Bar for Wheelchair	\$271

Invacare® Top End® Excelerator™ Stock Handcycle



HANDCYCLE SPECIFICATIONS:	
Weight:	50 pounds
Turning Radius:	12 feet for 360 degree circle, 8 feet on a 3-point turn
Transfer Height:	20" from ground with 2" cushion
Frame Style:	One piece main frame, upright position, 1-3/4" X .065" ovalized carbon steel tubing
Fork:	Carbon Steel
Propulsion / Shifting:	Internal front hub, twist activated shifter mounted on crank assembly
Cranks:	175mm (length) X 12" (width), adjustable height
Back:	Adjustable height from 14"-18"
Footrest:	Dual adjustable height and angle, includes straps
Hand Pedals:	Ovalized ergonomic with frictionless bearings
Brake:	Reverse Activated Drum and Locking Parking Brake mounted on fork
Wheelbase:	40" Axle to Axle
Camber:	3 degree
Wheels:	24" metal spoked wheels
Axles:	Quick release
Tires / Tubes:	24" X 1.75" treaded tires (40 PSI)
Positioning Straps:	Hook and loop seat restraint included
Cushions:	Seat cushion included
Safety / Visibility:	Flag and wheel reflectors included
Chain Guard:	Full length, shock resistant composite material except when Mountain Drive is chosen. With Mountain Drive selection a half moon shaped black composite chainguard is used.

DEFINITIONS:	
Assembly and Shipping:	Arrives factory assembled, rear wheels and some accessories will need to be installed. Shipping rate is \$150 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
Quad Twist Shifter Adaptor:	Three pronged adaptor which allows user to tap for shifting.
Tri-pin Hand Pedals:	3-Point contact hand pedal similar to hand controls used for vehicles. Allows for the hands to be kept in place for people with limited grip strength.
Quad Grips Hand Pedals:	Combination of Tri-Pin hand pedal and Quad Cuff hand pedal. Allows for the hands to be kept in place for people with limited grip strength.
Quad Cuff Hand Pedals:	Includes Power Plates on the top and the bottom of the handpedal to assist through both up and down stroke. Includes an adjustable metal cuff with a hook and loop strap to secure hands.
Mountain Drive:	Additional seven lower gears for help climbing hills. Works similar to a transmission. Not available with V crank option.
V - Cranks:	Option for people who need more clearance due to long legs or large midsections. Cranks are 180mm long and 20" wide. Not available with Mountain Drive option.
Upgradable Options and Accessories:	All Add-on options and accessories may be boxed separately and require self installation. Additionally, all add-on options may be shipped separately.

Warranty Info:

Invacare warrants the frames when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Invacare warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty