



The **Travel Pro** from Pride® brings exceptional value and easy transportation together in a compact, lightweight package. The stylish Travel Pro has a weight capacity of 275 lbs. and disassembles quickly and easily.



Pride®
Mobility Products Corp.
— LIVE YOUR BEST —



Specifications

Max. weight capacity	275 lbs.
Maximum speed ¹	Up to 4 mph
Ground clearance ²	1.375"
Turning radius ²	32.5"
Overall length ²	37"
Overall width ^{2,6}	19.5"
Tires	
Front	2" x 8" solid
Rear	2.5" x 8" solid
Anti-tips	Rear
Standard seating	Foldable, compact, vinyl-covered molded plastic
Drivetrain	24-volt DC motor, sealed mini-transaxle
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	Off-board, 2 amp
Per-charge range ^{1,5}	Up to 8 miles (w/ 2 x 12 AH)
Battery requirements ³	(2) 12V deep-cycle (standard)
Compatible batteries	12 AH
Weight	
Batteries (2 required)	12 AH (19.5 lbs. battery pack) ⁴
Weight w/o battery pack ²	3-wheel: 72 lbs.
Heaviest piece ²	3-wheel: 28 lbs. (front section)
Standard Colors	



Red



Blue

Features

- Easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Sleek sporty style
- Modular design for easy serviceability
- Front frame-mounted seat post offers maximum stability
- PDAC coded K0800 power operated vehicle



Silver Star Lifts & Ramps

These Silver Star® Lifts are compatible with the Travel Pro.

- Backpacker® Plus
- Backpacker® MV
- Power Tote Mini



Full Platform



Silver Boom 300

(US) 800-800-8586 • (Canada) 888-570-1113 • www.pridemobility.com

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³ AGM or Gel-Cell type required.

⁴ Battery weight may vary based on manufacturer.

⁵ This range was obtained with the 4-wheel model.

⁶ Overall scooter width at seat is 22.3".

*The information contained herein is correct at the time of publication.
We reserve the right to alter specifications without prior notice.*