

Etac Turner Pro

Intuitive and easy to use sit-to-stand aid

oetac[®]

Etac Turner Pro

An easy-to-use turn aid that encourages natural standing and movement for the user and minimizes the risk of injury for the carer

Encourages natural movement

Etac Turner Pro is a turn aid with an intuitive design that offers safe patient turning with standing support. This turn aid enables a good posture and minimises the risk of injury for the carer when transferring.

Etac Turner Pro encourages natural standing and sitting movement for the user. The only adjustments needed is the height of the leg support and the handle. After that, the transfer can be carried out safely while the carer and the user maintain eye contact and communication.

The orange handle provides an intense, high contrasting colour which can help users with impaired vision or dementia. The small-scale footplate makes it the perfect choice of transfer aid when transferring in narrow spaces.

Etac Turner Pro is compact and easy to grip, and due to its low weight of only 8 kg, it is easy to bring along.



Product	Item no
Etac Turner Pro	1609010
Weight	
Weight of unit	81
Maximum user weight (SWL)	200 I
Dimensions	
Total width	450 m
Total depth	460 m
Handle height	800 - 1200 m
Footplate Ø	400 m
Foot plate height	20 m
Leg support width	380 m

Leg support height Leg support adjustable height min - max

150 mm tht min - max 300 - 550 mm

Materia

Powder-coated aluminium, plastic, steel.

Washing instructions









Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution.

The product can be safely cleaned in a cabinet washer-disinfector at 85°C for 3 minutes.

3

Features



Compact

Etac Turner Pro is easy to transport due to its low weight of only 7.5 kg. It is compact, evenly balanced and easy to grip.



Angled leg support

Soft and comfortable leg support with lateral support that follows the movement from sitting to standing. Easily adjusted to the required height.



Contrasting coloured handle

The orange handle provides a strong, high contrasting colour which can help users with impaired vision or dementia. It also offers many grip possibilities for both carer and



Functional and easy to move

The base plate in cast aluminum is very easy to turn and it has an integrated anti-slip function. It is also very easy to move indoors on different floors, including carpets thanks to transport wheels.



Etac Turner Pro encourages a natural standing and sitting movement for the user and provides a naturally good posture that minimizes the risk of injury to the carer during the transfer.





Founded in 1973, Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our ambition is to provide solutions that optimise quality of life for the individual, family members and caregivers.

For the latest news and continuously updated product information – visit: www.etac.com

Etac AB
Kista Science Tower
SE-164 51 Kista, Sweden
Tel +46 371 - 58 73 30
info@etac.se www.etac.com

