

TiLite TR Series 3 Order Form

Effective February 1, 2015 (New Items in Red)

TITANIUM RIGID FRAME

Page 1 of 8

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:		Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	State: Zip:
Ship Date:	Order #:	Phone:	Fax:

Please Email or Fax to the TiLite Customer Service Team at orders@tilite.com or 866.586.2413 / 509.586.2413



TRFS1 ☐ **TiLite TR Series 3** (Minimal Adjustability. 265 Pound Weight Limit.).....\$3,495

TIFIT®: TAILORED INDIVIDUALIZED FRAME

SEAT WIDTH A Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

TRSW1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....STD
 TRSW2 ☐ 20" (HCPCS code only applies to 20" Rear Seat Width A.).....\$450 E2201

SEAT DEPTH B Measure from front of back post to the front edge of seat sling. Front edge of seat sling will be .5" behind bend at front of chair frame, unless Glide™ Forks, or a 90° Front Angle is ordered, in which case front edge of seat sling begins approximately 1" behind the bend.

TRSD1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....STD
 TRSD2 ☐ 20".....\$475 E2203
 TRSD3 ☐ Custom Frame Depth: _____" Frame Depth Longer or Shorter Than
 Seat Depth (A Custom Frame Depth of 20" is \$475.).....NCO

FRONT SEAT HEIGHT C Measure from floor to top of seat tube at beginning of bend.

TRFSH1 ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5"
☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21".....STD

REAR SEAT HEIGHT D Measure from floor to top of seat tube at back post.

REAR SEAT HEIGHT RANGES

Rear Wheel Size	Rear Seat Height D		
	No Rear Suspension	Frog Legs® Rear Suspension	
20"	12.5" - 21"	15.0" - 21"	TRRSH1 <input type="checkbox"/> 12.5" <input type="checkbox"/> 13" <input type="checkbox"/> 13.5" <input type="checkbox"/> 14"
22"	13.5" - 21"	16.0" - 21"	<input type="checkbox"/> 14.5" <input type="checkbox"/> 15" <input type="checkbox"/> 15.5" <input type="checkbox"/> 16"
24"	14.5" - 21"	16.5" - 21"	<input type="checkbox"/> 16.5" <input type="checkbox"/> 17" <input type="checkbox"/> 17.5" <input type="checkbox"/> 18"
25"	15.0" - 21"	17.0" - 21"	<input type="checkbox"/> 18.5" <input type="checkbox"/> 19" <input type="checkbox"/> 19.5" <input type="checkbox"/> 20"
26"	15.5" - 21"	17.5" - 21"	<input type="checkbox"/> 20.5" <input type="checkbox"/> 21".....STD

SEAT TO FOOTREST B/FOOTRESTS Measure from front edge of seat sling to top rear of footrest.

Standard Footrests Measurement must be at least 2" shorter than Front Seat Height C.

☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17"
☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19"
 TRFTR3 ☐ Angle Adjustable FootrestSTD
 TRFTR1 ☐ Titanium Footrest w/Flat ABS Cover.....\$175
 TRFTR2 ☐ Titanium Open Loop Footrest\$225
 TRFTR12 ☐ Flip-Back Footrest (Not available with ABS Cover
 with 2" Sides or Front Roller.).....\$195 K0040, K0037




High Mount Footrest

☐ 7"-9" ☐ 9"-11" ☐ 11"-13" Set at _____ (Use 1/2" increments)
 TRFTR4 ☐ Angle Adjustable High Mount Footrest\$85 K0040, K0037
 TRFTR14 ☐ ABS Cover w/2" Sides (Not available with Flip-Back Footrest.).....\$95 K0041

SEAT BACK TYPE/PUSH HANDLES/SEAT BACK HEIGHT E Measure from top of back post to top of seat tube at rear of frame.

SEAT BACK TYPE

TRSB1 ☐ Folding Aluminum Adjustable Height.....NCO
 TRSB2 ☐ Folding Titanium Adjustable Height with New Rounded Rigidizer Bar
 (Titanium back will match Frame Finish. Not available with Bolt-On Push Handles.).....\$225
 TRSB15 ☐ Non-Folding, Fixed Height Titanium with New Rounded Rigidizer Bar
 (Titanium back will match Frame Finish. Not available with Push Handles.).....\$245



THE ULTIMATE RIDE.®

tilite.com
 2701 West Court Street | Pasco, Washington 99301
 509.586.6117 Tel | 800.545.2266 Toll-Free
 509.586.2413 Fax | 866.586.2413 Toll-Free Fax

THE TITANIUM EXPERTS.®

TiLite® OF0007_TR_Series3_0215RevU

TiLite TR Series 3 Order Form

Effective February 1, 2015 (New Items in Red)

TITANIUM RIGID FRAME

Page 2 of 8


TiFIT®: TAILORED INDIVIDUALIZED FRAME

TISHAFT™ BACK RELEASE BARS (Not available with Non-Folding, Fixed Height Titanium Back.)

- TRSB3 ☐ AluminumSTD
 TRSB4 ☐ Titanium.....\$65
 TRSB5 ☐ Carbon Fiber.....\$65

PUSH HANDLES

- TRSB6 ☐ Integrated Push Handles.....\$85
 TRSB7 ☐ Integrated Fold-Down Push Handles (Requires Tension Adjustable by Straps Back Upholstery.).....\$100
 TRSB8 ☐ Bolt-On Push Handles (Includes 4" Deep Backrest Rigidizer Bar. Not available with Ti Back.).....\$175
 TRSB9 ☐ Bolt-On Height-Adjustable Push Handles (Includes 4" Deep Backrest Rigidizer Bar. Not available with Ti Back.).....\$295



FOLDING SEAT BACK HEIGHT ① (TRSB1 and TRSB2)

Push Handles (Selected Above)	Short	Medium	Tall	Set at: (Use 1/2" increments)
None / Bolt-On	<input type="checkbox"/> 8.5"- 11.5"	<input type="checkbox"/> 12"- 16.5"	<input type="checkbox"/> 16.5"- 20.5"	Set at: _____
Integrated	<input type="checkbox"/> 10.5"- 12.5"	<input type="checkbox"/> 13.5"- 16.5"	<input type="checkbox"/> 17.5"- 20.5"	Set at: _____
Integrated Fold-Down	<input type="checkbox"/> 13"- 15.5"	<input type="checkbox"/> 14.0"- 18.5"	<input type="checkbox"/> 18.5"- 21"	Set at: _____

NON-FOLDING, FIXED-HEIGHT ② TITANIUM SEAT BACK with New Rounded Rigidizer Bar (TRSB15) (Titanium back will match Frame Finish. Required for fixed Side Guards. Not available with Swing Away Tubular or Swing Away - Flip Back Tubular Armrests and Integrated Push Handles, Integrated Fold-Down Push Handles or Bolt-On Height Adjustable Push Handles.)

- ☐ 5" ☐ 5.5" ☐ 6" ☐ 6.5" ☐ 7" ☐ 7.5" ☐ 8" ☐ 8.5" ☐ 9" ☐ 9.5" ☐ 10"
☐ 10.5" ☐ 11" ☐ 11.5" ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5"
☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"

AVAILABLE ON ALL BACKS

- TRSB11 ☐ Backrest that Locks when Folded (Not available with Ergonomic Seat.).....NCO
 TRSB12 ☐ 4" Deep Backrest Rigidizer Bar (Available with Aluminum and Titanium backrest option. 4" Deep Rigidizer Bar is NOT a rounded style rigidizer bar.).....NCO
 TRSB13 ☐ Custom Height Backrest Rigidizer Bar (The Standard Rigidizer Bar height is 5" above the Seat Sling for a short, 6" for a medium and 10" for a tall. Custom Height Rigidizer Bar must be at least 4-1/2" above Seat Sling. Maximum Height is 3" less than Lowest Seat Back Height without push handles.).....\$150
 Set at: _____ (Use 1/4" increments)

SEAT BACK ANGLE ③ Measure angle from front of back post to floor.


- TRBSA2 **Folding Back**
☐ 80° ☐ 83° ☐ 86° ☐ 89° ☐ 92° ☐ 95° ☐ 98° ☐ 101°.....NCO
 TRBSA3 **Non-Folding Back**
☐ 80° ☐ 81° ☐ 82° ☐ 83° ☐ 84° ☐ 85° ☐ 86° ☐ 87° ☐ 88° ☐ 89° ☐ 90° ☐ 91°
☐ 92° ☐ 93° ☐ 94° ☐ 95° ☐ 96° ☐ 97° ☐ 98° ☐ 99° ☐ 100° ☐ 101°.....NCO

CENTER OF GRAVITY ④ Measure from front of back post at the seat tube to the center of the rear axle. If no center of gravity is selected, chair will ship with the center of gravity in the rearward-most position.

- TRCOG1 ☐ 1" ☐ 1.25" ☐ 1.5" ☐ 1.75" ☐ 2" ☐ 2.25" ☐ 2.5" ☐ 2.75" ☐ 3" ☐ 3.25" ☐ 3.5" ☐ 3.75"
☐ 4" ☐ 4.25" ☐ 4.5" ☐ 4.75" ☐ 5" ☐ 5.25" ☐ 5.5" ☐ 5.75" ☐ 6".....STD
 TRAMP1 ☐ Amputee Adapter (Black anodize only.)\$175 E0959

FRONT SEAT WIDTH (SEAT TAPER) ⑤ Measure from inside of front frame tube to inside of opposite front frame tube. (Required with V-Front End ⑥. Not Available with Standard Front End ⑦.)

- TRFSW1 ☐ No Taper (2" Narrower than Rear Seat Width ⑧.).....STD
 TRFSW2 ☐ 1" Taper (3" Narrower than Rear Seat Width ⑧. Includes Tension Adjustable by Straps Seat Upholstery.).....\$250
 TRFSW4 ☐ 2" Taper (4" Narrower than Rear Seat Width ⑧. Includes Tension Adjustable by Straps Seat Upholstery.).....\$250
 TRFSW6 ☐ 3" Taper (5" Narrower than Rear Seat Width ⑧. Includes Tension Adjustable by Straps Seat Upholstery. Not available with 12" Seat Width.).....\$250



TiLite TR Series 3 Order Form

Effective February 1, 2015 (New Items in Red)

TITANIUM RIGID FRAME

Page 3 of 8

TiFIT®: TAILORED INDIVIDUALIZED FRAME

FOOTREST WIDTH ❶ or ❷ Select either a Standard Front End ❶ if no Seat Taper ❶ is selected on page 2 or a V-Front End ❷ if a Seat Taper ❶ is selected on page 2.

TRFE1 STANDARD FRONT END ❶ Find your Rear Seat Width ❸ in the Chart below and select one of the options to the right. Measure from inside of front frame tube to inside of opposite front frame tube.**STD**



REAR SEAT WIDTH ❸ (Outside Measurement)	FOOTREST WIDTH ❶		
	No Taper (2" Narrower than Rear Seat Width)	2" Taper (4" Narrower than Rear Seat Width)	4" Taper (6" Narrower than Rear Seat Width)
12"	<input type="checkbox"/> 10"	<input type="checkbox"/> 8"	
13"	<input type="checkbox"/> 11"	<input type="checkbox"/> 9"	
14"	<input type="checkbox"/> 12"	<input type="checkbox"/> 10"	<input type="checkbox"/> 8"
15"	<input type="checkbox"/> 13"	<input type="checkbox"/> 11"	<input type="checkbox"/> 9"
16"	<input type="checkbox"/> 14"	<input type="checkbox"/> 12"	<input type="checkbox"/> 10"
17"	<input type="checkbox"/> 15"	<input type="checkbox"/> 13"	<input type="checkbox"/> 11"
18"	<input type="checkbox"/> 16"	<input type="checkbox"/> 14"	<input type="checkbox"/> 12"
19"	<input type="checkbox"/> 17"	<input type="checkbox"/> 15"	<input type="checkbox"/> 13"
20"	<input type="checkbox"/> 18"	<input type="checkbox"/> 16"	<input type="checkbox"/> 14"

OR

TRFE2 V-FRONT END ❷ Measure from inside of front frame tube to inside of opposite front frame tube 2-1/2" above footrest. (Titanium Footrests have 1" of adjustment up and down.)....**\$250**

☐ **NO V-Front Taper**

☐ 8" ☐ 8.5" ☐ 9" ☐ 9.5" ☐ 10" ☐ 10.5" ☐ 11" ☐ 11.5"
☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5"
☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5"



REAR WHEEL SPACING ❹/PLUGS Measure from outside of seat back post to inside of rear tire.
.5" Rear Wheel Spacing not available with 0° or 2° Camber ❺.

TRRWS1 ☐ .5" ☐ .75" ☐ 1" ☐ 1.25" ☐ 1.5" ☐ 1.75" ☐ 2".....**STD**

TRCBR12 ☐ Adjustable Rear Wheel Spacing Camber Plugs (Not available with Power Adaptable Frame.).....**STD**

TRCBR9 ☐ Additional Spacers to Increase Rear Wheel Spacing.....**NCO**

TRCBR11 ☐ Zero Play Aluminum Camber Plugs (Not available with Power Adaptable Frame or more than 6° of Camber ❺.).....**NCO**

CAMBER ❺/CAMBER TUBES

TRCBR1 ☐ 0° ☐ 2° ☐ 4° ☐ 6° ☐ 8° ☐ 12° (More than 6° of Camber ❺ is not available with Zero Play Aluminum Camber Plugs or Carbon Fiber Camber Tube.).....**STD**

TRCBR15 ☐ Aluminum Camber Tube.....**STD**

TRCBR3 ☐ Titanium Camber Tube (Not available with Power Adaptable Frame.).....**\$175**


TRCBR6 ☐ Carbon Fiber Camber Tube (Not available with Power Adaptable Frame or more than 6° of Camber ❺.).....**\$175**

TRCBR13 ☐ Camber for Additional Camber Tube ☐ 0° ☐ 2° ☐ 4° ☐ 6° ☐ 8° ☐ 12° (More than 2° of Camber ❺ is not available with the Power Adaptable Frame. More than 6° of Camber ❺ is not available with Zero Play Aluminum Camber Plugs or Carbon Fiber Camber Tube.).....**STD**

TRCBR5 ☐ Additional Aluminum Camber Tube.....**\$100**

TRCBR4 ☐ Additional Titanium Camber Tube (Not available with Power Adaptable Frame.).....**\$225**

TRCBR7 ☐ Additional Carbon Fiber Camber Tube (Not available with Power Adaptable Frame or more than 6° of Camber ❺.).....**\$225**



FRONT ANGLE ❶ Measure from frame front to floor.

TRFA1 ☐ 60° ☐ 65° ☐ 70° ☐ 75° ☐ 80° ☐ 85° (Not available with 6" Front Wheels.) ☐ 90° (Not available with 5" or 6" Front Wheels.)...**STD**

ERGONOMIC SEAT ❶ Measure from front of back post to beginning of upward bend of seat tube. (Not available with Tension Adjustable Bolt On Seat Upholstery or G10 Solid Seat.)

TRERG1 ☐ 4" ☐ 5" ☐ 6" ☐ 7" ☐ 8" ☐ 9" ☐ 10"**\$300**

FROG LEGS® REAR SUSPENSION (Not available with Power Adaptable Frame, User Friendly Flip-up, Titanium with Single Swivel Caster or Quick Release Titanium with Single Swivel Caster Anti-Tips.)

TRFRG2 ☐ Rider's Weight Under 180 lbs. (Black anodize only. Not available with Fendered Sideguards.).....**\$595**

TRFRG3 ☐ Rider's Weight Over 180 lbs. (Black anodize only. Not available with Fendered Sideguards.).....**\$595**

TiLite TR Series 3 Order Form

Effective February 1, 2015 (New Items in Red)

TITANIUM RIGID FRAME

Page 4 of 8

FRAME STYLE

TILITE SPECIAL EDITIONS If selecting one of the TiLite Special Editions you must still fill out the entire order form. Edition options will override inconsistent selections. Options cannot be added or deleted.

TRP1 ☐ SuperLite Edition

TRFTR2 Titanium Open Loop Footrest\$225

Seat Back Type

TRSB2 Folding Titanium Back.....\$225

TRSB5 Carbon Fiber TiShaft™ Back Release Bar (Not available with Non-Folding, Fixed Height Titanium Back.).....\$65

OR

TRSB15 Non-Folding, Fixed Height Titanium Back.....\$245

TRCBR6 Carbon Fiber Camber Tube (Not available with a Power Adaptable Frame. When a Power Adaptable Frame is selected, an Aluminum Camber Tube will be used.).....\$175

TRFK2 TiLite Slipstream Single-Sided Fork.....\$150

TRTFS1 Titanium Fork Stem.....\$50

Retail: \$845 - \$890

SuperLite Edition Price: \$755

TRP5 ☐ Black Edition (Black Edition includes: Black Lava Frame Color, Black Anodize Package, LiteSpeed Front Wheels (Must Choose

Front Wheel Size, Carbon Fiber TiShaft™ Back Release Bar and TiLite Slipstream Single-Sided Fork.)\$950

Exclusive To The Black Edition

☐ Black on Black Tattoo

Black Lava Footrest

Carbon Fiber Backrest Stripe with Black TiLite Logo

Stainless Quick Release Axles with Black TiLite Tattoo



TRPWR1 ☐ Power Adaptable Frame for Frank Mobility System (No charge when 20" Rear Seat Width **A** and 20" Seat Depth **B** are selected. Requires LiteSpeed Billet Aluminum Front Wheels and 0 degree camber. Not available with Carbon Fiber or Titanium Camber Tubes, Slipstream Forks.).....\$300

TRRRF1 ☐ Reinforced Frame (No charge when 20" Rear Seat Width **A** and 20" Seat Depth **B** are selected.)\$300

FRAME FINISHES

Natural Titanium

TRFF1 ☐ SatinSTD

TRSCBT1 ☐ Satin Finish Buffing Pads.....\$15

TRFF2 ☐ Polished\$650

TRFF3 ☐ Tattooed Titanium™ (Tattoos on front frame tubes. Frame finished in Satin. Click on the images below to view the Tattoos online.) ...\$295

TRFF31 ☐ The King of Rock

TRFF32 ☐ Aloha



TRFF33 ☐ Snow Leopard

TRFF34 ☐ Wood Grain



High Gloss Paint

TRFF8 ☐ Super Wet Black

TRFF10 ☐ TiLite Red

TRFF29 ☐ Neon Pink

TRFF28 ☐ Neon Green

Metallic Paint

TRFF14 ☐ Ocean Blue Metallic

TRFF15 ☐ Tangerine Metallic

TRFF16 ☐ Salsa Red Metallic

TRFF30 ☐ Acid Green Metallic

Pearlescent Paint

TRFF17 ☐ Midnight Blue Pearl

TRFF18 ☐ White Pearl

TRFF27 ☐ Titanium Grey

Matte Paint

TRFF19 ☐ Black Lava

TRFF20 ☐ Deep Blue Ice

TRFF25 ☐ Electric Plum.....STD

Candy Paint

TRFF22 ☐ Ruby Red

TRFF23 ☐ Azure Blue

COLOR ANODIZE PACKAGE

(Color Anodize Package includes: Aluminum Camber Tube (Color Anodize not available with Power Adaptable Camber Tube), Camber Mount, Footrest Clamp, Caster Mount, Caster Fork, and Caster Cap. If no anodize color is selected, the default will be black.)

TRANO1 ☐ Black.....STD

TRANO5 ☐ Sapphire TRANO3 ☐ Red Pepper TRANO6 ☐ Burnt Orange TRANO7 ☐ Sub-Lime TRANO8 ☐ Ultraviolet.....\$150

FRONT WHEELS

TRFW1 ☐ 3"x.75" Micro-Caster (Not available with Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.).....STD

TRFW2 ☐ 3"x.75" Light-Up Micro-Caster (Not available with Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.).....\$70

TRFW24 ☐ 3"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub.....\$165

TRFW14 ☐ 4"x.75" Plastic Wheel w/Poly Tire (Not available with Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.).....NCO

TRFW3 ☐ 4"x.75" Light-Up Micro-Caster (Not available with Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.).....\$70

TRFW25 ☐ 4"x.75" LiteSpeed Billet Aluminum Wheel w/Poly Tire ☐ Silver Hub OR ☐ Black Hub\$125

TRFW26 ☐ 4"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame.).....\$80

TRFW27 ☐ 4"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub.....\$175

TRFW6 ☐ 5"x1" Plastic Wheel w/Poly Tire (Not available with 90° Front Angle **1**, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.).....NCO

TRFW28 ☐ 5"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with 90° Front Angle **1**, or Power Adaptable Frame.).....\$80

TRFW29 ☐ 5"x.75" LiteSpeed Billet Aluminum Wheel w/Poly Tire ☐ Silver Hub OR ☐ Black Hub (Not available with 90°

Front Angle **1**.).....\$145

TRFW30 ☐ 5"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub (Not available with 90°

Front Angle **1**.).....\$185

TRFW9 ☐ 6"x1" Plastic Wheel w/Poly Tire (Not available with 85° or 90° Front Angle **1**, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.).....NCO





TiLite Series 3 Order Form

Effective February 1, 2015 (New Items in Red)

TITANIUM RIGID FRAME

Page 5 of 8










CUSTOM CONFIGURATION OPTIONS

TRFW34 <input type="checkbox"/> 6"x1" LiteSpeed Billet Aluminum Wheel w/Poly Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub (Not available with 85° or 90° Front Angle  , Slipstream Fork, Glide™ Single-Sided Fork.).....	\$165
TRFW31 <input type="checkbox"/> 6"x1.25" LiteSpeed Billet Aluminum Wheel w/Pneumatic Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub (Not available with 85° or 90° Front Angle  , Slipstream Fork, Glide™ Single-Sided Fork.).....	\$165
TRFW32 <input type="checkbox"/> 6"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with 85° or 90° Front Angle  , Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.).....	\$85
TRFW33 <input type="checkbox"/> 6"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub (Not available with 85° or 90° Front Angle  , Slipstream Fork, Glide™ Single-Sided Fork.).....	\$195

FRONT FORKS

TRFK1 <input type="checkbox"/> TiLite Standard Fork.....	STD	
TRFK2 <input type="checkbox"/> TiLite Slipstream Single-Sided Fork (Not available with Power Adaptable Frame.).....	\$150	
TRGL1 <input type="checkbox"/> Out-Front Glide™ Suspension Fork (Black anodize only.).....	\$425	E1015
Rider's Weight: <input type="checkbox"/> 125 lbs. <input type="checkbox"/> 125-180 lbs. <input type="checkbox"/> 180-265 lbs.		
TRGL2 <input type="checkbox"/> Out-Front Glide™ Single-Sided Suspension Fork (Black anodize only. Not available with 6" Front Wheels.).....	\$545	E1015
Rider's Weight: <input type="checkbox"/> 125 lbs. <input type="checkbox"/> 125-180 lbs. <input type="checkbox"/> 180-265 lbs.		
TRTFS1 <input type="checkbox"/> Titanium Fork Stem.....	\$50	

REAR WHEELS (See Rear Seat Height Ranges Chart on Page 1 for Rear Seat Height limitations.)

TRRW14 TiLite SHADOW™.....	STD	
<input type="checkbox"/> 20"(451)		
<input type="checkbox"/> 22"(501)		
<input type="checkbox"/> 24"(540)		
<input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth )		
<input type="checkbox"/> 26"(590) (Requires 15" minimum Seat Depth )		
TRRW13 Golz TWIN-STAR Exchange.....	\$995	
<input type="checkbox"/> 24"(540)		
<input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth )		
<input type="checkbox"/> 26"(590) (Requires 15" minimum Seat Depth )		
TRRW7 Spinerger SPOX <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Yellow <input type="checkbox"/> Blue <input type="checkbox"/> White	\$795	
<input type="checkbox"/> 22"(501) - Std. Hub		
<input type="checkbox"/> 24"(540) - Std. Hub		
<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth )		
<input type="checkbox"/> 26"(590) - Sport Hub (Requires 15" minimum Seat Depth  , Not available with Quad Release Axles.)		
TRRW8 Spinerger LX <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Yellow <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Go Green <input type="checkbox"/> Bright Orange <input type="checkbox"/> Flamingo Pink.....	\$995	
<input type="checkbox"/> 22"(501) - Std. Hub		
<input type="checkbox"/> 24"(540) - Std. Hub		
<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth )		
<input type="checkbox"/> 26"(590) - Sport Hub (Requires 15" minimum Seat Depth  , Not available with Quad Release Axles.)		
TRRW15 Topolino CarbonCore (Not available with Solid Treaded Tires.).....	\$1,195	
<input type="checkbox"/> 24"(540)		
<input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth )		
TRRW11 Mag (Not available with Short Tabs, Quad Release Axles or Spoke Guards.).....	\$50	
<input type="checkbox"/> 22"(501)		
<input type="checkbox"/> 24"(540)		

AXLES

TRAXL1 <input type="checkbox"/> Stainless Quick Release.....	STD	
TRAXL2 <input type="checkbox"/> Titanium Quick Release.....	\$225	
TRAXL3 <input type="checkbox"/> Quad Release (Not available with Mag or 26"(590) Spinerger Spox or LX Rear Wheels.).....	\$135	

PNEUMATIC REAR TIRES

TRRTR1 <input type="checkbox"/> Treaded (Not available with 25"(559) Rear Wheels.).....	STD	
TRRTR2 <input type="checkbox"/> Airless Inserts (Not available with Topolino Rear Wheels.).....	\$80	E2213
TRRTR3 <input type="checkbox"/> 25"(559) Primo Treaded (Standard with 25"(559) Rear Wheels.).....	STD	
TRRTR4 <input type="checkbox"/> Primo (Not available with Mag Rear Wheels.).....	NCO	
TRRTR16 <input type="checkbox"/> 20"(451) Kenda® Grey High Pressure (110psi).....	\$95	
TRRTR5 24"(540), 25"(559) and 26"(590) Kenda® High Pressure (110psi) <input type="checkbox"/> Grey <input type="checkbox"/> Red <input type="checkbox"/> Blue	\$95	
TRRTR6 <input type="checkbox"/> 24"(540) Grey Knobby (Available with Shadow Rear Wheels Only.).....	\$95	
TRRTR7 <input type="checkbox"/> 26"(590) Black Knobby (Available with Shadow, Spinerger SPOX, or Spinerger LX Rear Wheels Only.).....	\$95	
TRRTR8 <input type="checkbox"/> 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution (Not available with Mag Rear Wheels.).....	\$125	
TRRTR9 <input type="checkbox"/> 24"(540) and 25"(559) Schwalbe® Slick Tires (Not available with Mag Rear Wheels.).....	\$125	

SOLID REAR TIRES

TRRTR10 <input type="checkbox"/> 20"(451) Street Tread.....	\$65	
TRRTR11 <input type="checkbox"/> 22"(501) Pyramid.....	\$75	
TRRTR18 <input type="checkbox"/> 24"(540) and 25"(559) Herringbone Tread (Not available with Topolino Rear Wheels.).....	\$115	
TRRTR15 <input type="checkbox"/> 24"(540) and 25"(559) Treaded (Not available with Topolino Rear Wheels.).....	\$95	

TiLite TR Series 3 Order Form

Effective February 1, 2015 (New Items in Red)

TITANIUM RIGID FRAME

Page 6 of 8

CUSTOM CONFIGURATION OPTIONS

HANDRIMS

Handrims by OUT-FRONT®

TRRIM1 <input type="checkbox"/> Aluminum – Silver Anodized.....	STD
TRRIM21 <input type="checkbox"/> Aluminum – Black Anodized (Not available with Mag or 20"(451) Rear Wheels.).....	NCO
TRRIM12 <input type="checkbox"/> Titanium – Satin Finish (Not available with Mag, 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.).....	\$400
The Natural-Fit® (Original, Larger Profile) (Not available with 22"(501) Mag Rear Wheels.)	
TRRIM14 <input type="checkbox"/> Standard Grip (Not available with Short Tabs.).....	\$350
TRRIM15 <input type="checkbox"/> Super Grip (Not available with Short Tabs.).....	\$395
TRRIM16 <input type="checkbox"/> Without Thumb Grip.....	\$325
The Natural-Fit® LT (Reduced Weight, Smaller Profile) (Not available with 22"(501) Mag Rear Wheels.)	
TRRIM17 <input type="checkbox"/> Standard Grip (Not available with Short Tabs.).....	\$380
TRRIM18 <input type="checkbox"/> Super Grip (Not available with Short Tabs.).....	\$425
TRRIM19 <input type="checkbox"/> Without Thumb Grip.....	\$355
TRRIM10 <input type="checkbox"/> The Surge (All-In-One Oval with Gription™ Strip, Larger Profile) (Not available with 22"(501) Mag Rear Wheels.).....	\$350
TRRIM11 <input type="checkbox"/> The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile) (Not available with 22"(501) Mag Rear Wheels.).....	\$380
TRRIM20 <input type="checkbox"/> Q-Grip (Exclusive High-Friction, Anti-Slip Coating) (Not available with 20"(451), or Mag Rear Wheels.).....	\$295

Additional Handrims

TRRIM9 <input type="checkbox"/> Spinergy FlexRim (Available only on 24"(540) and 25"(559) Spinergy LX Rear Wheels. Please allow an additional 2-3 weeks delivery.).....	\$600
TRRIM5 <input type="checkbox"/> Plastic Coated – Black.....	\$125
TRRIM6 <input type="checkbox"/> Projections - 8 Vertical Projections (Not available with 20"(451), 22"(501), 26"(590) or Mag Rear Wheels.).....	\$325

Tab Length

TRRIM7 <input type="checkbox"/> Long Tabs.....	STD
TRRIM8 <input type="checkbox"/> Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.).....	\$20

WHEEL LOCKS

Wheel Locks by OUT-FRONT®

TRWLK7 <input type="checkbox"/> Compact Composite Push to Lock (Available with .75" to 1.5" Rear Wheel Spacing ① only.).....	NCO
TRWLK8 <input type="checkbox"/> Compact Composite Pull to Lock (Available with .75" to 1.5" Rear Wheel Spacing ① only.).....	NCO
TRWLK10 <input type="checkbox"/> Extension Handles for Compact Composite Push/Pull to Lock.....	\$75 E0961
TRWLK9 <input type="checkbox"/> Composite Scissor Lock.....	\$125

Additional Wheel Locks

TRWLK5 <input type="checkbox"/> Black Aluminum Scissor Lock.....	\$150
TRWLK3 Uni-Lock Set At: <input type="checkbox"/> Push OR <input type="checkbox"/> Pull.....	\$145
TRWLK4 <input type="checkbox"/> Extension Handles for Uni-Lock.....	\$75 E0961

BACK UPHOLSTERY

TRBUP2 <input type="checkbox"/> Velcro Adjustable (Not available with Fold-Down Push Handles.).....	STD
Tension Adjustable by Straps (Required with Fold-Down Push Handles.)	
TRBUP3 <input type="checkbox"/> Padded Nylon.....	\$300 E2611
TRBUP5 <input type="checkbox"/> AIR Breathable Upholstery (Not available with Ultrasuede or NiteLite Backrest Stripe.).....	\$300 E2611
TRUPC3 Ultrasuede Backrest Stripe <input type="checkbox"/> Black <input type="checkbox"/> Sapphire <input type="checkbox"/> Red Pepper <input type="checkbox"/> Burnt Orange <input type="checkbox"/> Sub-Lime	
<input type="checkbox"/> Ultraviolet <input type="checkbox"/> Espresso.....	\$90
TRUPC3 <input type="checkbox"/> TiLite NiteLite Reflective Backrest Stripe.....	\$90
TRBUP4 <input type="checkbox"/> Ride Designs® Corbac™ Adjustable Contour Back Support (Black only. Not available with Fold-Down Push Handles. See Chart on Page 7.).....	\$455 E2611
TRBUP0 <input type="checkbox"/> Omit Back Upholstery.....	NCO

SEAT UPHOLSTERY

TRSUP1 <input type="checkbox"/> Tension Adjustable Bolt On (Not available with Ergonomic or Tapered Seat ①. Black only.).....	STD
TRSUP7 iSleeve™ Under Seat Smartphone Holder (Not available with 12" or 13" Rear Seat Width ① or Seat Depth ②.).....	STD
TRSUP2 <input type="checkbox"/> Tension Adjustable by Straps (Standard with Ergonomic Seat ① and Front Seat Width (Seat Taper) ① of 1", 2" or 3". Black upholstery only.).....	\$150
TRSUP3 <input type="checkbox"/> G10 Solid Seat (Not available with 10" or 11" Rear Seat Width ① or Ergonomic Seat ①.).....	\$245 E2618
TRSUP0 <input type="checkbox"/> Omit Seat Upholstery.....	NCO

SEAT CUSHIONS

TRCSH1 Foam Cushion <input type="checkbox"/> 2" <input type="checkbox"/> 3" <input type="checkbox"/> 4".....	\$99
TRCSH2 <input type="checkbox"/> StimuLITE® Sport (2" Black 12"x12" to 20"x20". Please allow an additional 2-3 weeks delivery.).....	\$525 E2603
TRCSH3 <input type="checkbox"/> Corbee® (2" Black 6"x6" to 16"x16". Please allow an additional 2-3 weeks delivery.).....	\$525 E2603
TRCSH4 <input type="checkbox"/> ROHO® Low Profile® Quadro® Select (Black only. Please allow an additional 2-3 weeks delivery.).....	\$525 E2624
TRCSH6 <input type="checkbox"/> VARILITE® Evolution™ (Black only. See Chart below.).....	\$525 E2622
TRCSH7 <input type="checkbox"/> VARILITE® Zoid™ (Black only. See Chart below.).....	\$525 E2622

VARILITE® SEAT CUSHION AVAILABILITY CHART

Seat Width ①	VARILITE® Evolution™		VARILITE® Zoid™	
	Seat Depth ②			
12"	14"	8", 12", 14"		
13"	NONE	NONE		
14"	14", 16", 18", 20"	14", 16", 18"		
15"	15", 17"	15", 17"		
16"	14", 16", 18", 20"	16", 18", 20"		
17"	17"	17"		
18"	16", 18", 20"	16", 18", 20"		
19"	NONE	NONE		
20"	16", 18", 20"	16", 18", 20"		

CORBAC™ SEAT WIDTH/BACK HEIGHT AVAILABILITY CHART

Rear Seat Width ①	Seat Back Height ②
14"	13" - 15"
15"	13" - 16"
16"	14" - 16"
17"	14" - 16"
18"	14" - 16"

TiLite Series 3 Order Form

Effective February 1, 2015 (New Items in Red)

TITANIUM RIGID FRAME

Page 7 of 8

CUSTOM CONFIGURATION OPTIONS

SIDE GUARDS

TRSDG1 ☐ Cloth (Not available with Ride Designs Corbac Back Support.).....\$85

Aluminum

TRSDG2 ☐ Rigid Removable ☐ Small **OR** ☐ Std.....\$175

TRSDG3 ☐ Quad Removable w/Hole ☐ Small **OR** ☐ Std.....\$185

TRSDG4 ☐ Fold Down Rigid ☐ Small **OR** ☐ Std (Not available with Ergonomic Seat.).....\$295

TRSDG6 ☐ Fendered Adjustable (Available with 24" and 25" Rear Wheels. Requires Rear Seat Height **D** of 14.5" to 18.5")......\$295

Carbon Fiber

TRSDG7 ☐ Rigid Removable ☐ Small **OR** ☐ Std.....\$270

TRSDG8 ☐ Quad Removable w/Hole ☐ Small **OR** ☐ Std.....\$280

TRSDG9 ☐ Fold Down Rigid ☐ Small **OR** ☐ Std (Not available with Ergonomic Seat.).....\$400

TRSDG10 ☐ Fixed (Requires Non-Folding, Fixed-Height Titanium Seat Back.).....\$375

TRSDG11 ☐ Fendered Adjustable (Available with 24" and 25" Rear Wheels. Requires Rear Seat Height **D** of 14.5" to 18.5")......\$495

TRSDG12 ☐ Fendered Removable.....\$395

Available on Aluminum and Carbon Fiber Non-Fendered Side Guards

TRSDG13 ☐ Side Guards Custom Cut to Profile of Rear Wheels (Not available with Cloth.).....\$275

ARMRESTS (Desk Arms come with built-in Side Guards.)

TRARM1 Swing Away Tubular ☐ 9" to 12" **OR** ☐ 11" to 14" (Not Available with Non-Folding, Fixed Height Titanium Seat Back or Height Adjustable Push Handles.).....\$195 **E0973**

TRARM4 Swing Away - Flip Back Tubular ☐ 6" to 9" **OR** ☐ 9" to 12" (Not Available with Non-Folding, Fixed Height Titanium Seat Back.).....\$275 **E0973**

TRARM2 Removable Height Adjustable Desk Arm with Rigid Side Guard (Requires 14" minimum Seat Depth **B**. Not available with Ergonomic Seat.).....\$295 **E0973**

☐ 9" to 12" **OR** ☐ 12" to 15"

☐ Standard Length Pad (10") **OR** ☐ Full Length Pad (14")

TRARM3 Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 15" minimum Seat Depth **B**. Not Available with Ergonomic Seat.).....\$325 **E0973**

☐ 9" to 12" **OR** ☐ 12" to 15"

☐ Standard Length Pad (10") **OR** ☐ Full Length Pad (14")

Available on all Armrests

TRARM5 ☐ Covered in Black Ultrasuede.....\$95

SEAT BELTS

TRBLT1 ☐ Velcro Adjustable, Black Nylon Webbing.....\$45 **E0978**

TRBLT2 ☐ Bodypoint® 2" Push-Button Buckle (Auto-Style.).....\$85 **E0978**

TRBLT4 ☐ Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" or 20" Rear Seat Widths **A**.).....\$165 **E0978**

TRBLT5 ☐ Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" or 20" Rear Seat Widths **A**.).....\$155 **E0978**

TRBLT6 ☐ Bodypoint® 1.5" Wide Rehab Latch Buckle Padded (Airline-Style.) (Not available with 19" or 20" Rear Seat Widths **A**.).....\$160 **E0978**

TRBLT7 ☐ Bodypoint® 1.5" Wide Plastic Side Release Padded (Not available with 19" or 20" Rear Seat Widths **A**.).....\$120 **E0978**

TRBLT9 ☐ Bodypoint® Evoflex Pelvic Stabilizer.....\$210 **E0978**

TRBLT10 ☐ Bodypoint® Optional Band Clamp Mounting Hardware.....\$20 **E0978**

CALF STRAP

TRCLF2 Velcro Adjustable.....\$35 **K0038**

☐ Xtra Small (Recommended for 6" - 9" Footrest Widths.)

☐ Small (Recommended for 9" - 12" Footrest Widths.)

☐ Medium (Recommended for 11" - 15" Footrest Widths.)

☐ Large (Recommended for 14" - 17" Footrest Widths.)

☐ Xtra Large (Recommended for 16"- 19" Footrest Widths.)

TRCLF3 Bodypoint® Padded Velcro Adjustable.....\$90 **K0038**

☐ Medium (Recommended for 11"- 15" Footrest Widths.)

☐ Large (Recommended for 16"- 19" Footrest Widths.)

TRCLF1 ☐ Buckle Adjustable, Black Nylon Webbing.....\$25 **K0038**

BACK PACK

TRBCK1 ☐ Small - Black Nylon.....\$50

TRBCK2 ☐ Large - Black Nylon.....\$65

SEAT POUCH

TRPCH1 ☐ Small - Black Nylon.....\$40

TRPCH2 ☐ Large - Black Nylon.....\$45

SPOKE GUARDS by SPOZEE.COM

TRGRD1 ☐ Clear Lexan with TiMotion Logo (Not available with Mag Rear Wheels.).....\$115 **K0065**

TRGRD2 ☐ Scooter Frog (Requires Shadow Rear Wheels.).....\$275 **K0065**

TRGRD3 ☐ Color Spectrum (ROYGBIV) on Black (Requires Shadow Rear Wheels.).....\$275 **K0065**

TRGRD4 ☐ Purple Party (Requires Shadow Rear Wheels.).....\$275 **K0065**

ANTI-TIPS

TRTIP1 ☐ Rear Aluminum - Black.....\$140 **E0971**

TRTIP2 ☐ User Friendly Flip-up (Not available with Frog Legs® Rear Suspension.).....\$295 **E0971**

TiLite TR Series 3 Order Form

Effective February 1, 2015 (New Items in Red)

TITANIUM RIGID FRAME

Page 8 of 8


CUSTOM CONFIGURATION OPTIONS

ACCESSORIES


TRUEST1 Upper Extremity Support Tray <input type="checkbox"/> 12" to 14" OR <input type="checkbox"/> 13" to 15" Seat Widths A (Requires Armrests.).....	\$250	E0950
TRGS1 <input type="checkbox"/> GRIP Solutions Lap Tray.....	\$85	
TRLUG1 <input type="checkbox"/> Luggage Carriers.....	\$125	
TRMPCT2 Neoprene Impact Guards with Ultrasuede Stripe <input type="checkbox"/> Black <input type="checkbox"/> Sapphire <input type="checkbox"/> Red Pepper <input type="checkbox"/> Burnt Orange <input type="checkbox"/> Sub-Lime <input type="checkbox"/> Ultraviolet <input type="checkbox"/> Espresso.....	\$90	
TRMPCT2 <input type="checkbox"/> TiLite NiteLite Neoprene Impact Guards with Reflective Stripe.....	\$90	
TRCAT1 <input type="checkbox"/> Catch-All.....	\$75	
TRCH1 <input type="checkbox"/> Crutch Holder.....	\$195	E2207
TRTIP5 <input type="checkbox"/> Front Roller (Requires that Seat to Footrest Measurement B be at least 3" shorter than Front Seat Height C . Not available w/ Titanium Flip-Back Footrest.).....	\$95	
TRHS1 <input type="checkbox"/> TiLite 9 Piece SAE Round Handle Folding Hex Key Set.....	\$25	
TRAFP1 4" Logo Patch <input type="checkbox"/> Army <input type="checkbox"/> Navy <input type="checkbox"/> Air Force <input type="checkbox"/> Marines.....	\$90	
TRAFP2 Name Tape Patch <input type="checkbox"/> Army <input type="checkbox"/> Navy <input type="checkbox"/> Air Force <input type="checkbox"/> Marines.....	\$90	
TRAFP3 <input type="checkbox"/> Combat Wounded Purple Heart Patch.....	\$90	

FRAME NOTES

The rear frame will be square when the Center of Gravity **H** measurement is less than 3.5".



A. Round Rear Frame



B. Square Rear Frame

DESIGNS UNLIMITED

TiLite's TR Series 3 offers more than 470 billion different frame combinations! But when that's not enough, check out the Designs Unlimited Studio at tilite.com. There you will find drawings and pricing for hundreds of creative solutions to customers' customization requests. We call these "D.U.'s." If you don't find the D.U. you're looking for, contact the TiLite Customer Service Team at 800.545.2266 or customerservice@tilite.com and we'll design a D.U. just for you. Designs Unlimited customization may add an additional 2 to 3 weeks to the ship time of your chair.

GENERAL

TiLite manufactures a wide variety of wheelchairs to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. TiLite wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at tilite.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, The Titanium Experts, and TiShaft are registered trademarks of TiSport, LLC and are licensed to TiLite, LLC. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. Ride Designs® and Corbac® are registered trademarks of Aspen Seating, LLC. ROHO®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. Supracor®, Stimulite® and Corbee® are registered trademarks of Supracor, Inc. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. Varilite® is a registered trademark of Cascade Designs, Inc.

QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to TiLite dealers and distributors). A TiLite Quote or a TiLite Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR