

# Top End® Force™ 3 Disc Brake Enabled Handcycle



PROVIDER / CUSTOMER INFORMATION:		
Provider Information:		End User Information:
Date:	Contact Name:	Athlete Name:
Acct #:	Ship to Address:	
Contact:	City:	State:
Phone #:	Zip/Postal Code:	Phone:
Fax/Email:	Email:	
<input type="radio"/> Ship to Provider		<input type="radio"/> Ship to End User
ORDER REQUEST TYPE:		
<input type="radio"/> Quote		<input type="radio"/> Order
		PO # <input type="text"/>

BASE MODEL:		
<input checked="" type="radio"/> 10040	Force 3 Disc Brake Enabled Handcycle	\$5,499
Assembly: <i>(see definitions)</i>		
<input checked="" type="radio"/> 10048	Assembled	\$0
Shipping:		
Shipping Fee <i>(Continental US Only - See definitions)</i>		\$325

DISC BRAKE:		
<input checked="" type="radio"/> 65081	Front Disc Brakes	\$0

BACKREST TYPE:		
<input checked="" type="radio"/> 15172	Straight, Upright position, Adjustable Angle Backrest	\$0

FORK TYPE:		
<input checked="" type="radio"/> 75026	Wide Fork	\$0

DRIVE CONTROL:		
CRANK ASSEMBLY AND FORK:		
<input checked="" type="radio"/> 65179	Length and Height Adjustable Fork	\$0
CRANK ARM TYPE:		
<input checked="" type="radio"/> 65048	V-Cranks, 165mm Long x 14" Wide	\$0

HAND PEDALS:		
<input checked="" type="radio"/> 65137	Medium 3D Printed W/Power Plates	

SHIFTING:		
<input checked="" type="radio"/> 65230	Single Chainring with 10sp Shifter	\$0
CHAINRING:		
<input checked="" type="radio"/> 65104	Chainring, Single 34T	\$0

CASSETTE:		
<input checked="" type="radio"/> 65161	11-42T Front Wheel Cassette	\$0
BRAKING:		
<input checked="" type="radio"/> 65090	Locking Parking Brake	\$0

FOOTREST:		
<input checked="" type="radio"/> 60009	4-3/4" Footrest Width	\$0

WHEELBASE: <i>(Front to Rear Axle)</i>		
<input checked="" type="radio"/> 75078	60" Wheelbase	\$0
CAMBER ANGLE:		
<input checked="" type="radio"/> 45025	1 Degree Camber	\$0



WHEEL SETS:		
WHEEL SIZE:		
<input checked="" type="radio"/> 50171	Wheel Size, 26" Front, 20" Rear	\$0
WHEEL TYPE:		
<input checked="" type="radio"/> 50196	High Performance	\$0
TIRES:		
<input checked="" type="radio"/> 50248	Tires, 26" (559) & 20" (406)	\$0
<input checked="" type="radio"/> 50402	Threaded Axles	\$0

FRAME COLOR:		
<input checked="" type="radio"/> 55001	Red	Simple Red \$0

UPHOLSTERY:		
MATERIAL & COLOR:		
<input checked="" type="radio"/> 30001	Nylon, Black Upholstery	\$0
SEAT UPHOLSTERY:		
<input checked="" type="radio"/> 25060	Hook and Loop Positioning Strap	\$0
BACK UPHOLSTERY:		
<input checked="" type="radio"/> 25064	Back Cushion	\$0
SEAT CUSHION:		
<input checked="" type="radio"/> 25063	Seat Cushion	\$0

Top End Sports LLC.

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Top End unless otherwise noted.

Specifications and prices are subject to change without notification.

©2023 Top End Sports LLC.

60106218~M

Page 1 of 2

# Top End Force™ 3 Disc



## Brake Enabled Handcycle

ACCESSORIES:		
STANDARD:		
● 25061	Safety Flag	\$0
● 25072	Chain Guard	\$0
● 25067	Drafting Bumper	\$0
SAFETY:		
● 25041	Mirror, Footrest Mounted	\$0
HYDRATION:		
● 25102	Water Bottle and Cage	\$0

FORK COLOR		
○ TE55061	Fork Color: Agricultural Green	\$120
○ TE55042	Fork Color: Black Prism	\$120
○ TE55038	Fork Color: Bubblegum Pink	\$120
○ TE55035	Fork Color: Electric Blueberry	\$120
○ TE55031	Fork Color: Grasshopper Green	\$120
○ TE55063	Fork Color: Lilac Purple	\$120
○ TE55033	Fork Color: Midnight Blue	\$120
○ TE55062	Fork Color: Neon Green	\$120
○ TE55059	Fork Color: Neon Orange	\$120
○ TE55065	Fork Color: Neon Pink	\$120
○ TE55058	Fork Color: Neon Red	\$120
○ TE55060	Fork Color: Neon Yellow	\$120
○ TE55026	Fork Color: Raspberry Red	\$120
○ TE55034	Fork Color: Royal Blue	\$120
○ TE55030	Fork Color: Safety Yellow	\$120
○ TE55043	Fork Color: Silver Metallic	\$120
○ TE55028	Fork Color: Simple Orange	\$120
○ TE55024	Fork Color: Simple Red	\$120
○ TE55046	Fork Color: Stealth Grey	\$120
○ TE55029	Fork Color: Sunkissed Orange	\$120
○ TE55025	Fork Color: Tangerine Red	\$120
○ TE55041	Fork Color: Wet Black	\$120
○ TE55039	Fork Color: White Snow	\$120

HANDCYCLE SPECIFICATIONS:	
<b>Weight:</b>	40 Pounds
<b>Turning Radius:</b>	18 feet for 360 degree circle, 12 feet with a 3-point turn
<b>Transfer Height/ Ground Clearance:</b>	10" from ground with 2" cushion 5.5" Ground Clearance
<b>Frame Style:</b>	Wrap around ovalized lower frame, 7000 series aluminum
<b>Fork:</b>	Removable with adjustable height of 22.5" to 24.5
<b>Propulsion / Shifting:</b>	V-Cranks with 1x10 Speed Drivetrain
<b>Cranks:</b>	165mm x 14" Wide V-Cranks
<b>Back:</b>	Angle adjustable straight back 35-70 degrees
<b>Footrest:</b>	Dual adjustable length, 4-3/4" wide, includes straps
<b>Hand Pedals:</b>	Ergonomic 3D Printed handpedals with power plates
<b>Brake:</b>	Primary hand brake is a disc brake operated by a lever mounted on the right handpedal. Bike also has a locking parking rim caliper brake that is operated by a flip lever mounted on the fork.
<b>Wheelbase:</b>	60" Axle to axle
<b>Camber:</b>	1 Degree
<b>Wheels:</b>	559 metal spoked wheels
<b>Axles:</b>	Threaded
<b>Tires/Tubes</b>	26" (559) Front and 20" (406), Cruiser tires
<b>Positioning Straps:</b>	Hook and loop seat restraint included
<b>Cushions:</b>	Seat and back cushions included
<b>Safety / Visibility:</b>	Flag and wheel reflectors included

DEFINITIONS:	
<b>Force-3:</b>	The Force-3 has an updated and improved design. The frame is 17" wide with 15" seat and back
<b>Shipping:</b>	Shipping rate is \$325 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
<b>V Cranks:</b>	This crank option is great for customers whom require more clearance due to long legs or larger midsections. V Cranks also provides a more efficient cranking position by allowing users to use larger muscle groups.
<b>Straight Back:</b>	The straight back measures 22.5" tall.
<b>Upgrade Accessories:</b>	All Add-on options and accessories are boxed separately and require self installation.
<b>Force-3 Disc Brake Versions:</b>	The Force-3 Disc Brake version is equipped with a front disc brake for use as the primary hand brake. The parking brake is a cable driven caliper activated by a flip lever and is not intended to operate as standard brake.

### Warranty Info:

Top End warrants the frame and fork when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Top End warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty

# Top End Force™ 3 Disc Brake Enabled Handcycle



FRAME COLOR		
○ TE55105	Frame Color: Bomber Orange	\$0
○ TE55095	Frame Color: Carolina Blue	\$0
○ TE55093	Frame Color: Fighter Jet Grey	\$0
○ TE55090	Frame Color: Las Vegas Bronze	\$0
○ TE55092	Frame Color: Light Pink	\$0
○ TE55091	Frame Color: Mint Green	\$0
○ TE55094	Frame Color: Tiffany Blue	\$0
○ TE55053	Frame Color: Agricultural Green	\$0
○ TE55019	Frame Color: Black Prism	\$0
○ TE55015	Frame Color: Bubblegum Pink	\$0
○ TE55070	Frame Color: Buffed Aluminum - NC	\$0
○ TE55012	Frame Color: Electric Blueberry	\$0
○ TE55102	Frame Color: Flat Black	\$0
○ TE55008	Frame Color: Grasshopper Green	\$0
○ TE55055	Frame Color: Lilac Purple	\$0
○ TE55010	Frame Color: Midnight Blue	\$0
○ TE55054	Frame Color: Neon Green	\$0
○ TE55051	Frame Color: Neon Orange	\$0
○ TE55057	Frame Color: Neon Pink	\$0
○ TE55050	Frame Color: Neon Red	\$0
○ TE55052	Frame Color: Neon Yellow	\$0
○ TE55003	Frame Color: Raspberry Red	\$0
○ TE55011	Frame Color: Royal Blue	\$0
○ TE55007	Frame Color: Safety Yellow	\$0

FRAME COLOR		
○ TE55020	Frame Color: Silver Metallic	\$0
○ TE55005	Frame Color: Simple Orange	\$0
○ TE55001	Frame Color: Simple Red	\$0
○ TE55017	Frame Color: Stealth Black	\$0
○ TE55023	Frame Color: Stealth Grey	\$0
○ TE55006	Frame Color: Sunkissed Orange	\$0
○ TE55002	Frame Color: Tangerine Red	\$0
○ TE55018	Frame Color: Wet Black	\$0
○ TE55016	Frame Color: White Snow	\$0
○ TE55103	Frame Color: Wrinkle Black	\$0
○ TE25102	Water Bottle and Cage	\$0



Note: All Special requests are subject to approval, additional charges, and additional lead time

Special Request 1)	
Special Request 2)	
Special Request 3)	
Special Request 4)	
Special Request 5)	