

Top End® Force G™ Handcycle

PROVIDER / CUSTOMER INFORMATION:

Provider Information:		End User Information:	
Date:		Contact Name:	Athlete Name:
Acct #:		Ship to Address:	
Contact:		City:	State:
Phone #:		Zip/Postal Code:	Phone:
Fax/Email:		Email:	
<input type="radio"/> Ship to Provider		<input type="radio"/> Ship to End User	

ORDER REQUEST TYPE:

<input type="radio"/> Quote	<input type="radio"/> Order	PO # <input type="text"/>
-----------------------------	-----------------------------	---------------------------

BASE MODEL:

<input checked="" type="radio"/> 10012	Top End Force G Handcycle	\$6,999
--	---------------------------	---------

Assembly: (see definitions)

<input checked="" type="radio"/> 10048	Assembled	\$0
--	-----------	-----

Shipping:

Shipping Fee (Continental US Only - See definitions)	\$325
---	-------

SEAT WIDTH:

<input type="radio"/> 15015	15" Seat Width	\$0
<input type="radio"/> 15021	18" Seat Width	\$0

FOOTREST:

<input checked="" type="radio"/> 60009	4-3/4" Footrest Width	\$0
--	-----------------------	-----

FORK OPTIONS:

LOWER FORK STYLE:

<input checked="" type="radio"/> 75026	Wide Fork <i>Disc brake mounts will be welded on the frame</i>	\$0
--	---	-----

FORK ADJUSTABILITY:

<input checked="" type="radio"/> 65183	Adjustable Height Adjustable Length	\$0
--	--	-----

HAND PEDALS: (See Diagram Below)

HAND PEDAL SIZE:

<input checked="" type="radio"/> 65137	Medium 3D Printed (1-3/4") w/Power Plates (15 Degrees only)	\$0
--	--	-----



CRANK OPTIONS: (See Diagram Below)

CRANK STYLE:


<input type="radio"/> 65055	V Cranks	\$0
<input type="radio"/> 65056	S Cranks	\$0

CRANK LENGTH:


<input checked="" type="radio"/> 65040	180mm Long	\$0
--	------------	-----

CRANK WIDTH:

<input type="radio"/> 65021	20" Wide	\$0
<input type="radio"/> 65025	22" Wide	\$0



Hand Pedal
Size



Top End®

Force G™ Handcycle

CHAINRINGS:

<input type="radio"/> 65118	Single, 36T <i>Required with Rival 1x shifting package</i>	\$0
<input type="radio"/> 65059	Double, 48-34T <i>Required with Rival 2x shifting package</i>	\$0

SHIFTER OPTIONS: (see diagram)

<input type="radio"/> 65095	Shimano 105 2x Trigger Shifting Package <i>Shifter & brake mounted on right pedal only</i> Shimano 105 Shifter Shimano 105 Derailleurs Avid Brake Lever 11-34T 11 Speed Cassette Cable Driven chainring derailleur	\$0
<input type="radio"/> 65088	SRAM eTap 1x Eagle GX Electronic Shifting Package <i>Shifter & brake mounted on right pedal</i> SRAM eTap Shifter SRAM eTap Eagle GX 12 Speed Derailleur SRAM Time Trail Brake Lever 11-50T 12 speed cassette	\$1,415
<input type="radio"/> 65094	SRAM GX 1x Trigger Shifting Package <i>Shifter & brake mounted on right pedal only</i> SRAM GX Shifter SRAM GX Cable 12 Speed Derailleur Avid Brake Lever 11-50T 12 speed Cassette	\$0

BRAKING: (see definitions) (see exclusions)

<input type="radio"/> 65084	Rim Brakes, Tektro	\$0
<input type="radio"/> 65082	Front Disk Brake	\$361

CAMBER ANGLE:

<input checked="" type="radio"/> 45004	9 Degree Camber	\$0
--	-----------------	-----

WHEELBASE: (Front to Rear Axle) (See Definitions)

<input checked="" type="radio"/> 75019	Extended	\$0
--	----------	-----

WHEEL SETS:

(See Definitions: Includes Front and Rear Wheels, Tires and Tubes)

WHEEL SET SIZE:

<input checked="" type="radio"/> 50170	26" Front and Rear Wheels	\$0
--	---------------------------	-----

WHEEL SET TYPE:

<input checked="" type="radio"/> 50173	High Performance, 26" Front High Performance, 26" Rear	\$0
--	---	-----

ACCESSORIES:

STANDARD:

<input checked="" type="radio"/> 25067	Drafting Bumper	\$0
<input checked="" type="radio"/> 25072	Chain Guard <i>(half-moon, composite)</i>	\$0
<input checked="" type="radio"/> 25061	Safety Flag	\$0
<input checked="" type="radio"/> 25099	Mirror, Footrest Mounted	\$0

UPGRADABLE:

POSITIONING:

<input type="radio"/> 60016	Leg Support, Amputee	\$90
<input type="radio"/> 25039	Leg Guard (Tire/Chain)	\$90

HYDRATION:

<input checked="" type="radio"/> 25021	Water Bottle and Cage	\$60
--	-----------------------	------

Speed Chart for Shifting Packages

Shifting Package TE65094

	Cadence (Pedal RPM)			
Gear Combination	70	80	90	100
34x11	15.6	17.9	20.1	22.3
34x34	5.1	5.8	6.5	7.2
48x11	22.1	25.2	28.4	31.5
48x34	7.1	8.2	9.2	10.2

Shifting Package TE65094 and TE65088

	Cadence (Pedal RPM)			
Gear Combination	70.0	80.0	90.0	100.0
36x11	16.5	18.9	21.3	23.6
36x50	3.6	4.2	4.7	5.2

Top End®

Force G™ Handcycle

UPHOLSTERY:			
MATERIAL:			
<input checked="" type="radio"/>	30001	Nylon, Black Upholstery	\$0
BACK UPHOLSTERY:			
<input checked="" type="radio"/>	25064	Back Cushion	\$0
SEAT UPHOLSTERY:			
<input checked="" type="radio"/>	25063	Seat Cushion	\$0
<input checked="" type="radio"/>	25060	Hook and Loop Positioning Strap	\$0

FORK COLOR:			
<input type="radio"/>	55067	N/A	Same as Frame Color Not applicable with buffed frame, must select buffed fork
<input type="radio"/>	55041	Black	Wet Black
<input type="radio"/>	55042	Black	Black Prism
<input type="radio"/>	55040	Black	Stealth Black
<input type="radio"/>	55033	Blue	Midnight Blue
<input type="radio"/>	55034	Blue	Royal Blue
<input type="radio"/>	55035	Blue	Electric Blueberry
<input type="radio"/>	55046	Gray	Stealth Grey
<input type="radio"/>	55061	Green	Agriculture Green
<input type="radio"/>	55031	Green	Grasshopper Green
<input type="radio"/>	55062	Green	Neon Green
<input type="radio"/>	55029	Orange	Sunkissed Orange
<input type="radio"/>	55028	Orange	Simple Orange
<input type="radio"/>	55059	Orange	Neon Orange
<input type="radio"/>	55038	Pink	Bubblegum Pink
<input type="radio"/>	55065	Pink	Neon Pink
<input type="radio"/>	55063	Purple	Lilac Purple
<input type="radio"/>	55024	Red	Simple Red
<input type="radio"/>	55025	Red	Tangerine Red
<input type="radio"/>	55026	Red	Raspberry Red
<input type="radio"/>	55058	Red	Neon Red
<input type="radio"/>	55043	Silver	Silver Metallic
<input type="radio"/>	55039	White	White Snow
<input type="radio"/>	55030	Yellow	Safety Yellow
<input type="radio"/>	55060	Yellow	Neon Yellow
<input type="radio"/>	55106	Orange	Bomber Orange
<input type="radio"/>	55101	Blue	Carolina Blue
<input type="radio"/>	55099	Grey	Fighter Jet Grey
<input type="radio"/>	55068	Black	Flat Black
<input type="radio"/>	55096	Bronze	Last Vegas Bronze
<input type="radio"/>	55098	Pink	Light Pink
<input type="radio"/>	55097	Green	Mint Green

FRAME COLOR:			
<input type="radio"/>	55018	Black	Wet Black
<input type="radio"/>	55019	Black	Black Prism
<input type="radio"/>	55017	Black	Stealth Black
<input type="radio"/>	55010	Blue	Midnight Blue
<input type="radio"/>	55011	Blue	Royal Blue
<input type="radio"/>	55012	Blue	Electric Blueberry
<input type="radio"/>	55023	Gray	Stealth Grey
<input type="radio"/>	55053	Green	Agriculture Green
<input type="radio"/>	55008	Green	Grasshopper Green
<input type="radio"/>	55054	Green	Neon Green
<input type="radio"/>	55006	Orange	Sunkissed Orange
<input type="radio"/>	55005	Orange	Simple Orange
<input type="radio"/>	55051	Orange	Neon Orange
<input type="radio"/>	55015	Pink	Bubblegum Pink
<input type="radio"/>	55057	Pink	Neon Pink
<input type="radio"/>	55055	Purple	Lilac Purple
<input type="radio"/>	55001	Red	Simple Red
<input type="radio"/>	55002	Red	Tangerine Red
<input type="radio"/>	55003	Red	Raspberry Red
<input type="radio"/>	55050	Red	Neon Red
<input type="radio"/>	55020	Silver	Silver Metallic
<input type="radio"/>	55016	White	White Snow
<input type="radio"/>	55007	Yellow	Safety Yellow
<input type="radio"/>	55052	Yellow	Neon Yellow
<input type="radio"/>	55066	Raw	Buffed Aluminum

Note: All Special requests are subject to approval, additional charges, and additional lead time

Special Request 1)	
Special Request 2)	
Special Request 3)	
Special Request 4)	
Special Request 5)	

Top End®

Force G™ Handcycle

FORK COLOR

<input type="radio"/> 55100	Blue	Tiffany Blue	\$120
<input type="radio"/> 55104	Black	Wrinkle Black	\$0
FRAME COLOR			
<input type="radio"/> 55105	Orange	Bomber Orange	\$0
<input type="radio"/> 55095	Blue	Carolina Blue	\$0
<input type="radio"/> 55093	Grey	Fighter Jet Grey	\$0
<input type="radio"/> 55090	Bronze	Las Vegas Bronze	\$0
<input type="radio"/> 55092	Pink	Light Pink	\$0
<input type="radio"/> 55091	Green	Mint Green	\$0
<input type="radio"/> 55094	Blue	Tiffany Blue	\$0
<input type="radio"/> 55102	Black	Flat Black	\$0
<input type="radio"/> 55056	Purple	Purple Death	\$0
<input type="radio"/> 55103	Black	Wrinkle Black	\$0
<input type="radio"/> 25102		Water Bottle and Cage	\$0



Top End®

Force G™ Handcycle

HANDCYCLE SPECIFICATIONS:	
Weight:	Approximately 35 pounds
Turning Radius:	18 feet for 360 degree circle, 12 feet on a 3-point turn
Transfer Height:	8" from ground with 2" cushion
Frame Style:	Ovalized 6061 T6 aluminum aircraft tubing with internal skeletal structure for stiffness
Fork:	Removable, travel ready, with adjustable fork height and length
Propulsion / Shifting:	Up to 22 speeds with choice of shifting components.
Cranks:	Choice of style, length and width.
Back:	Adjustable angle from 55 to 90 degrees
Footrest:	Dual adjustable length, includes straps
Hand Pedals:	Ovalized ergonomic with frictionless bearings
Brake:	Caliper brakes standard, optional front disk brake
Overall Length:	86" overall length
Camber:	9 degree
Wheels:	26" wheels with optional upgrades
Axles:	Threaded
Tires / Tubes:	Standard tires included.
Positioning Straps:	Hook and loop positioning strap included
Cushions:	Seat and back cushion included
Safety / Visibility:	Flag, mirror and reflectors included
Chain Guard:	Black plastic half-moon shaped chain guard included.

DEFINITIONS:	
Assembled vs Unassembled and shipping:	If TE10048 (Assembled) is selected the bike will be mostly assembled upon delivery. Shipping rate is \$325 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
V Cranks:	V cranks position the hand pedals at a more narrow width than S cranks, but provide less clearance with the user's torso and legs.
S Cranks:	S Cranks are recommended for positioning the hand pedals at shoulder width or larger. S Cranks provide the most clearance with the user's torso and legs.
Wheelbase:	Wheelbase is measured from the middle of the rear axle to the middle of the front hub. - Extended length = 60"
Omit Wheels:	When an omit wheels option is selected, the following will NOT be included. Front and rear wheels, tires tubes, and axles.
Tektro Braking:	Braking packages include a caliper and cables for the hand brake lever included in the shifting package and flip lever parking brake assembly.
Disc Braking:	Braking package includes a disc brake rotor, caliper, cables and a Tektro flip lever parking brake assembly.
Omit Shifting:	When selected, the shifter, chainring derailleur, cassette derailleur, cassette, friction shifter for brake caliper, chain and shifter cables will NOT be included. The hand pedals and fork will be made to accept trigger shifter components mounted on the right hand pedal. A hand brake lever, caliper and cables as well as a flip-lever parking brake will be included. The brand/style is determined by the braking option selected.

Warranty Info:

Top End Sports LLC warrants the frames and forks when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Top End warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty