



USB
Charger
Port



LC-525 shown in
Durasoft Timber fabric

MADE IN
USA

INFINITY COLLECTION

The Infinity Collection from Pride® Power Lift Recliners features infinite-position sleep recline lift chairs designed to deliver comfort, style and performance. Available in many different fabric selections, the Infinity Collection models deliver maximum relaxation and will match any décor.

AMERICA'S #1
POWER LIFT RECLINER

POWER
RECLINERS

Protect your
investment with
the Pride® Lift
Chair Care Kit



1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)
pridemobility.com prideliftchairs.com

Pride®
Mobility Products Corp.
— LIVE YOUR BEST —

INFINITY COLLECTION

AMERICA'S #1
POWER LIFT RECLINER

POWER
LIFT
RECLINERS



LC-525 shown in
Durasoft Ember fabric



SPECIFICATIONS

	Petite Wide	Medium	Large
	LC-525PW ^{2,4}	LC-525M ^{2,4}	LC-525L ^{2,4}
	Infinite-Position Zero gravity Trendelenburg Chaise Lounger 5'2" & Below	Infinite-Position Zero gravity Trendelenburg Chaise Lounger 5'4"- 6'0"	Infinite-Position Zero gravity Trendelenburg Chaise Lounger 5'9"- 6'2"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	19"	20"	21"
Seat Depth	19"	20"	22"
Seat Width	24"	20"	22"
Top of Back to Seat	25.5"	26.5"	27.5"
Back Style	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back
Fully Padded Chaise	Yes	Yes	Yes
Heat & Massage Option	Yes	Yes	Yes
Total Comfort Seating	No	No	No
Distance From Wall	28"	29"	30"
Overall Chair Width	38"	33"	35"
Chair Weight	130 lbs.	130 lbs.	128 lbs.
Warranty	7-Year Prorated	7-Year Prorated	7-Year Prorated

INFINITY FABRIC SELECTION

QS

DURASOFT




LEXIS STA-KLEEN®

sta•kleen.
PERFORMANCE FABRIC




ULTRALEATHER™

Uf
ULTRALEATHER



ARIA³ CRYPTON[®]

super fabric



Look for the Quick Ship
Indicator for the shortest
lead times available.

BENEFITS

- Industry-first hand control with USB charger port
(not available with Heat and Massage)
- Standard head & arm covers
- Standard large dual pockets
- Removable back
- Wooden frame composed of engineered furniture grade laminate/hardwood for added strength and durability
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Fire retardant (California specifications; bulletin 117-2013)
- Integrated battery backup in case of power failure¹
- Transformer features self-diagnostic electronics with 9' cord
- Lifetime warranty on steel lift frame
- 7-year prorated warranty

OPTIONS*

- Deluxe Heat & Massage
- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Non-skid leg levelers
- Footrest extension available

Additional lead time may
apply for some options.

*Additional charges apply for all options.

1. Batteries not included.
2. Available with footrest extension.
3. Crypton® Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.
4. Available with Deluxe Heat & Massage only. Standard Heat & Massage not available.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.
The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)
pridemobility.com prideliftchairs.com

Pride
Mobility Products Corp.
— LIVE YOUR BEST —