




Protect your investment with the Pride® Power Lift Recliner Care Kit



LC-835 shown in Durasoft Deep Sky fabric

MADE IN
USA

Specialty COLLECTION

The Specialty Collection from Pride® Power Lift Recliners features the LC-835, a 2-position chaise lounger that offers the comfort and convenience of a lift chair in a space-saving design. The LC-835 has a stylish look and can be placed just five inches from the wall, making it the ultimate lift chair for smaller living spaces.

AMERICA'S **#1**
POWER LIFT RECLINER

POWER
RECLINERS

1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)
pridemobility.com prideliftchairs.com

Pride
Mobility Products Corp.
— LIVE YOUR BEST —

Specialty COLLECTION

AMERICA'S #1
POWER LIFT RECLINER

POWER
LIFT
RECLINERS



LC-835 shown in
Durasoft Deep Sky fabric

LC-835 FABRIC SELECTION



DURASOFT



EMBER



OAT



DEEP SKY



TIMBER



Look for the Quick Ship
Indicator for the shortest
lead times available.

OPTIONS*

- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Non-skid leg levelers

**Additional lead time may
apply for some options.**

*Additional charges apply for all options.
1. Batteries not included.

LC-835 WALL HUGGER



2-POSITION TV RECLINE

The 2-position TV recline is achievable on the LC-835. This is a comfortable position for performing many relaxing activities including reading, watching TV, working on your tablet or for conversation with family and friends.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.
The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



Medium

LC-835

2-Position
Wall Hugger

SPECIFICATIONS

Weight Capacity	300 lbs.
Seat-to-Floor	21"
Seat Depth	21"
Seat Width	20"
Top of Back to Seat	25"
Back Style	Blown Fiber Pillow Back
Fully Padded Chaise	Yes
Heat & Massage Option	No
Total Comfort Seating	Yes
Distance From Wall	5"
Overall Chair Width	37"
Chair Weight	155 lbs.
Warranty	2-Year Limited

BENEFITS

- Industry-first Total Comfort Seating features a unique spring design and one-piece soft foam seat for outstanding relaxation
- Wooden frame composed of engineered furniture grade laminate/hardwood for added strength and durability
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Standard head & arm covers
- Fire retardant (California specifications; bulletin 117-2013)
- Standard, large dual pockets
- Removable back
- 300 lbs. weight capacity
- Integrated battery backup in case of power failure¹
- Transformer features self-diagnostic electronics with 9' cord
- Lifetime warranty on steel lift frame
- 2-year limited warranty

1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)
pridemobility.com prideliftchairs.com

Pride
Mobility Products Corp.
— LIVE YOUR BEST —