



# Zippie Voyage

May 2013

EIZ14

|  |  |
|--|--|
| Date: _____<br>Dealer Acct #: _____<br>Dealer: _____<br>Dealer Contact: _____<br>Dealer Address: _____<br>Dealer City: _____ ST: _____ ZIP: _____<br>Dealer Phone: ( ) _____ Fax: ( ) _____<br><br>Confirm Via: <input type="checkbox"/> Email <input type="checkbox"/> Fax<br><br>Mark For: _____ | <input type="checkbox"/> Quote <input type="checkbox"/> Order<br><br><p style="text-align: center;"><b>SHIPPING INFORMATION</b></p> PO#: _____<br>Ship To: _____<br>Attention: _____<br>Address: _____<br>Address: _____<br>Ship To City: _____ ST: _____ ZIP: _____<br>Ship To Phone: ( ) _____ Fax: ( ) _____<br>Confirmation Email: _____<br><br>DOB (Date of Birth): mn/dd/yr ____/____/____ |
|--|--|

Designed in conjunction with Baby Jogger allowing their patented quick fold technology which allows one to fold a base unit after removing the seat in one simple step and with a secure fold latch for easy storage. Each unit includes:

- Hand operated Wheel Lock
- Quick release Front and rear wheels
- 45 Tilt
- 12" rear tires & Lt. weight front wheels
- Telescoping stroller handle
- Fwd & Rwd facing seating capability
- Front caster suspension & swivel locks
- Adjustable seating Shell
- Under carriage storage
- Canopy w/ breathable window/mesh pockets
- (may be limited when some options chosen)

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.  
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

| ✓ | Option # | HCPCS | Definition | Price | ✓ | Option # | HCPCS | Definition | Price |
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|

## Voyage (Weight Capacity 60#)

### 1. FRAME STYLE

|                          |        |       |  |            |                          |        |       |  |            |
|--------------------------|--------|-------|--|------------|--------------------------|--------|-------|--|------------|
| <input type="checkbox"/> | 198BS7 | E1232 | Voyage with Tilt<br>Includes Solid Seat and Back | \$2,375.00 | <input type="checkbox"/> | 198BS8 | E1234 | Voyage with Tilt<br>Seat and Back not included | \$2,135.00 |
|--------------------------|--------|-------|--|------------|--------------------------|--------|-------|--|------------|

### 2. BASE OPTIONS (Vent tray options in accessory section)

|                          |        |       |                                   |          |                          |         |  |  |          |
|--------------------------|--------|-------|-----------------------------------|----------|--------------------------|---------|--|--|----------|
| <input type="checkbox"/> | 198KJ9 | E1014 | Recline Feature ( 90 - 135 range) | \$325.00 | <input type="checkbox"/> | 198KJ10 |  | WC-19 Transit Option<br>(includes H Harness) | \$315.00 |
|--------------------------|--------|-------|-----------------------------------|----------|--------------------------|---------|--|--|----------|

### 3. COLORS

#### Color Options - Solid

| Solid Colors             |         |             |    |
|--------------------------|---------|-------------|----|
| <input type="checkbox"/> | 198CL9  | Purple      | NC |
| <input type="checkbox"/> | 198CL20 | Pink Sorbet | NC |
| <input type="checkbox"/> | 198CL11 | Green       | NC |
| <input type="checkbox"/> | 198CL12 | Navy Blue   | NC |
| <input type="checkbox"/> | 198CL13 | Orange      | NC |
| <input type="checkbox"/> | 198CL17 | Turquoise   | NC |
| <input type="checkbox"/> | 198CL16 | Black       | NC |
| <input type="checkbox"/> | 198CL18 | Burgundy    | NC |

#### Custom Option (Only available with Prints)

|                          |         |   |         |
|--------------------------|---------|---|---------|
| <input type="checkbox"/> | 198C001 | Custom Color Selection: Choose &<br>Enter from Solid Color Options: | \$70.00 |
|--------------------------|---------|---|---------|

#### Canopy Pattern

|                          |         |               |    |
|--------------------------|---------|---------------|----|
| <input type="checkbox"/> | 198UF01 | Cross Over    | NC |
| <input type="checkbox"/> | 198UF02 | Center Stripe | NC |

#### Color Options - Print (Std. Trim listed behind each print)

|                          |         |                                  |    |
|--------------------------|---------|----------------------------------|----|
| <input type="checkbox"/> | 198CL19 | Leopard (Std w/ Black Trim)      | NC |
| <input type="checkbox"/> | 198CL14 | Butterflies (Std w/Purple Trim)  | NC |
| <input type="checkbox"/> | 198CL15 | Red Race Flags (Std w/ Red Trim) | NC |



The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.  
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

| ✓ | Option # | HCPCS | Definition | Price | ✓ | Option # | HCPCS | Definition | Price |
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|

**4. SEAT OPTIONS (Depth adjustable from 5" - 13" with 1" backrest)**

Choose section A (when ordering base with seat/back) or B. (when ordering base without seat/back)

**Section A (seat options when base with seat/back is ordered)**

**Cushion includes:**  
 - Contoured Medium Density Foam  
 - 6" to 13" Hip Width  
 - 3" Lateral Hip Blocks (Pr)  
 (Hip Blocks Not avail. on Planar w/o Hip Block cushion)  
 - Slight Medial Thigh Contour  
 - Fixed Ischial Ledge  
 - Built in Depth Growth  
 - Removable Neoprene Cover

|                          |         |       |                                    |          |
|--------------------------|---------|-------|------------------------------------|----------|
| <input type="checkbox"/> | 198CU28 | E2292 | Planar Seat w/o Hip Blocks         | NC       |
| <input type="checkbox"/> | 198CU41 | E2292 | Planar Seat w/ Hip Block (Pr)      | \$225.00 |
| <input type="checkbox"/> | 198CU42 | E2294 | Contoured Wedge Seat               | \$225.00 |
| <input type="checkbox"/> | 198CU43 | E2294 | Contoured Anti Thrust (Ledge) Seat | \$225.00 |

**Choose Depth (Based on 1" T Back)**

|                          |         |  |                             |    |
|--------------------------|---------|--|-----------------------------|----|
| <input type="checkbox"/> | 198SD78 |  | Short 5" - 8" Depth         | NC |
| <input type="checkbox"/> | 198SD79 |  | Medium 7" through 10" Depth | NC |
| <input type="checkbox"/> | 198SD80 |  | Long 9" through 13" Depth   | NC |

**ADJUSTABLE CONTOUR SEAT includes:**  
 -Adjustable Pre-Ischial Ledge  
 - Adjustable 1" High Foam Wedges for:  
 Lateral Hip and Thigh Support  
 Medial Knee and Thigh Support  
 - 3" Seat Depth Growth Kit  
 - Removable Neoprene Cover

|                          |         |       |                         |          |
|--------------------------|---------|-------|-------------------------|----------|
| <input type="checkbox"/> | 198ST65 | E2294 | Adjustable Contour Seat | \$320.00 |
|--------------------------|---------|-------|-------------------------|----------|

**Choose Foam**

|                          |        |  |                  |    |
|--------------------------|--------|--|------------------|----|
| <input type="checkbox"/> | 198FM1 |  | Soft Foam        | NC |
| <input type="checkbox"/> | 198FM2 |  | Medium/Soft Foam | NC |
| <input type="checkbox"/> | 198FM3 |  | Medium Foam      | NC |

**Choose Depth (Based on 1" T Back)**

|                          |         |  |                                 |    |
|--------------------------|---------|--|---------------------------------|----|
| <input type="checkbox"/> | 198SD78 |  | Short 5" - 8"                   | NC |
| <input type="checkbox"/> | 198SD79 |  | Medium 7" through 10"           | NC |
| <input type="checkbox"/> | 198SD80 |  | Long 9" through 13" (4" growth) | NC |

OR

**Section B. (seat options when base ordered without standard seat/back)**

**Cushion includes:**  
 - Contoured Medium Density Foam  
 - 6" to 13" Hip Width  
 - 3" Lateral Hip Blocks (Pr)  
 (Hip Blocks Not avail. on Planar w/o Hip Block cushion)  
 - Slight Medial Thigh Contour  
 - Fixed Ischial Ledge  
 - Built in Depth Growth  
 - Removable Neoprene Cover

|                          |         |       |                                    |          |
|--------------------------|---------|-------|------------------------------------|----------|
| <input type="checkbox"/> | 198CU39 | E2292 | Planar Seat w/o Hip Blocks         | \$225.00 |
| <input type="checkbox"/> | 198CU44 | E2292 | Planar Seat w/ Hip Blocks(Pr)      | \$335.00 |
| <input type="checkbox"/> | 198CU37 | E2294 | Contoured Wedge Seat               | \$335.00 |
| <input type="checkbox"/> | 198CU38 | E2294 | Contoured Anti Thrust (Ledge) Seat | \$335.00 |

**Choose Depth (Based on 1" T Back)**

|                          |         |  |                             |    |
|--------------------------|---------|--|-----------------------------|----|
| <input type="checkbox"/> | 198SD78 |  | Short 5" - 8" Depth         | NC |
| <input type="checkbox"/> | 198SD79 |  | Medium 7" through 10" Depth | NC |
| <input type="checkbox"/> | 198SD80 |  | Long 9" through 13" Depth   | NC |

**ADJUSTABLE CONTOUR SEAT includes:**  
 -Adjustable Pre-Ischial Ledge  
 - Adjustable 1" High Foam Wedges for:  
 Lateral Hip and Thigh Support  
 Medial Knee and Thigh Support  
 - 3" Seat Depth Growth Kit  
 - Removable Neoprene Cover

|                          |         |       |                         |          |
|--------------------------|---------|-------|-------------------------|----------|
| <input type="checkbox"/> | 198ST60 | E2294 | Adjustable Contour Seat | \$420.00 |
|--------------------------|---------|-------|-------------------------|----------|

**Choose Foam**

|                          |        |  |                  |    |
|--------------------------|--------|--|------------------|----|
| <input type="checkbox"/> | 198FM1 |  | Soft Foam        | NC |
| <input type="checkbox"/> | 198FM2 |  | Medium/Soft Foam | NC |
| <input type="checkbox"/> | 198FM3 |  | Medium Foam      | NC |

**Choose Depth (Based on 1" T Back)**

|                          |         |  |                       |    |
|--------------------------|---------|--|-----------------------|----|
| <input type="checkbox"/> | 198SD78 |  | Short 5" - 8"         | NC |
| <input type="checkbox"/> | 198SD79 |  | Medium 7" through 10" | NC |
| <input type="checkbox"/> | 198SD80 |  | Long 9" through 13"   | NC |

OR

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.  
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

| ✓ | Option # | HCPCS | Definition | Price | ✓ | Option # | HCPCS | Definition | Price |
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|

### 5. BACKRESTSHELL

**Desired Backrest Height Setting** (Bold indicates default, when not picked)

|                          |          |  |                        |    |                          |     |                          |     |                          |     |                          |            |                          |     |                          |     |                          |     |                          |     |
|--------------------------|----------|--|------------------------|----|--------------------------|-----|--------------------------|-----|--------------------------|-----|--------------------------|------------|--------------------------|-----|--------------------------|-----|--------------------------|-----|--------------------------|-----|
| <input type="checkbox"/> | 198BT111 |  | Short 7" to 13" Height | NC | <input type="checkbox"/> | 7"  | <input type="checkbox"/> | 8"  | <input type="checkbox"/> | 9"  | <input type="checkbox"/> | <b>10"</b> | <input type="checkbox"/> | 11" | <input type="checkbox"/> | 12" | <input type="checkbox"/> | 13" |                          |     |
| <input type="checkbox"/> | 198BT112 |  | Tall 10" to 17" Height | NC | <input type="checkbox"/> | 10" | <input type="checkbox"/> | 11" | <input type="checkbox"/> | 12" | <input type="checkbox"/> | <b>13"</b> | <input type="checkbox"/> | 14" | <input type="checkbox"/> | 15" | <input type="checkbox"/> | 16" | <input type="checkbox"/> | 17" |

### 6. BACK OPTIONS

Choose section A (when ordering base with seat/back) or B. (when ordering base without seat/back)

#### Section A (Backrest options when base with seat/back is ordered)

| " T " Back 1" Thick      |         |       |  | OR | " T " Back 2.5" Thick    |   |       |  |         |
|--------------------------|---------|-------|--|----|--------------------------|---|-------|--|---------|
| <input type="checkbox"/> | 198BT35 | E2291 | Short 8" to 10"                                    |    | NC                       | 2.5" "T" Back is used for seat depths as small as 3.5", backrest has removable foam for additional depth growth without adjusting seat pan. |       |  |         |
| <input type="checkbox"/> | 198BT28 | E2291 | Medium 10" to 12"                                  | NC | <input type="checkbox"/> | 198BT115  | E2291 | Short 8" to 10"                                    | \$70.00 |
| <input type="checkbox"/> | 198BT29 | E2291 | Tall 12" to 14"                                    | NC | <input type="checkbox"/> | 198BT116  | E2291 | Medium 10" to 12"                                  | \$70.00 |
| <input type="checkbox"/> | 198BT56 | E2291 | Extra Tall 14" to 16" (n/a w Short Backrest Shell) | NC | <input type="checkbox"/> | 198BT117  | E2291 | Tall 12" to 14"                                    | \$70.00 |
|                          |         |       |  |    | <input type="checkbox"/> | 198BT118  | E2291 | Extra Tall 14" to 16" (n/a w Short Backrest Shell) | \$70.00 |

#### Section B. (Backrest options when base ordered without standard seat/back)

| " T " Back 1" Thick      |          |       |  | OR       | " T " Back 2.5" Thick    |   |       |  |          |
|--------------------------|----------|-------|--|----------|--------------------------|---|-------|--|----------|
| <input type="checkbox"/> | 198BT103 | E2291 | Short 8" to 10"                                    |          | \$305.00                 | 2.5" "T" Back is used for seat depths as small as 3.5", backrest has removable foam for additional depth growth without adjusting seat pan. |       |  |          |
| <input type="checkbox"/> | 198BT104 | E2291 | Medium 10" to 12"                                  | \$305.00 | <input type="checkbox"/> | 198BT107  | E2291 | Short 8" to 10"                                    | \$375.00 |
| <input type="checkbox"/> | 198BT105 | E2291 | Tall 12" to 14"                                    | \$305.00 | <input type="checkbox"/> | 198BT108  | E2291 | Medium 10" to 12"                                  | \$375.00 |
| <input type="checkbox"/> | 198BT106 | E2291 | Extra Tall 14" to 16" (n/a w Short Backrest Shell) | \$305.00 | <input type="checkbox"/> | 198BT109  | E2291 | Tall 12" to 14"                                    | \$375.00 |
|                          |          |       |  |          | <input type="checkbox"/> | 198BT110  | E2291 | Extra Tall 14" to 16" (n/a w Short Backrest Shell) | \$375.00 |

### 7. HEAD SUPPORT OPTIONS (Hardware allows the headrest to be flushmounted to backrest when above shell)

When choosing a headrest, Hardware Bracket must be chosen as well

|                          |         |       |                             |          |                          |         |       |   |          |
|--------------------------|---------|-------|-----------------------------|----------|--------------------------|---------|-------|---|----------|
| <input type="checkbox"/> | 198HE40 | E0955 | Whitmyer 8" Plush           | \$225.00 | <input type="checkbox"/> | 198HE29 | E0955 | Whitmyer Sm Adjust-a-Plush  | \$270.00 |
| <input type="checkbox"/> | 198HE41 | E0955 | Whitmyer 10" Plush          | \$225.00 | <input type="checkbox"/> | 198HQ6  | K0108 | Whitmyer DFS (Forehead Strap)<br>Squared for Single pad, <b>must order</b><br>Whitmyer headrest | \$120.00 |
| <input type="checkbox"/> | 198HE42 | E0955 | Whitmyer SM Contour Cradle  | \$225.00 |                          |         |       |   |          |
| <input type="checkbox"/> | 198HE43 | E0955 | Whitmyer Med Contour Cradle | \$225.00 | <input type="checkbox"/> | 198HE31 | E1028 | Hardware Bracket Only (can be chosen w/o headrest)  | \$105.00 |

### 8. FOOTREST

|                          |         |       |                                |          |
|--------------------------|---------|-------|--------------------------------|----------|
| <input type="checkbox"/> | 198F312 | K0040 | Angle Adj. Footrest (Required) | \$155.00 |
|--------------------------|---------|-------|--------------------------------|----------|

Footrest adjustment range 1.5" - 12", Short Range, the footplate will be inverted. Choose desired setting, if not chosen Long range is default.

|                          |         |       |                                      |    |                          |         |  |  |    |
|--------------------------|---------|-------|--------------------------------------|----|--------------------------|---------|--|--|----|
| <input type="checkbox"/> | 198F313 | K0037 | <b>Preset in short range 1.5"-4"</b> | NC | <input type="checkbox"/> | 198F314 |  | <b>Preset in long range 4.5" - 12"</b> | NC |
|--------------------------|---------|-------|--------------------------------------|----|--------------------------|---------|--|--|----|

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.  
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

| ✓ | Option # | HCPCS | Definition | Price | ✓ | Option # | HCPCS | Definition | Price |
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|

### 9. SUPPORTS

Laterals Height Adjustment Range based on Backrest shell height chosen in Sec. 5. i.e. Short Height Lats adj. 3" - 12", Medium Height Lats Adj. 3" - 17"

For the supports, Voyage will be set up at factory, insert dimensions on last page spec drawing. Not required, will default to std. settings

**\*Priced as each, selecting these options indicates Qty 2. Hardware Bracket must be chosen with Medial Support.**

#### Lateral Hip Support

Only available with Planar seat without Hip Blocks and Adj. Contour seat.

|                          |          |       |                                  |          |
|--------------------------|----------|-------|----------------------------------|----------|
| <input type="checkbox"/> | 198KC4/6 | E0956 | SM (4" x 4", 3/8" thick)*        | \$100.00 |
| <input type="checkbox"/> | 198KC3/5 | E0956 | LG (5.38" x 4 3/4", 1/2" thick)* | \$100.00 |

#### Lateral Thigh Support

Only available with Planar seat without Hip Blocks and Adj. Contour seat.

|                          |          |       |                                  |          |
|--------------------------|----------|-------|----------------------------------|----------|
| <input type="checkbox"/> | 198KF4/6 | E0956 | SM (4" x 4", 3/8" thick)*        | \$100.00 |
| <input type="checkbox"/> | 198KF3/5 | E0956 | LG (5.38" x 4 3/4", 1/2" thick)* | \$100.00 |

#### Medial Support

|                          |         |       |                              |          |
|--------------------------|---------|-------|------------------------------|----------|
| <input type="checkbox"/> | 198KH12 | E0957 | SM (3"H x 2.5"W x 4.5"L)     | \$125.00 |
| <input type="checkbox"/> | 198KH13 | E0957 | LG ( 3 3/4"H x 4"W x 6"L)    | \$125.00 |
| <input type="checkbox"/> | 198HE77 | E1028 | Hardware Medial Support Brkt | \$225.00 |

#### Lateral Trunk Support

|                          |            |       |   |          |
|--------------------------|------------|-------|---|----------|
| <input type="checkbox"/> | 198KE25/24 |       | Fixed Angle Adjustable*                             | NC       |
|                          | 198KD15/19 | E0956 | <input type="checkbox"/> Large Straight Pad 3" x 6" | \$100.00 |
|                          | 198KD13/17 | E0956 | <input type="checkbox"/> Small Straight Pad 2" x 6" | \$100.00 |
| <input type="checkbox"/> | 198KE42/41 | E1028 | Swing Away*   | \$225.00 |
|                          | 198KD15/19 | E0956 | <input type="checkbox"/> Large Straight Pad 3" x 6" | \$100.00 |
|                          | 198KD16/20 | E0956 | <input type="checkbox"/> Large Curved Pad 3" x 6"   | \$100.00 |
|                          | 198KD13/17 | E0956 | <input type="checkbox"/> Small Straight Pad 2" x 6" | \$100.00 |
|                          | 198KD14/18 | E0956 | <input type="checkbox"/> Small Curved Pad 2" x 6"   | \$100.00 |

### 10. POSITIONING AND TRUNK SUPPORTS

#### Pelvic Positioning Belt

Must use when ordering Anterior Trunk Supports.

|                          |         |       |                            |         |
|--------------------------|---------|-------|----------------------------|---------|
| <input type="checkbox"/> | 198A14  | E0978 | 1" Padded                  | \$45.00 |
| <input type="checkbox"/> | 198A15  | E0978 | 1.5" Padded                | \$45.00 |
| <input type="checkbox"/> | 198A263 | E0978 | 1" Padded Auto Style       | \$55.00 |
| <input type="checkbox"/> | 198A256 | E0978 | 1.5" Padded Auto Style     | \$55.00 |
| <input type="checkbox"/> | 198A818 | E0978 | Comfort Fit APB5145 XS Pad | \$95.00 |
| <input type="checkbox"/> | 198A819 | E0978 | Comfort Fit APB5146 Sm Pad | \$95.00 |

#### Comfort Fit ATS Anterior Trunk Support (Lil Kiddos)

|                          |         |       |                              |          |
|--------------------------|---------|-------|------------------------------|----------|
| <input type="checkbox"/> | 198A820 | E0960 | Com Fit ATS9107 xxs 9.5 pad  | \$152.00 |
| <input type="checkbox"/> | 198A821 | E0960 | Com Fit ATS9100 xs 12.5" pad | \$152.00 |
| <input type="checkbox"/> | 198A822 | E0960 | Com Fit ATS9101 sm 15"pad    | \$152.00 |

#### Anterior Trunk Butterfly Support

Dimensions are seat to shoulder height.

|                          |        |       |                                   |          |
|--------------------------|--------|-------|-----------------------------------|----------|
| <input type="checkbox"/> | 198K11 | E0960 | XX Small Butterfly Pad (< 10")    | \$110.00 |
| <input type="checkbox"/> | 198K12 | E0960 | X Small Butterfly Pad (10" - 12") | \$110.00 |
| <input type="checkbox"/> | 198K13 | E0960 | Small Butterfly Pad (12" - 14")   | \$110.00 |
| <input type="checkbox"/> | 198K14 | E0960 | Medium Butterfly Pad (14" - 16")  | \$110.00 |
| <input type="checkbox"/> | 198K15 | E0960 | Large Butterfly Pad (>16")        | \$110.00 |

#### Trunk Support Options

|                          |        |       |                     |         |
|--------------------------|--------|-------|---------------------|---------|
| <input type="checkbox"/> | 198SS1 | E0960 | Shoulder Strap Pads | \$25.00 |
|--------------------------|--------|-------|---------------------|---------|

#### Anterior Trunk H Strap Support (included with Transit)

|                          |        |       |                          |          |
|--------------------------|--------|-------|--------------------------|----------|
| <input type="checkbox"/> | 198SS2 | E0960 | XX Small H Strap 5 Point | \$165.00 |
|--------------------------|--------|-------|--------------------------|----------|

### 11. WHEELS AND CASTERS

Standard with Quick release rear and front wheels, and front wheel suspension with caster swivel locks

|                          |         |  |  |     |
|--------------------------|---------|--|--|-----|
| <input type="checkbox"/> | 198WQ07 |  | 12" Pneumatic w/ Airless Inserts with 7" casters | STD |
|--------------------------|---------|--|--|-----|

Customer Service: 800-333-4000 Fax: 800-300-7502 www.sunrisemedical.com

Specifications are subject to change without notice. Copyright © 2013 Sunrise Medical 120448 REV. B

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.  
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

| ✓ | Option # | HCPCS | Definition | Price | ✓ | Option # | HCPCS | Definition | Price |
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|

**12. ACCESSORIES (See Available Pictures Below)**

| Shoe Holders             |        |                |                      |          |
|--------------------------|--------|----------------|----------------------|----------|
| <input type="checkbox"/> | 198A68 | E0952<br>E0951 | 6" Long Shoe Holders | \$120.00 |
| <input type="checkbox"/> | 198A69 | E0952<br>E0951 | 8" Long Shoe Holders | \$120.00 |

| Foot Straps              |        |       |                             |         |
|--------------------------|--------|-------|-----------------------------|---------|
| <input type="checkbox"/> | 198A79 | E0952 | 16" Foot Straps with Covers | \$35.00 |
| <input type="checkbox"/> | 198A80 | E0952 | 18" Foot straps with Covers | \$35.00 |

| Upper Extremity Support Tray (Angle Adjustable Hardware Included) |        |       |            |          |
|---|--------|-------|------------|----------|
| (Tray Mount adjusts from 4.5"-7.5" with 1" tall cushion )         |        |       |            |          |
| <input type="checkbox"/>  | 198A91 | E0950 | Light Grey | \$145.00 |
| <input type="checkbox"/>  | 198A92 | E0950 | Clear Poly | \$165.00 |

| Misc Accessories         |         |  |                                     |          |
|--------------------------|---------|--|-------------------------------------|----------|
| <input type="checkbox"/> | 198LF5  |  | Cup Holder/ Bottle/Feeding          | \$45.00  |
| <input type="checkbox"/> | 198LF7  |  | Utility Hook (i.e.suction canister) | \$25.00  |
| <input type="checkbox"/> | 198A462 |  | Rain Cape                           | \$175.00 |

| Floor Sitter/ Feeder Base |        |       |                                 |          |
|---------------------------|--------|-------|---------------------------------|----------|
| <input type="checkbox"/>  | 198A78 | E1399 | Folding Booster Base (Ht. Adj.) | \$625.00 |

| Accessories              |         |                |   |          |
|--------------------------|---------|----------------|---|----------|
| <input type="checkbox"/> | 198LF1  | E1029          | Dovetail Vertical Frame Mounted Hardware for LTV ventilators (limit 15#) includes 1 Utility hook  | \$340.00 |
| <input type="checkbox"/> | 198LF11 | E1029          | Dovetail Vertical Frame Mounted Hardware for Trilogy ventilators (limit 15#) includes 1 Utility hook                                    | \$340.00 |
| <input type="checkbox"/> | 198LF2  | E1029          | Battery and Vent/Utility Underframe Tray (limit 50#, ideal for U1 or Grp 22 Batteries) includes 1 Utility hook                          | \$400.00 |
| <input type="checkbox"/> | 198LF3  | E1029          | Front Mount Utility platform (limit 50#)  | \$355.00 |
| <input type="checkbox"/> | 198LF4  | E1029          | Rear Add on Accessory platform (limit 45# or one Grp 24 battery) includes 1 Utility hook. Must be removed or folded before folding base | \$135.00 |
| <input type="checkbox"/> | 198A74  | K0105          | IV Support Pole   | \$145.00 |
| <input type="checkbox"/> | 198A713 | E2208          | O2 Tank Holder  | \$175.00 |
| <input type="checkbox"/> | 198A714 | K0105<br>E2208 | IV/O2 Combination   | \$245.00 |


|   |   |   |  |   |
|---|---|---|--|---|
|  |  |  |  |  |
| IV Support Pole   | Cup Holder  | Front Mount Utility Platform  | Utility Hook   | Upper Extremity Support Tray  |

|   |   |   |  |   |
|---|---|---|--|---|
|  |  |  |  |  |
| Folding Booster Base  | O2 Tank Holder  | Battery and Vent/Utility Underframe Tray  | Dovetail Vertical Vent Mount (LTV/Trilogy)   | Rear Add On Accessory Platform  |

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.  
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

| ✓ | Option # | HCPCS | Definition | Price | ✓ | Option # | HCPCS | Definition | Price |
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|

**Desired Set Up Specifications (Not Required)**



Child's Name \_\_\_\_\_  
 Age: \_\_\_\_\_  
 Diagnosis: \_\_\_\_\_  
 Height: \_\_\_\_\_  
 Weight: \_\_\_\_\_

- |                                    |                                    |
|------------------------------------|------------------------------------|
| 1. Top of Head to Seat: _____      | 7. Width of Head: _____            |
| 2. Top of Shoulder to Seat*: _____ | 8. Shoulder Width: _____           |
| 3. Under Arm to Seat: _____        | 9. Chest Width: _____              |
| 4. Trunk Depth: _____              | 10. Hip Width: _____               |
| 5. Actual Sitting Depth*: _____    | 11. Left Lower Leg Length*: _____  |
| 6. Foot Length: _____              | 12. Right Lower Leg Length*: _____ |

\* Critical Dimensions