

Seatar

A chair for indoor use
Weight 11.6 kg
Self-propelled.

The chair is manufactured
using injection technology
with advanced plastic
substances.

All the metal parts are made
from stainless steel
(including the bearings).

Hermetic brake locking.

Folding push up handgrips.

Adjustable footrest.

An option for a chamber
pot.

An option for a rigid
suitcase.



This chair enables
safe passage from
chair to chair and is
convenient for use.



Folding and
assembling can be
performed by the
user while sitting as
well.

The small folded
dimensions enable
the chair's availability
for any purpose.



WheelAble

Be Spontaneous



The **Zriha Group**, 6 Shrir St, Koren Ind. Zone, Maalot Tel: +972-3-906-6906 info@seataro.com

Seataro

WheelAble

Be Spontaneous



INNOVATION TAKES YOU TO NEW PLACES

www.seataro.com

A FOLDING CHAIR FOR INACCESSIBLE TOILETS AND SHOWERS