



| Item # | Seat                | Color                                                                                             | UOM  |
|--------|---------------------|---------------------------------------------------------------------------------------------------|------|
| FLEX   | 16.75" Stadium Seat | Electric Blue  | 1/ea |



## ZooMe Flex Folding Travel Scooter

### Travel Scooter

- Ultra-lightweight for superior portability.
- Quick and Easy "Flex" folding design for transport and storage in one easy motion.
- High-gloss shrouds for a protective finish.
- Easy throttle control.
- Delta tiller for users with limited dexterity.
- Auto-folding seat with padded, adjustable armrests and fold down back.
- Flat-free, non-marking tires for worry-free travel.
- Easy-to-adjust tiller angle places the scooter's controls in the ideal position.
- Anti-tip wheels for additional safety.

### specifications

|                    |                             |                       |                                                                            |
|--------------------|-----------------------------|-----------------------|----------------------------------------------------------------------------|
| Weight Capacity    | 275 lbs.                    | Brakes                | Electromagnetic                                                            |
| Top Speed          | 4 mph                       | Freewheel Mode        | No                                                                         |
| Estimated Range    | 13 miles                    | Front Wheels          | 6" x 1.5"                                                                  |
| Turning Radius     | 53"                         | Rear Wheels           | 7" x 2.35"                                                                 |
| Climbing Angle     | 6°                          | Anti-Tip Wheels       | 2" x .5"                                                                   |
| Ground Clearance   | 2"                          | Base Weight           | 53.6 lbs.                                                                  |
| Overall Dimensions | 37" (L) x 18" (W) x 36" (H) | Battery Weight (Pair) | 5 lbs.                                                                     |
| Seat Dimensions    | 16.75" (W) x 15" (D)        | Heaviest Piece Weight | 53.6 lbs                                                                   |
| Motor              | 24V x 120W                  | Total Weight          | 58.6 lbs.                                                                  |
| Controller         | PG S Drive 45A              | Warranty              | 14 months on Electronics, 6 months on Batteries, Limited Lifetime on Frame |
| Batteries          | 25.9V x 10AH Lithium        |                       |                                                                            |
| Battery Charger    | 2A Offboard                 |                       |                                                                            |