

Soft-Touch® Booster Car Seat™



Soft contoured support for a comfortable ride


- Portable, lightweight and durable
- FMVSS213 compliant
- FAA compliant
- Vehicle LATCH system required
- Easy to wipe clean
- Attachment Straps for a variety of uses
- Small Weight Capacity: 51 lbs. to 90 lbs.
- Large Weight Capacity: 81 lbs. to 130 lbs.

Approved for:  





Gianna in the Small:
Height (45")
Weight (51 lbs)
Age (6 y)

The Special Tomato Booster Car Seat... taking you where you want to go!




Color Options and Features



Soft-Touch® Booster Car Seat	Product Number
Dark Gray Small	71010309
Lilac Small	71010302
Aqua Small	71010304
Dark Gray Large	71010409
Lilac Large	71010402
Aqua Large	71010404

Seat Specifications	Small	Large
Inside Seat Depth	11.5" (29 cm)	13.5" (34 cm)
Inside Seat Width	12" (30.5 cm)	14" (35.5 cm)
Seat Back Height	29" (74 cm)	35" (89 cm)
Harness Shoulder Height	14.5-18.5" (37- 47 cm)	17.75-22.25" (45 – 56.5 cm)
Overall Seat Depth	16" (41 cm)	19" (48 cm)
Overall Seat Height*	34" (86 cm)	40" (102 cm)
Overall Seat Width	16" (41 cm)	18" (46 cm)
Weight Capacity (Min)	51 lbs (23 kg)	81 lbs (37 kg)
Weight Capacity (Max)	90 lbs (41 kg)	130 lbs (59 kg)
Weight of Seat	13 lbs (6 kg)	15 lbs (7 kg)



*Overall seat height: Customer **MUST** measure from the seat surface to the headliner to be sure the large will fit into your vehicle.

Optional
Seat Base
Choices!

