



SPECIFICATIONS

	EC Transport	EC 2000	EC 2000HD	EC 3000	EC 4000	EC 4000 Recliner
Seat Width	17"/19"	16"/18"/20"	22"/24"	16"/18"/20"	16"/18"/20"	16"/18"/20"
Seat Depth	16"	16"	18"/18"	16"	16"/16"/18"	17"
Overall Width						
Fixed Arm	20.5"/22.5"	22.5"/24.5"	N/A	N/A	N/A	N/A
Removable Arm	N/A	24.5"/26.5"/29.5"	30.5"/32.5"	23"/25"	23"/25"/27"	25"
Overall Depth						
Fixed Arm, Fixed F/R	N/A	41.5"	N/A	N/A	N/A	N/A
Fixed Arm w/o riggings	N/A	32"	N/A	N/A	N/A	N/A
Removable w/o riggings	24"	29"	29"/32"/32"	30.5"	30.5"	34.5"
Seat to Floor Height	18.5"	19.75"	19.75"	19.75"	19.75"	19.75"
Back Style	Fold Down	Fixed	Fixed	Fixed	Adjustable	Recline
Back Height	18.5"	16"	16"	16"	16" - 18"	22"
Arm Styles	Fixed Full	Fixed Full/ Removable Full Removable Desk	Removable Full Removable Desk	Removable Desk	Flip-Back Desk Flip Back Full	Removable Desk (on S/A) Ht. Adj. Full (on ELR)
Front Riggings	Swing-away	Fixed Swing-away Elevating	Swing-away Elevating	Swing-away Elevating	Swing-away Elevating	Swing-away Elevating
Rear Axle	Single Position	Dual Position	Single Posit.	Dual Position	Dual Position	Dual Position
Rear Wheels	8x1 in. Rubber	Composite Mag	Composite Mag	Composite Mag	Composite Mag	Composite Mag
Hand Rims	N/A	Composite	Chromed Steel	Composite	Aluminum	Aluminum
Upholstery	Nylon	Vinyl	Vinyl	Nylon	Nylon	Nylon
Product Weight (w/o riggings)	18.4 lbs. (Aluminum) 24 lbs. (Steel)	41 lbs.	57 lbs.	35.5 lbs.	33.5 lbs.	38 lbs.
Weight Capacity	250 lbs.	250/250/300 lbs.	450/450 lbs.	250 lbs.	250 lbs.	250 lbs.
Warranty						
Frame	3 Years	3 Years	3 Years	3 Years	Lifetime Limited	Lifetime Limited
Tires	None	None	None	None	None	None
Upholstery	None	None	None	None	None	None
Other Components	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year
Reimbursement Code ¹	E1038	K0001	K0006/K0007	K0003	K0004	K0004 & K0028
Eligibility For No Charge Upgrade			K0006	K0001/K0002	K0001/K0002/K0003	

1. Meets requirements for HCPCS code

 Breezy

BREEZY EC WHEELCHAIRS

The best combination of
durability, selection and value.





Breezy EC 2000

A quality wheelchair with options offering the flexibility to accommodate a variety of users.

Standard Features

- Dual axle allows the seat to be easily lowered to hemi height
- Chrome-plated, carbon steel frame is built to last, great for rental fleets
- Composite mag wheels are molded with contours to enhance strength and durability
- Full size chart pocket
- Industry-compatible front riggings
- 250 lb. weight capacity (300 lb. on 20" seat width)

Options for Versatility

- Removable desk-length or full-length or fixed full-length armrests
- Swing-away footrests, elevating legrests

Models

216RADPS	16" width, removable padded desk arms, swing-away footrests
216RADPE	16" width, removable padded desk arms, elevating legrests
218FAFPS	18" width, fixed padded full arms, swing-away footrests
218RADPS	18" width, removable padded desk arms, swing-away footrests
218RADPE	18" width, removable padded desk arms, elevating legrests
218RAFPS	18" width, removable padded full arms, swing-away footrests
218RAFPE	18" width, removable padded full arms, elevating legrests
220RADPS	20" width, padded desk arms, swing-away footrests
220RADPE	20" width, padded desk arms, elevating legrests
220RAFPS	20" width, padded full arms, swing-away footrests
220RAFPE	20" width, padded full arms, elevating legrests



Breezy EC 2000HD Heavy-Duty

Extra wide seat widths and reinforced construction accommodate special needs.

Features

- Heavy-duty reinforced chrome steel frame supports weights up to 450 lbs.
- Seat widths of 22" or 24" accommodate larger users.
- Durable, double-embossed vinyl upholstery
- Swing-away footrests or elevating legrests with padded calf
- Oversized chart pocket
- Choice of aluminum swing-away footrests or elevating legrests and padded full or desk armrests

Models

2HD22RADPS	22" width, padded desk arms, swing-away footrests
2HD22RADPE	22" width, padded desk arms, elevating legrests
2HD24RADPS	24" width, padded desk arms, swing-away footrests
2HD24RADPE	24" width, padded desk arms, elevating legrests



Breezy EC® Wheelchairs

As the world's leading manufacturer of custom mobility products, Sunrise Medical knows wheelchairs.

We also know that for basic patient transportation needs your requirements are different. You need wheelchairs that have been **designed with durability** in mind to provide years of trouble-free service at a cost that maximizes the return on your investment.

You also require choices—**choices** in style, size, features and configurations to accommodate all of your customers' needs, as well as a choice of suppliers.

Seven Breezy EC models, covering the spectrum from basic to high-strength, lightweight chairs establish the Breezy EC family as the best combination of **quality, selection and value**.



Breezy EC 4000 High-Strength, Lightweight

Durability in a lightweight design and unique options make this the most versatile EC model.

Features

- Frame design eliminates seat guides and allows for easy installation of custom backs and accessories
- Flip-back arms with an easy release mechanism
- Weighs 33.5 lbs without riggings
- Meets requirements for K0004 Medicare funding
- Height adjustable seat back (16"/17"/18")

Models

416RADPS	16" width, padded desk arms, swing-away footrests
418RADPE	18" width, padded desk arms, elevating legrests
418RADPS	18" width, padded desk arms, swing-away footrests
420RADPS	20" width, padded desk arms, swing-away footrests



Breezy EC 3000 Lightweight

A value-priced lightweight design to meet today's funding requirements.

Features

- Weighs 35 lbs without riggings meeting requirements for K0003 funding
- Scratch-resistant composite handrim
- Durable chrome-plated, carbon steel frame
- Dual axle allows seat to be easily adjusted to hemi height
- Removable desk arms with padded armrests
- Comfortable padded nylon upholstery
- Industry-compatible swing-away footrests or elevating legrests

Models

316RADPS	16" width, swing-away footrests
316RADPE	16" width, elevating legrests
318RADPS	18" width, swing-away footrests
318RADPE	18" width, elevating legrests
320RADPS	20" width, swing-away footrests
320RADPE	20" width, elevating legrests



Breezy EC Transport

Compact and lightweight, this economical model offers the most convenience for transporting.

Standard Features

- Painted black carbon steel offers clean lines in a durable frame
- Folding back allows chair to fold into a compact unit for easier transport and storage
- Padded nylon upholstery is comfortable and durable
- 8" wheels in front and back provide a smooth ride
- Full-length arms and a nylon seat belt are standard
- Contoured handgrips and easy-to-access wheel locks offer convenience to the caregiver
- Weighs only 24 lbs. without footrests
- Swing-away footrests with heel loops allow for easy transfers

Models

ECTC17SB	17" seat width, black steel frame
ECTC19SB	19" seat width, black steel frame

Lightweight Features

- Lightweight model offers even greater maneuverability weighing only 18.4 lbs.
- Red anodized aluminum frame
- All the features of the standard model

Models

ECTC17ALR	17" seat width, red aluminum frame
ECTC19ALR	19" seat width, red aluminum frame



Standard



Lightweight

Breezy EC 4000 Recliner

Reclining backrest offers patient comfort and re-positioning.

Features

- Full back and headrest support
- Reclines from 90 to 160 degrees
- Removable 10" headrest
- Swing-away footrest models have removable desk armrests
- Elevating legrest models have height-adjustable full armrests
- High-strength, lightweight frame is easily maneuverable
- Dual-axle can be easily lowered to hemi seat height

Models

4R18RADPS	18" width, reclining back, padded desk arms, swing-away footrest
4R16HAFPE	16" width, reclining back, height adjustable padded full arms, elevating legrests
4R18HAFPE	18" width, reclining back, height adjustable padded full arms, elevating legrests
4R20HAFPE	20" width, reclining back, height adjustable padded full arms, elevating legrests



Breezy EC Accessories



	Breezy EC Transport	Breezy EC 2000	Breezy EC 2000HD(20-24")	Breezy EC 2000HD(26-30")	Breezy EC 3000	Breezy EC 4000	Breezy EC Recliner
Overhead Anti-Theft Device	N/A	GEY-OATH	GEY-OATH	N/A	GEY-OATH	GEY-OATH	N/A
Rear Anti Tipppers	N/A	GEY-ATL	GEY-ATL	GEYB-AT	GEY-ATL	GEY-4AT	Included
Anti-Fold Device	N/A	GEY-ATH	GEY-ATH	N/A	GEY-ATH	GEY-ATH	GEY-ATH
Oxygen Tank Holder ^{2 3}	N/A	GEY-O2	GEY-O2	GEYB-O2	GEY-O2	GEY-O2	N/A
Oxygen Tank/IV Holder ^{2 3}	N/A	GEY-O2IV	N/A	N/A	N/A	N/A	N/A
Telescoping IV Pole	N/A	GEY-IVP	GEY-IVP	GEYB-IVP	GEY-IVP	GEY-IVP	N/A
Wheel Lock Extension Handles	N/A	GEY-EXT	GEY-EXT	N/A	GEY-EXT	GEY-EXT	GEY-EXT
Amputee Axle Plate ⁴	N/A	GEY-AMP	GEY-HDAMP	N/A	N/A	GEY-4AMP	N/A
Amputee Axle Plate Wheel Lock	N/A	GEY-WDG	GEY-WDG	N/A	N/A	GEY-WDG	N/A
Chart Pocket	N/A	Included	Included	N/A	GEY-CP	GE-CP	GE-CP
Footrests	ECTCAR-FTREST (steel) ECTCSB-FTREST (alu)	GEY-FRC	GEY-FRA	GEYB-FRA	GEY-FRC	GEY-FRC	GEY-FRC
Footrests w/Heel Loop	N/A	GEY-FRCHL	N/A	N/A	GEY-FRCHL	GEY-FRCHL	GEY-FRCHL
Elevating Legrest	N/A	GEY-ELRC	GEY-ELRA	GEYB-ELRA	GEY-ELRC	GEY-ELRC	GEY-ELRC
Low Profile ELR	N/A	GEY-ELRCLP	N/A	N/A	N/A	GEY-ELRCLP	GEY-ELRCLP
Articulating Legrests	N/A	GEY-ALR	GEY-ALR	N/A	GEY-ALR	GEY-ALR	GEY-ALR
Ht. Adj Padded Desk Armrest	N/A	GEY-HADP	GEY-HADP	N/A	GEY-HADP	GEY-4HADP	GEY-HADP
Ht. Adj Padded Full Armrest	N/A	GEY-HAFP	GEY-HAFP	N/A	GEY-HAFP	GEY-4HAFP	GEY-HAFP
2" Seat Extension	N/A	GEY-2SREXT* * = 16, 18 or 20	N/A	N/A	N/A	GEY-4SREXT* * = 16 or 18	N/A
Seat Belt (Automotive Style)	Included	GEY-SBA	GEY-SBA	N/A	GEY-SBA	GEY-SBA	GEY-SBA
Seat Belt (Velcro Style)		GEY-SBV	GEY-SBV	N/A	GEY-SBV	GEY-SBV	GEY-SBV
Super Hemi ⁵ (15.75")STFH	N/A	N/A	N/A		N/A	GEY-SHEMI	N/A

2. Must use with rear anti-tippers
3. Mounts on left only
4. Must use with GEY-WDG
5. Includes 20" wheels, 6" casters & fork, GEY-WDG wheel-locking extension



Dual axle
Allows seat to be easily lowered to hemi height. Available on 2000, 3000 and 4000 models



Chart pocket
Included on 2000 and 2000HD models.



Quick-release flip-back arms
Make transfers easy. Available on 4000 models.



Height adjustable back
Included on 4000 models.



JAY® Seating Options


Make sure your customers get the comfort they deserve.

** Medicare beneficiaries that qualify for a wheelchair also qualify for a general use seat cushion (E2601/E2602) and a general use wheelchair back cushion (E2611/E2612). The use of these products provides better sitting tolerance and comfort.*

JAY® Basic Cushion

The Jay Basic provides a soft and durable support for clients with mild postural needs who are at low risk for skin breakdown. The Jay Basic Cushion has been code verified by SADMERC to meet the requirements of E2601/E2602.

- Beveled to fit between the seat rails on wheelchairs with seat upholstery.
- Durable molded foam provides mild contouring for increased sitting tolerance and postural stability
- Provides stability to reduce sliding and shearing
- Available for chair widths 14” to 24”.



JAY® Basic Back

The Jay Basic Back has been redesigned to easily adapt to standard wheelchairs. Choose strap-style hardware for easy attachment - just strap around wheelchair back canes. The clamp-style hardware provides 0° and 8° of angle adjustability and snaps easily into place.

- Designed for clients with mild postural needs
- Provides mild posterior pelvic-stabilization and thoracic support
- Three height options provide a more personalized fit
- Lightweight for improved propelling efficiency and easy transport
- Lightweight, contoured, aluminum shell and soft foam insert



* An order for each item billed must be signed and dated by the treating physician, kept on file by the supplier, and made available to the DMERC upon request. It is advised that the order for the seat and back cushions be obtained at the same time as the order for the wheelchair to avoid claim processing errors.