

Child Restraint System **CARROT 3**



- **Adaptive Positioning**

High density foam inserts for optimal positioning and maximum comfort

- **Grows With Your Child**

30-108 lbs / 37" – 60"

- **"Free Angle" Recline 180°**

Allows for easy diaper changes and conforms to vehicle's seat

- **LATCH Equipped**

Lower Anchors and Tethers for Children

Accessories:

Seat Angle Adjustment Wedge



Helps stabilize sitting posture

Lateral Support



Maintain midline positioning of the trunk for increased postural control

Modular Extensions



Extensive growth where needed – length & height

Support Tray



Provides added support to the upper body (Crash Tested)

Lightweight • Compact Folding • Convenient

 **Convoid**
Now You're Going Places



888-Convoid (266-8243)
www.convoid.com

Small Support Kit



Chest Support



Head Support

Accessories



Support Tray
(Crash Tested)



Seat Angle
Adjustment Wedge



Air Mesh
Cushion



Tether Strap
(Recommended)

Foot & Leg Positioning



Rodded Footrest
(Small, Medium, Large)


Pelvis Positioning



Rodded Seat Extension
(Large Curved)



Rodded Seat Extension
(Short 90°)



Seat Extension
(2" or 4")



Pelvis Lock
Harness
(S, M or L)

Trunk Positioning



Backrest Extension
(2" or 4")



Firm Lateral Trunk
Support



Shoulder Protector
Wings



3D Body
Balance Pads


Upholstery Colors



Red



Blue



Black

Carrot 3 Base Model Specifications:

Seat Width:	11.5"
Seat Depth Adjustment:	9.8" - 11.3"
Backrest Height Adjustment:	11.8" - 12.8"
Head Support Height:	8.5"
Recline Adjustment:	90 - 180°
Weight:	19.1 lbs.

Carrot 3 Specifications with Extensions:

Seat Depth Adjustment:	9.8" - 22.3"
Backrest Height Adjustment:	11.8" - 22.8"
Shoulder Wing Height:	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"
Lateral Supports Width Adjustment:	10" - 12"
Lateral Supports Height Adjustment:	2"

* Measurements do not include head support & growth extensions

^ Base model includes 1.5" of seat depth growth **A** & 1" of backrest growth **B**

Base Model Measurements: * ^

Seat Depth:	9.8" - 11.3"
Backrest Height:	11.8" - 12.8"
Head Support Height:	8.5"

* Measurements do not include head support & growth extensions

^ Base model includes 1.5" of seat depth growth **A** & 1" of backrest growth **B**

Growth Extensions:

1a Rodded seat extension - large	7"
1b Rodded seat extension - small	2"
2 Seat extension - 4"	4"
3 Seat extension - 2"	2"
4 Backrest extension - 2"	2"
5 Backrest extension - 4"	4"
6 Shoulder protector wings	6"

Note: The same type of 2" & 4" extensions cannot be used simultaneously.

User Measurement Guide:

Minimum User Weight:	30 lbs.
Maximum User Weight:	108 lbs.
Minimum User Height:	37"
Maximum User Height:	60"
Recommended User Age:	Ages 3-15

