


COMPARE TO MASON® GEL-U-SEAT™,
BLUE CHIP® GEL PRO DLX™, INVACARE®
COMFORTMATE, INVACARE® - INFINITY LO
CONTOUR & FLA® POSTURA® GELFOAM™



14886

COMPARE TO MASON® -
TITANIUM, JAY® CARE, PRIDE® -
SYNERGY, BLUE CHIP - AMARA 400
& INVACARE® - INFINITY
GENTLE CONTOUR



14882

Gel "E" 3" Wheelchair Cushions with Gel Bladder

- 14886 18" (W) x 16" (D) x 3" (H), 275 Wt Cap., 1/cs
(Formerly 14890)
14891 16" (W) x 16" (D) x 3" (H), 275 Wt Cap., 1/cs
14894 20" (W) x 16" (D) x 3" (H),
275 Wt Cap., 1/cs

HCPCS code: E2603

- 14902 22" (W) x 18" (D) x 3" (H), 325 Wt Cap., 1/cs
14903 24" (W) x 18" (D) x 3" (H), 350 Wt Cap., 1/cs

HCPCS code: E2604

FEATURES:

- Pressure reducing cushion designed to assist in the prevention, treatment and management of pressure ulcers.
- Viscous Gel provides optimized pressure reduction, support and comfort.
- Fire Retardant polyurethane foam shell.
- Removable and replaceable cover has water resistant "Top" and "Non-Slip" vinyl base.
- 18 Month Limited Warranty on cover and cushion.

Molded - Skin Protection & Positioning Wheelchair Cushions with Polymer Filled Bladders & Abductor

- 14882 18" (W) x 16" (D) x 3" (H), 250 Wt. Cap., 1/cs
14883 20" (W) x 16" (D) x 3" (H), 250 Wt. Cap., 1/cs

HCPCS code: E2607

- 14884 22" (W) x 18" (D) x 3" (H), 350 Wt. Cap., 1/cs

HCPCS code: E2608

FEATURES:

- Unique construction of high resilient (HR) molded foam provides optimum seat stabilization, positioning, pressure reduction and correct seating posture.
- Unique Polymer bladders beneath the Ischium for pressure reduction.
- Soft, durable, sealed inner liner provides protection from incontinence.
- Removable and replaceable cover is washable, flame retardant, reduces pressure, reduces friction & shearing, is urethane coated for water resistance and vapor permeability, is constructed with a non-slip bottom
- 1 Year Limited Warranty on cover and cushion.

SPECIFICATIONS

Item #	Width	Depth	Height	Product Weight	Carton Shipping Weight	Weight Limit
14886	18"	16"	3"	5.7 lbs.	6.7 lbs.	275 lbs.
14891	16"	16"	3"	4.4 lbs.	5.4 lbs.	275 lbs.
14894	20"	16"	3"	4.4 lbs.	5.4 lbs.	275 lbs.
14902	22"	18"	3"	12.5 lbs.	14.0 lbs.	325 lbs.
14903	24"	18"	3"	14.0 lbs.	15.5 lbs.	350 lbs.
14882	18"	16"	3"	3.9 lbs.	4.9 lbs.	250 lbs.
14883	20"	16"	3"	4.4 lbs.	5.4 lbs.	250 lbs.
14884	22"	18"	3"	8.1 lbs.	9.6 lbs.	350 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.