

HTR Accessories

Substitute XX in part numbers below with number corresponding to your color choice: 25-Mauve, 26-Grotto and 27-Blueberry. Part numbers for the accessories include the hardware kits to attach to the chairs



Headrest Extension
(Option No. HRA1)
(Part No. 460069XX)
Dimensions 6" x 15"
Adjustment Range 3"
Hardware included



4" Abduction Pommel
(Option No. HRA2)
(Part No. 460070XX)
Hardware included
Not available on
HTR3000/5000



2" Side Bolsters (Pair)
Option No. HRA10
Part No. 450148XX
Dimensions 6" x 15"
Reduces seat width



Lateral Trunk Support (Pair)
(Option No. HRA20)
(Part No.:
1072292 - Mauve
1072299 - Grotto
1072300 - Blueberry)
Dimensions 4.5" x 6"
Adjustment Range 4"
Hardware included



Seat Positioning Strap
- Auto Style
(Option No. HRA30)
(Part No. 450149)
Hardware included



Adjustable Height IV Pole
(Option No. HRA35)
(Part No. 1072293)
Adjustment Range 35" to 61"
Hardware included



Standard Lap Tray
(Option No. HRA42)
(Part No. 460087 - 16")
(Part No. 460088 - 18")
Dimensions:
16" wide - 23" x 18"
18" wide - 23" x 20"
Standard on HTR3500
and HTR5500



Vinyl Tray Pad 1"
(Option No. HRA4CA)
(Part No. 460071XX - 16")
(Part No. 460072XX - 18")
Secures with Hook and
Loop Straps



Contoured Seat
Option No. HRA43
Part No. 450088-XX - 16"
Part No. 450089-XX - 18"
Standard on HTR3500 and
HTR5500



Headwings (Pair)
(Option No. HRA44)
(Part No.:
1072294 - Mauve
1072301 - Grotto
1072302 - Blueberry)
Dimensions 7" x 11.5"
Adjustment Range 4"
Hardware included
Standard on HTR3500
and HTR5500



Winged Headrest
(Option No. HRA46)
(Part No. 450127-XX)
Secures with hook and
loop straps



Contoured Headrest
(Option No. HRA47)
(Part No. 450094-XX - 16")
(Part No. 450095-XX - 18")
Secures with hook and
loop straps



Front Anti-Tippers (pair)
(Option No. HRA60)
(Part No. 460073)
Add 9" to chair depth



Footboard Pad
(Option No.:
HRA68-16 - 16")
HRA68-18 - 18")
(Part No. 450121 - 16")
(Part No. 450122 - 18")
Secures with hook and
loop straps



12" Wheel with Flat
Free Inserts
(Option No. HRA12)
Optional on HTR5000
and HTR5500
Width reduction of 4"



Swing Away Footrests
(Option No. HRA92)
5.5" Horizontal adjustment
5" Footboard adjustment
15" Total additional depth



Articulating Elevating
Leg Rest
(Option No. HRA95)
5.5" Horizontal adjustment
13" Total footboard
adjustment
29" Total additional depth



One-Piece Elevating Footboard
(Option No. HRA93)
5.5" Horizontal adjustment
7.5" Footboard adjustment
23" Total additional depth
Standard on HTR3500
and HTR5500