



| Item # | Seat | UOM |
|-----------------|--------------------|------|
| 2800ECBU-RCL | 18" Captain's Seat | 1/cs |
| 2800ECBU-RCL-20 | 20" Captain's Seat | 1/cs |



HPCPS
K0823



Image EC

Standard Power Wheelchair
Mid-Wheel Drive

- Includes interchangeable color shrouds in red and blue.
- Front suspension system for extra-smooth ride.
- Available with full-size Captain's Seat in 18" or 20" widths.
- Height-adjustable swivel seat features fold-down backrest.
- Seat back reclines and features adjustable headrest.
- Adjustable-length controller mount.
- Armrests are padded, adjustable and removable.
- Optional swing-away controller arm (Item #SF1009) sold separately.
- Adjustable, flip-up footplate.
- Positioning belt included.
- Optional elevating leg rests (requires Item #ECLER-P, Item #SFELR & Item #LK3JELR) sold separately.
- Flat-free, non-marking tires.
- Easy freewheel operation.
- Easy access to batteries from the rear of the chair. (Figure A)

| specifications | | | |
|----------------------|---------------------------------------|--|--------------------|
| HPCPS | K0823 | Batteries (pair) | 12V x 35AH (U1) |
| Weight Capacity | 300 lbs. | Charger | 4A Offboard |
| Top Speed | 4 mph | Brakes | Electromagnetic |
| Maximum Range | 15 miles | Caster Wheels | 6" x 2" Solid |
| Turning Radius | 27.5" | Drive Wheels | 10" x 3" Flat-free |
| Climbing Angle | 6° | Base Weight | 88 lbs. |
| Ground Clearance | 2" | Battery Weight (pair) | 53 lbs. |
| Dimensions | 42" (L) x 24" (W) | Seat Weight | 51 lbs. |
| Floor-to-Seat Height | 20" - 23" | Total Weight | 192 lbs. |
| Seat Dimensions | 18" (W) x 18" (D) 20" (W) x 18" (D) | Warranty on Frame | Lifetime |
| Controller | PG VR2 50A | Warranty on Electronic Controller/Drive Train Components | 14 Months |
| Motors (2) | 24V x 200W x 3900 rpm | Warranty of Batteries | 6 Months |