

The Jazzy® 614 HD features patented Mid-Wheel 6® technology and Active-Trac® ATX Suspension for an outstanding level of performance, maneuverability and tight turning radius, all in one heavy-duty package. All-around performance indoors and out over most surfaces makes the Jazzy 614 HD an excellent heavy-duty power chair choice.

ACTIVE-TRAC® ATX SUSPENSION





6" CASTERS



EXTENDED SEATING SIZES

AMERICA'S FAVORITE POWER CHAIR

Learn more online.









go**jazzy**.com

Specifications

Weight Capacity	450 lbs.
Maximum Speed	Up to 4.25 mph
Ground Clearance	2.25"
Turning Radius	22.75"
Overall Length	35.5"
Overall Width	25.25"
Drive Wheels	14" solid
Front Wheels	6" solid casters
Rear Wheels	6" solid casters
Suspension	Active-Trac® (ATX)
Drivetrain	Two-motor, in-line, Mid-Wheel 6®
Braking System	Regenerative and electro-mechanical
Standard Electronics	50A, Dynamic Shark controller
Battery Charger	Off-board, 5 amp
Per-Charge Range	Up to 11 miles
Battery Requirements	(2) 12 volt, NF-22
Battery Weight	37 lbs. each
Base Weight	122 lbs.
Standard Seat Weight	52 lbs.
Standard Seat	20" x 18-20", 22" x 20"-22"
Max. Seat Dimensions	24" x 24"
Warranty	Lifetime limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

Colors







Viper Blue

Models

Jazzy 614HD 2HD-SS	Code: K0824
Jazzy 614HD 2HD-C	Code: K0825

^{*}Additional charges may apply.

Features and Benefits

- Depth adjustable, limited recline high back seat available
- · Universal, fixed, in-line joystick mount
- · Flip-up height adjustable foot platform
- Lap belt (50", 60", 70", 80", 90")

Seating Options*

- Deluxe contour, limited recline high back available sizes: up to 24"x24"
- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 16"x16" or 20"x22"
- Synergy[®] Seat

Options

- · Elevating leg rests
- Heel loops
- · Residual limb support
- Swing- away footrests
- · Swing-away joystick

Accessories

- · Cane/crutch holder
- Cell phone holder
- · Cup holder
- Flag holder
- · Military patches
- · Oxygen tank holder
- Rear basket
- Saddle bag
- · Walker holder
- · Weather cover
- · XLR USB charger

Recommended Lifts

- Backpacker® Series
- Outlander Full Platform
- · Commander 450



Find a compatible lift for your vehicle online: PRIDEMOBILITY.COM/LIFT-FINDER



(US) 800-800-8586 • (Canada) 888-570-1113 • pridemobility.com •









The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.