



SUPERIOR PERFORMANCE
INDOORS AND OUT

Jazzy® Select™ 6

The Jazzy® Select™ 6 provides maximum front stability, combined with the advantages of in-line motor technology. Active-Trac® Suspension and 6" front and rear casters for superior performance indoors and out. An impressive selection of standard convenience features make it the easiest power chair to use.

#1 Physician Approved

JAZZY® SELECT™ 6
shown in Jazzy Red with PG
GC controller and foot platform


New contemporary
design high-back seat



Charger port
relocated under
controller



New redesigned foot platform



New swivel handle

Pride
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

JZ SEL 6 SS 04/07



The Jazzy® Select™ 6 series models include:

- Jazzy Select 6-C
- Jazzy Select 6-SS

FEATURES

- Active-Trac® Suspension
- 6" casters on front and rear for maximum stability
- In-line motor technology for increased efficiency, torque & range
- 40 amp, PG GC controller
- Articulating rear caster beam
- Battery removal through side panel
- Easy access charger port and circuit breaker
- Service friendly access to controller
- Height and depth-adjustable foot platform

OPTIONS

- High-back with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.

Jazzy® Select™ 6	
Weight capacity:	300 lbs.
Maximum speed:	Up to 4 mph
Ground clearance:	2"
Turning radius:	20.5"
Overall length:	35.625" without foot riggings
Overall width:	23.5"
Wheels:	Drive: 10" solid
Front casters:	6" solid
Rear casters:	6" solid
Suspension:	Active-Trac®
Drive train:	Two-motor, Mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	40 amp, PG GC controller
Specialty controls:	N/A
Battery charger:	Off-board, 5 amp
Per-charge range:	Up to 15 miles
Battery requirements:	(2) 12 volt, deep cycle
Battery size:	U-1
Transport weight	
Battery weight:	24.5 lbs. each
Base weight:	100 lbs.
Max. Seat dimensions:	20" x 20"
Warranty:	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-months limited warranty on drive motors

Colors:



Pride SilverStar Lifts ★

- Full Platform
- Power chair lift with straps
- Backpacker™ Plus



Backpacker™ Plus



Full Platform

These lifts are compatible with the Jazzy® Select™ 6.