



Item #	Product Description	UOM
KINGCOBRA422CS-PGV	22" Sport-Style Captain's Seat	1/ea



King Cobra PGV Executive Scooter

- Personal Golf Vehicle with Hard Top Canopy and Golf Bag Attachment.
- Intuitive Drive System for Ultimate Stability.
- Front and Rear Suspension for the Smoothest Ride Possible.
- 6" Wide Turf Tires to Handle All Types of Terrain.
- Built-in USB Port for Charging your Mobile Devices.
- Built-in Cup Holder to Conveniently Carry Your Drinks.
- Storage Compartment Built into the Tiller.
- One-touch, infinite tiller adjustment allows the user to easily position the tiller into the ideal driving position.
- Ergonomic Delta Control Handles for Ease of Steering and Throttle Control.
- 22" W x 20" D Sport Style Captain's Seat is height adjustable and features a fold-down, reclining backrest.

specifications

Weight Capacity	450 lbs.	Charger	8A Offboard
Top Speed	10 mph	Brakes	Electromagnetic & Hand
Maximum Range	35 miles	Front Wheels	14.5" x 6"
Turning Radius	83"	Rear Wheels	17" x 6"
Climbing Angle	8°	Color	Vivid Red
Ground Clearance	4"	Base Weight	209 lbs.
Dimensions	85" (L) x 29" (W)	Battery Weight (pair)	94 lbs.
Floor-to-Seat Height	27" - 30"	Seat Weight	60 lbs.
Seat Type	Sport Style Captain's	Total Weight	363 lbs.
Seat Dimensions	22" (W) x 20" (D)	Warranty on Frame	Limited Lifetime
Controller	Dynamic Rhino II 180A	Warranty on Electronics	14 Months
Motor	24V x 1300W x 5200 rpm	Warranty on Batteries	6 Months
Batteries (pair)	12V x 75AH		