WWW.DRIVEMEDICAL.COM

EASILY CONVERTS FROM A TRANSPORT CHAIR TO A SELF PROPELLED WHEELCHAIR.



Poly-Fly High Strength, Light Weight Wheelchair/Flyweight Transport Chair Combo

TR16 16" Seat Width, Silver Frame and Black Upholstery, 1/cs
TR18 18" Seat Width, Silver Frame and Black Upholstery, 1/cs
TR20 20" Seat Width, Silver Frame and Black Upholstery, 1/cs

LELR-TF Swing-away Elevated Legrest, Tool Free,

(Optional front riggings), 1pr/cs

FEATURES:

- Can be used as standard self propelled wheelchair or a transport chair all in one.
- Quick release 24" wheels can be removed to transition self propelled chair to a transport chair.
- All aluminum frame.
- Comes with two sets of aluminum wheel locks: One for use with wheelchair, one for use with transport chair.
- Fold down back. Comes with Deluxe Back Release.
- Nylon upholstery.
- Comes standard with swing away foot rests with composite footplates and heel loops.
- 8" casters
- Limited Lifetime Warranty.

EASILY FOLDED FOR STORAGE!

QUICK RELEASE ON 24" WHEEL.





SPECIFICATIONS

or Edit Idaniano												
	Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/ Riggings	Product Weight	Carton Shipping Weight	Weight Limit
	TR16	23.5"	12"	37"	19.5"	16"	16"	16"	42"	33 lbs.	38 lbs.	250 lbs.
	TR18	24.5"	12"	37"	19.5"	16"	18"	16"	42"	34 lbs.	39 lbs.	250 lbs.
Ī	TR20	25.5"	13"	39.5"	19.5"	16"	20 "	16"	41 "	36 lbs.	41 lbs.	250 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

^{*}Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.