

# LXI/LX



Quickie's traditional folding chair with value pricing, a broad range of options,durable design and style.





improving people's lives

## Nothing Beats a Quickie<sup>™</sup>

Both Quickie® LX I and LX offer multiple features and options in lightweight manual folding wheelchair designs. These dependable, versatile chairs provide outstanding ease-of-use and performance at very affordable prices.



### **Compatibility:**

- Several common parts to Quickie ultra lightweight folding family chairs
- Helps reduce inventories and interchange parts when necessary

• Compliments Quickie's strong tradition of rehab

New Angle Adjustable Backrest:

the pivot point of the spine

focused chair options

Value Options Standard:

• Camber-angle adjustable axle plate

• 15 frame colors • Swing-In/out hangers



Angle







#### **NEW Options:**

- Transit Option
- Spinergy® wheels

SPECIFICATIONS	LXI (K0005)	LX (K0004)
Frame Style	Folding – Traditional	Folding – Traditional
Transport Weight	16.2 lbs.	16.2 lbs.
Overall Weight	28 lbs.	28 lbs.
User Weight Capacity	265 lbs.	265 lbs.
Seat Widths	12" - 20" (20" upcharge)	12" - 20" (20" upcharge)
Seat Depth	12" - 20"	12" - 20"
Back Angle Adjustment	Yes	Yes
Front Seat To Floor Height	16" - 20.5"	16" - 20.5"
Rear Seat to Floor Height	15" - 19.5"	15" - 19.5"
Footrest/Frame Angle	Swing-away	Swing-away
WARRANTY		
Warranty	Limited/Lifetime	Limited/Lifetime





For additional product specifications please visit: www.sunrisemedical.com

 Mountain wheels • Case Logic® side pockets

- Camouflage frame colors
- Natural Fit® handrims
- Military service patches.