

**SEAT ROTATES
360° AND LOCKS
EVERY 90°!**



Bath Stool with Padded Rotating Seat

RTL12061 4 retail boxes/cs

Available
1 Pack

FEATURES:

- The rotating padded seat reduces twisting and reaching while in the shower.
- The rotating top of the seat can be removed for easy cleaning.
- Comes with removable tray for storing personal items.
- Ideal for individuals who have balance or mobility issues.
- Seat Height can be adjusted from 16" to 23".
- Weight capacity 400 lbs.
- Limited Lifetime Warranty.



**SUPPORT
COLLAR
REDUCES
RATTLE!**



**RTL12004KD
NEW PINCH FREE COVER
ALLOWS FOR PUSH PINS
TO BE DEPRESSED
WITHOUT PINCHING
FINGER. THE AUDIBLE
"CLICK" WILL BE HEARD
WHEN SECURE!**

K.D. Bath Stool

RTL12004KD Knocked Down, White,
(Has replaced 12004KD-1), 1/cs
12004KDRB-1 Knocked Down Retail, Blue,
(Has replaced 12004B-2), 1/cs
12004KDRT-1 Knocked Down Retail, Teal,
(Has replaced 12004T-2), 1/cs
12004KDRY-1* Knocked Down Retail, Yellow, 1/cs
12002T Replacement Tips for 12204, 4/bx

Available
1 Pack

FEATURES:

- Aluminum frame is lightweight, durable and corrosion proof.
- Impact resistant, composite seat
- Height adjustable legs. is crack-proof and tarnish resistant.
- Height adjustable legs.
- Support collar reduces rattle.
- Limited Lifetime Warranty.



12004KDRY-1



12004KDRT-1



12004KDRB-1

SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
RTL12061	15"	14"	15"	14"	23"	16"	7.5 lbs.	33 lbs.	400 lbs.
RTL12004KD	13.5"	13"	13.5"	13"	17"	13.5"	3.5 lbs.	5.5 lbs.	300 lbs
12004KDRB, T, Y-1	13.5"	13"	13.5"	13"	17"	13.5"	3.5 lbs.	5.5 lbs.	300 lbs

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

To order call toll free: 877.224.0946

BATH SAFETY

drive