



Shown: Double Entry Scooter Lift

SCOOTER EXTERIOR LIFTS

Add a whole new dimension to your mobility with a Silver Star Scooter Lift by Pride, the reliable solution to taking your scooter on the road. An all-steel frame and protective coating mean all Silver Star Lifts are built to be durable, while quality construction means you can depend on your lift for the long haul.

- Single Entry, Double Entry and Quad XL lifts available
- 3-year transferable warranty
- 325 lbs. lifting capacity

SilverStarLifts ★

A Division of **Pride** Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • (800) 800-8586 • www.silverstarmobility.com • www.pridemobility.com

SCOOTER EXTERIOR LIFTS



Double Entry Scooter Lift



Single Entry Scooter Lift



Quad XL Scooter Lift

DOUBLE ENTRY SCOOTER LIFT

The Double Entry Lift allows you to enter and exit the lift platform from either side with your 3 or 4-wheel Pride® Scooter.

Base weight lift	106 lbs. ¹
Maximum scooter weight	325 lbs.
Maximum product width	27"
Maximum travel distance from the ground	20"
Motor	Sealed 12 volt DC
Optional battery pack	12 volt, on-board
Optional charger	100 volt
Hitch class	Class III or IV ³
Warranty	3-year transferable

SINGLE ENTRY SCOOTER LIFT

Silver Star Single Entry Scooter Lift by Pride® enables you to easily lift your Pride 3-wheel scooter.

Base weight lift	112 lbs. ¹
Maximum scooter weight	325 lbs.
Maximum product width	27"
Maximum travel distance from the ground	20"
Motor	Sealed 12 volt DC
Optional battery pack	12 volt, on-board
Optional charger	100 volt
Hitch class	Class II or III ³
Warranty	3-year transferable

QUAD XL SCOOTER LIFT

The Silver Star Quad XL Scooter Lift by Pride® is sturdy, longer scooter lift for use with Pride's larger products, the Hurricane® and Legend® XL.

Base weight lift	112 lbs. ¹
Maximum scooter weight	325 lbs.
Maximum product width	27"
Maximum travel distance from the ground	20"
Motor	Sealed 12 volt DC
Optional battery pack	12 volt, on-board
Optional charger	100 volt
Hitch class	Class II or III ³
Warranty	3-year transferable

¹ Weight does not include any accessories

² Vehicle dependent

³ 190 lb. maximum scooter weight with Class II hitch

Note: Silver Star Lifts, while designed specifically for Pride Mobility products, are capable of lifting other scooter brands.

The Double Entry and Single Entry Scooter Lifts can easily lift the following Pride® Scooters ²

- Celebrity® X
- Go-Go® Ultra
- Legend®
- Maxima
- Rally®
- Revo®
- Sonic®
- Sundancer®
- Victory®

The Quad XL Scooter Lift can easily lift the following Pride® Scooters ²

- Hurricane®
- Legend® XL

SilverStarLifts

A Division of **Pride** Mobility Products Corp.

182 Susquehanna Ave. • Exeter, PA 18643 • (800) 800-8586 • www.silverstarmobility.com • www.pridemobility.com