

Super Light, Folding Transport Chair with Carry Bag

SL18 Silver Frame and Blue Upholstery, 1/cs
 839* Replacement Carry Bag for SL18, 1/cs
 835N Replacement Carry Bag for SL18
 (Will replace 839), 1/ea

FEATURES:

- Padded, flip back armrests.
- Weighs only 19 lbs. Easy to push or transport.
- Aluminum frame is lightweight and strong.
- Composite wheels are lightweight and maintenance free.
- Folds like a conventional folding chair for added convenience.
- Comes with carry bag (item # 839).
- Fold down footrests.
- Standard rear wheel locks.
- Attractive frame with durable, lightweight nylon upholstery.
- Limited Lifetime Warranty.



SL18

TOOL FREE
REMOVABLE,
FRONT CASTERS

WEIGHT
CAPACITY
250 LBS

FOLDS LIKE A
CONVENTIONAL
FOLDING CHAIR



835N

COMES WITH CARRY BAG



SPECIFICATIONS

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/ Riggings	Product Weight	Carton Shipping Weight	Weight Limit
SL18	26"	26"	37"	20"	17"	18"	17"	30"	19 lbs.	24 lbs.	250 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

transport chairs

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM