



Order form

# ROHO® AGILITY® Max Contour back support



Dealer Name\*:

Account #\*:

City, State, Zip\*:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference\*:

Client Weight\*:

**\*Required field**

Special considerations that need to be addressed here (e.g., diagnosis):

**ROHO® Back supports and accessories** | Back supports include attachment hardware, cover, pump, and operating instructions. Refer to reference pages at end of order form for additional information and images.

## AGILITY® Standard Series Max Contour Back System

HCPCS coding: E2615 Positioning back, posterior lateral

HCPCS coding: E2620 Positioning back, planar with lateral supports

Length	20"	\$957	\$957	\$957	\$957	\$957	\$957	\$957
	18"	\$957	\$957	\$957	\$957	\$957	\$957	\$957
	16"	\$957	\$957	\$957	\$957	\$957	\$957	\$957
		14"	15"	16"	17"	18"	19"	20"
		<b>Width</b>						

Fixed hardware only available with 16" length back systems.

Direct mount hardware compatible with Permobil power wheelchairs w/ compatible Corpus seating.

### Air and lumbar support, hardware options | Choose one.

AG6-WWLLAL-1	<b>Quick release hardware</b>	E2615
AG6-WWLLAL-2	<b>Fixed hardware</b> <i>Only available with 16" length back systems</i>	E2615
AG6-WWLLAL-P	<b>Direct mount hardware</b>	E2615

## Air only, hardware options | Choose one.

AG6-WWLLA-1	<b>Quick release hardware</b>	E2620
AG6-WWLLA-2	<b>Fixed hardware</b> <i>Only available with 16" length back systems</i>	E2620
AG6-WWLLA-P	<b>Direct mount hardware</b>	E2620

## Head support options

Must select size of pad and hardware

### Therafin Headrest and hardware options for ROHO AGILITY Standard Series Back Support Systems

AG-HRPAD-SM	<b>Small (5" height x 8" width)</b>	\$199	E0955
AG-HRPAD-MED	<b>Medium (5" height x 10" width)</b>	\$199	E0955
AG-HRPAD-LG	<b>Large (5" height x 14" width)</b>	\$199	E0955
AG-HRPAD-BRKT	<b>Multi-Axis headrest hardware</b> <i>Required</i>	\$244	E1028
AG-HR-BRKT	<b>Universal headrest adaptor plate</b>	\$107	
AG-HR-UT	<b>Headrest bracket and unitrack to fit 15" - 20" width</b>	\$21	
AG-HR-UT14	<b>Headrest bracket and unitrack to fit 14" width</b>	\$21	

## Other options

### ROHO Privacy shield

PS14	<b>10" width x 14" height</b> <i>Fits wheelchair sizes 13" -15"</i>	\$65
PS16	<b>13" width x 14" height</b> <i>Fits wheelchair sizes 15.5" width -17"</i>	\$65
PS18	<b>15" width x 14" height</b> <i>Fits wheelchair sizes 17.5" -19"</i>	\$65
PS20	<b>17" width x 14" height</b> <i>Fits wheelchair sizes 19.5" - 21"</i>	\$65

### Other

AG-LUMBAR	<b>Lumbar pad, AGILITY back</b>	\$23
AG-WRENCH	<b>Allen wrench, 4mm AGILITY back</b>	\$12
HP45	<b>ROHO Inflation Pump</b>	\$10

## Additional Order Instructions (for Permobil):

## Notes & Comments (not for Permobil):

## Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code and the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

## Order Acknowledgement:

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

*NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*

## References

### Section 1 | ROHO AGILITY Standard Series sizing matrix

**A – Back Shell Length**

**A<sup>2</sup> – Back Shell and Foam Length**


**B – Back Cane Width Range**

(there will be a 1.5cm loss of seat depth if an AGILITY Aluminum is fitted to the smallest wheelchair width within each size range)

**C – Inner Width**

**D – Flat Center Section of Shell**

**E – Length of Wings**



### AGILITY Standard Series Max Contour sizing matrix

Number	A	A	A <sup>2</sup>	A <sup>2</sup>	B	B	C	C	D	D	E	E
AG6-1416	40.5	16"	48.5	19"	33-35.5	13" - 14"	28	11"	15	6"	20.5	8"
AG6-1418	45.5	18"	53.5	21"	33-35.5	13" - 14"	28	11"	15	6"	20.5	8"
AG6-1420	51	20"	58.5	23"	33-35.5	13" - 14"	28	11"	15	6"	20.5	8"
AG6-1516	40.5	16"	48.5	19"	35.5-38	14" - 15"	32	12.5"	20.5	8"	20.5	8"
AG6-1518	45.5	18"	53.5	21"	35.5-38	14" - 15"	32	12.5"	20.5	8"	20.5	8"
AG6-1520	51	20"	58.5	23"	35.5-38	14" - 15"	32	12.5"	20.5	8"	20.5	8"
AG6-1616	40.5	16"	48.5	19"	38-40.5	15" - 16"	33	13"	23	9"	20.5	8"
AG6-1618	45.5	18"	53.5	21"	38-40.5	15" - 16"	33	13"	23	9"	20.5	8"
AG6-1620	51	20"	58.5	23"	38-40.5	15" - 16"	33	13"	23	9"	20.5	8"
AG6-1716	40.5	16"	48.5	19"	40.5-43	16" - 17"	37	14.5"	25.5	10"	23	9"
AG6-1718	45.5	18"	53.5	21"	40.5-43	16" - 17"	37	14.5"	25.5	10"	23	9"
AG6-1720	51	20"	58.5	23"	40.5-43	16" - 17"	37	14.5"	25.5	10"	23	9"
AG6-1816	40.5	16"	48.5	19"	43-45.5	17" - 18"	38	15"	29	11.5"	23	9"
AG6-1818	45.5	18"	53.5	21"	43-45.5	17" - 18"	38	15"	29	11.5"	23	9"
AG6-1820	51	20"	58.5	23"	43-45.5	17" - 18"	38	15"	29	11.5"	23	9"
AG6-1916	40.5	16"	48.5	19"	45.5-48.5	18" - 19"	42	16.5"	30.5	12"	23	9"
AG6-1918	45.5	18"	53.5	21"	45.5-48.5	18" - 19"	42	16.5"	30.5	12"	23	9"
AG6-1920	51	20"	58.5	23"	45.5-48.5	18" - 19"	42	16.5"	30.5	12"	23	9"
AG6-2016	40.5	16"	48.5	19"	48.5-51	19" - 20"	43	17"	32	12.5"	23	9"
AG6-2018	45.5	18"	53.5	21"	48.5-51	19" - 20"	43	17"	32	12.5"	23	9"
AG6-2020	51	20"	58.5	23"	48.5-51	19" - 20"	43	17"	32	12.5"	23	9"

### Section 1 | Product images

**Direct mount hardware**



**ROHO Privacy Shield**

