

# Invacare® Top End® Terminator™ Everyday Ultra Lightweight Active Performance Rigid Wheelchair

## Price List and Order Form

Effective January 1, 2006 Revised June 1, 2015

For accuracy of ordering, please fax order to Top End at 727-522-1007 or e-mail to TopEndOrders@invacare.com

For questions call 800-532-8677 or 727-522-8677

To order Literature: 800-828-6282

www.topendwheelchair.com



|                               |       |        |       |                    |       |
|-------------------------------|-------|--------|-------|--------------------|-------|
| Provider Name                 | _____ | Date   | _____ | Top End Job Number | _____ |
| Account #                     | _____ | P.O. # | _____ | Phone #            | _____ |
| Address                       | _____ | City   | _____ | State              | _____ |
| Zip Code                      | _____ |        |       |                    |       |
| Provider Primary Contact Name | _____ | Email  | _____ | Phone No           | _____ |
| Customer Name                 | _____ | Email  | _____ | Phone No           | _____ |
| Customer Address              | _____ | City   | _____ | State              | _____ |
| Zip Code                      | _____ |        |       |                    |       |

Indicate: ☐ Order or ☐ CAD or ☐ Quote Ship order to: ☐ Provider ☐ Customer Order will ship to Provider unless otherwise noted.

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

☐ RACE\* 101 CAD (see definition on page 3) ..... 125.00

PO required \_\_\_\_\_ E-mail address: \_\_\_\_\_

SEAT WIDTH: ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18"  
☐ Other \_\_\_\_\_ (20" or wider require HD Option)

SEAT DEPTH: ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18"

### FRONT SEAT-TO-FLOOR HEIGHT:

☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21"

### REAR SEAT-TO-FLOOR HEIGHT: <sup>8</sup>

☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5"

☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21"

### Knee-to-Heel Range

#### Front Seat-to-Floor

|     |                   |     |                   |
|-----|-------------------|-----|-------------------|
| 17" | 12 1/4" – 16 1/4" | 20" | 15 1/4" – 19 1/4" |
| 18" | 13 1/4" – 17 1/4" | 21" | 16 1/4" – 20 1/4" |
| 19" | 14 1/4" – 18 1/4" |     |                   |

### FRAME TYPE:

- ☐ TED Top End Terminator Everyday (K0005) ..... 2,490.00 <sup>1</sup>  
☐ TED61 Heavy Duty Upgrade (K0005) ..... 500.00 <sup>2</sup>

### FRONT FRAME STYLE: (CHOOSE ONE)

- Tapered Front End ..... Standard  
☐ Tapered Seat Width ..... No Charge <sup>3</sup>  
Taper Difference rear to front ☐ -1" ☐ -2" ☐ -3" ☐ -4" Other \_\_\_\_\_  
☐ Non Tapered V Front End ..... No Charge <sup>4</sup>

### FOOTREST OPTIONS FOR TAPERED FRONT END OR TAPERED SEAT WIDTH:

- One Piece Tubular Aluminum Footrest ..... Standard <sup>5</sup>  
☐ TED10 Plastic Footrest Cover (tubular footrest only) ..... 75.00  
☐ TED89 Carbon Fiber Footrest Cover (tubular footrest only) ..... 125.00  
☐ Adjustable Angle One Piece Footplate ..... No Charge  
☐ TED191 Adjustable Angle Footplate w/Risers ..... 95.00 <sup>5</sup>  
Footrest Height from Floor \_\_\_\_\_

### FOOTREST OPTIONS FOR NON TAPERED V FRONT END

- Adjustable Angle One Piece Footplate ..... Standard <sup>6</sup>  
☐ TED191 Adjustable Angle Footplate w/Risers ..... 95.00 <sup>5</sup>  
Footrest Height from Floor \_\_\_\_\_

### LOWER REAR FRAME


- Curved ..... Standard  
☐ TEDTi5 Square ..... 100.00

### SUSPENSION OPTIONS:

- ☐ TEDTi4 Rear Suspension Option (E1015) ..... 350.00  
<sup>7</sup>Weight of User \_\_\_\_\_

### FOOTNOTES

- <sup>1</sup> Terminator Everyday built of 6061T6 Aluminum in 1" tubing only.  
<sup>2</sup> 350 lb. weight limit for TED. Restrictions: Only available with folding back, square rear frame, single camber system (A). Solid tires, swing-away arms, one-sided forks, anodized forks and suspension not available. TED Heavy Duty will be built of 4130 Chrome-moly steel.  
<sup>3</sup> Tapered seat width tapers rear to front up to a maximum difference of 4". The front seat width measures the same as the outside of the footrest or straight so there is no tapered or V front end. Not available with Ergo seat or screw-on seat upholstery.  
<sup>4</sup> Not available w/tapered seat width or tubular footrest.  
<sup>5</sup> One piece tubular footrest is open and does not have a cover. Either a plastic or carbon fiber footrest cover is available. Tubular footrest not available with V style front end. For short leg lengths see knee-to-heel range chart on page 2 and consider option TED191.  
<sup>6</sup> Non-tapered V style front end requires adjustable angle one-piece footplate.



### ERGO SEAT OPTION

- ☐ TED104 Ergo Seat w/tapered front end ..... 250.00 <sup>9</sup>  
☐ TED110 Ergo Seat w/V style front end ..... 350.00

### FOOTNOTES

- <sup>7</sup> Front suspension forks not available at this time but can be ordered from Frog Legs and installed by provider. Not compatible with 85 or 90 degree front frame angle, seat depths less than 16" or one-sided forks. Maximum weight capacity is 250 lbs. Minimum rear seat height is 16" and for 26" or 25" knobby wheels is 17" with rear suspension option.  
<sup>8</sup> All rear seat heights are without suspension. 14.5" is lowest rear seat height for 24" wheels with Primo, Kenda or solid tires, 15" for 24" with Schwalbe and 1 3/8" tires and 25" wheels except w/knobby, 15.5" for 24" w/knobby tires and 16" for 25" knobby or 26" wheels.  
<sup>9</sup> Ergo seat option has the most benefit when seat height differentials are between 2" and 4". Not available with tapered seat width. A minimum rear seat height of 15" is required with 24" rear wheels.

# Invacare®Top End® Terminator™ Everyday


## Ultra Lightweight Active Performance Rigid Wheelchair

### REAR TO FRONT CASTER DISTANCE:

STD= Seat Depth plus 3" ☐ Other \_\_\_\_\_

### FRONT FRAME ANGLE: <sup>10</sup>

☐ 70° ☐ 75° ☐ 80° ☐ 85° ☐ 90°



### FOOT SPACE WIDTH: standard as per diagram above

Exception: Seat widths of 12", 13", and 14" will have footspace width of 8" ☐ Other \_\_\_\_\_

### BACK HEIGHT: (CHOOSE ONE) <sup>11</sup>

☐ Adjustable Height Back without Push Handles ..... Standard

☐ 9" – 11" ☐ 10" – 14" ☐ 12" – 16" ☐ 14" – 18" ☐ 16" – 20"

☐ Fixed Ht. Back without Push Handles (8" – 16") ..... No Charge

☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16"

### BACK TYPE: (FOLDING OR OPTIONAL FIXED)

**FOLDING BACK:** Fold-down back, Locking, Adjustable height/ angle ..... Standard <sup>12</sup>

#### SIDE GUARD OPTIONS FOR FOLDING BACK:

- ☐ TED18 Removable Side Guards (with hole) ..... 125.00
  - ☐ 6" tall (sm) ☐ 7.5" tall (lg)
- ☐ TED79 Removable Side Guards/Fender ..... 175.00
  - ☐ 6" tall (sm) ☐ 7.5" tall (lg)
- ☐ TED121 Fold-in Style Black Aluminum Side Guards ..... 125.00 <sup>13</sup>
  - ☐ 4.5" tall (sm) ☐ 5.75" tall (lg)
- ☐ TED142 Fold-in Style Carbon Fiber Side Guards ..... 325.00
  - ☐ 4.5" tall (sm) ☐ 5.75" tall (lg)

☐ **OPTIONAL FIXED BACK** ..... No Charge

Fixed Angle Welded at: ☐ 0° ☐ 3° ☐ 6° ☐ 9°

Backrest Crossbar Opening (3.5" std.) ..... Standard

#### SIDE GUARD OPTIONS FOR FIXED BACK:

- ☐ TED20 Fixed Aluminum Side Guards ..... 175.00 <sup>14</sup>
  - ☐ Even w/Wheels ☐ 1" Above Wheels
- ☐ TED113 Fixed Aluminum Side Guards with/Fender ..... 225.00 <sup>15</sup>
- ☐ TED83 Fixed Carbon Fiber Side Guards ..... 250.00 <sup>16</sup>
  - ☐ Even w/Wheels ☐ 1" Above Wheels
- ☐ TED119 Fixed Carbon Fiber Side Guards w/Fender ..... 300.00

### FOOTNOTES

<sup>10</sup> An 85 degree front frame angle not available with 6" front casters. 90 degree front frame angle not available with 5" or 6" front casters.

<sup>11</sup> 3/4" tube telescoping in 1" tube. Back height limited to 16" on fixed back style. Fixed height not available on folding back.

<sup>12</sup> Folding lock-down back is great for transferring the chair in/out of the car.

<sup>13</sup> Fold-in Style side guards flip back for easier transfers and are adjustable in height. Compatible with both high (push or pull to lock) or undermount (scissor) wheels locks. Not available with adjustable height push handles if back height is 10-14" or less.

<sup>14</sup> Fixed back only.

<sup>15</sup> Fixed back only. Not available w/tapered seat width.

<sup>16</sup> Fixed back only. Not available w/tapered seat width.

### BACK UPHOLSTERY AND ACCESSORIES:

#### Ballistic Nylon Adjustable Tension Upholstery with /Straps &

**Anchors** ..... Standard

- ☐ MESH Breathable Mesh Adjustable Tension Upholstery with/Straps ..... No Charge
- ☐ TED16 Fold Over Back Upholstery ..... 100.00 <sup>17</sup>
- ☐ TED30 Back Pack ..... 50.00
- ☐ TED82 Push Handles Adjustable Height Bolt-On ..... 200.00 <sup>18</sup>
- ☐ TED13 Push Handles, Non-Adjustable Height, Bolt-On ... 75.00
- ☐ TED90 Welded Integral Push Handles ..... 200.00 <sup>19</sup>
- ☐ TED97 Bolt-Through Integral Push Handles ..... 100.00 <sup>19</sup>

### FRONT WHEELS/FORK UPGRADE:

**Semi-Pneumatic** ☐ 4" ☐ 5" ☐ 6" ..... Standard <sup>20</sup>

☐ 3" Micro ..... No Charge

☐ TED165 3" Aluminum Soft Roll ..... 90.00

☐ TED65 Aluminum Soft Roll ☐ 4" ☐ 5" ☐ 6" ..... 90.00 <sup>20</sup>

☐ TED133 One Sided Forks ..... 150.00 <sup>21</sup>

| Wheel/Tire Chart – Applies to Standard & Spinergy® Wheels |     |     |     |     |
|-----------------------------------------------------------|-----|-----|-----|-----|
| Tire Choices                                              | 22" | 24" | 25" | 26" |
| Primo: 1" high pressure grey clincher                     | X   | X   | X   | X   |
| Kenda: 1" high pressure clincher                          |     | X   | X   | X   |
| Schwalbe 1" Black flat resistant                          |     | X   | X   | X   |
| 1 3/8" grey treaded                                       | X   | X   |     | X   |
| 1 3/8" w/flat free                                        | X   | X   |     |     |
| Solid 1" no-flat                                          |     | X   | X   |     |
| Knobby" 2" wide grey                                      |     | X   | X   |     |

### REAR WHEEL OPTIONS: (SEE WHEEL/TIRE CHART ABOVE)

- ☐ 22" Sun® Spoke ..... No Charge
- ☐ 24" Sun Spoke ..... No Charge
- ☐ 25" Sun Spoke ..... No Charge
- ☐ 26" Sun Spoke ..... No Charge
- ☐ TED67 22" Spinergy® Spox Wheels ..... 695.00 <sup>22</sup>
- ☐ TED67 24" Spinergy Spox Wheels ..... 695.00 <sup>22</sup>
- ☐ TED134 24" Spinergy LX Wheels \*(12 spoke) ..... 895.00
- ☐ TED67 25" Spinergy Spox Wheels ..... 695.00
- ☐ TED135 25" Spinergy LX Wheels\*(12 spoke) ..... 895.00
- ☐ TED67 26" Spinergy Spox Wheels ..... 695.00
- ☐ TED150 26" Spinergy LX Wheels\*(12 spoke) ..... 895.00

**\*Spinergy Spoke Color: (If color is not chosen, spokes will be black.)**

☐ SPCBK – Black ☐ SPCY – Yellow ☐ SPCBW - White

☐ SPCR -- Red ☐ SPCBL – Blue

### TIRE OPTIONS: (SEE WHEEL/TIRE CHART ABOVE)

**Primo** ..... Standard

☐ Pneumatic – 1 3/8" threaded ..... No Charge

☐ TED123 Schwalbe Flat resistant: black ..... 90.00

☐ TED73 Pneumatic – FF for 22 or 24 x 1 3/8" Tires ..... 125.00 <sup>23</sup>

☐ TED60 Knobby Tires (K0067) ..... 75.00

☐ TED7 Solid/Flat Free Tires ☐ Grey ☐ Black ..... 75.00 <sup>24</sup>

☐ TED101 Kenda® Tires ☐ Red ☐ Yellow ..... 75.00

### FOOTNOTES

<sup>17</sup> TER195, TER197 and UPH1/UPH2 upholstery not available with fold over back upholstery.

<sup>18</sup> Not available with Swingaway tubular padded arms or back heights less than 11" (fixed height TED13 push handles or any of the integral push handles will accommodate Swingaway arms). Push handles are 16" tall and bolted to the back post and have 10" of adjustment up or down.

<sup>19</sup> Welded integral push handles only available on Fixed Back with fixed height. Not available with foldover back upholstery.

<sup>20</sup> 6" casters not available with seat depths of 14" or less, 85 degree front frame angle, 5" and 6" casters not available with 90 degree front frame angle.

Aluminum soft roll casters are 1.4" wide with soft roll tires. For suspension fork option, see page 1. The caster housing is flared 2" on each side for 6" casters and 1.5" for 4" and 5" casters.

<sup>21</sup> Not available with 6" casters or suspension. Black anodized color only.

<sup>22</sup> Flat free inserts not available on Spinergy wheels.

<sup>23</sup> Not available with 25" wheels.

<sup>24</sup> Not available w/heavy duty option.

# Invacare®Top End® Terminator™ Everyday Ultra Lightweight Active Performance Rigid Wheelchair

## AXLES:

### Quick Release.....Standard

|                                |                         |        |
|--------------------------------|-------------------------|--------|
| <input type="checkbox"/> TED27 | Quad Release Axles..... | 125.00 |
| <input type="checkbox"/> TED76 | Titanium Axles.....     | 150.00 |

## HANDRIM OPTIONS:

### Aluminum Anodized Handrims ..... Standard

|                                 |                                            |                      |
|---------------------------------|--------------------------------------------|----------------------|
| <input type="checkbox"/> TED75  | Titanium Handrims.....                     | 600.00 <sup>25</sup> |
| <input type="checkbox"/> TED8   | Coated handrims (K0059).....               | 100.00               |
| <input type="checkbox"/> TED99  | Natural Fit® Standard Grip Handrims.....   | 295.00 25            |
| <input type="checkbox"/> TED98  | Natural Fit Super Grip Handrims.....       | 340.00 25            |
| <input type="checkbox"/> TED199 | Natural Fit LT Standard Grip Handrims..... | 325.00 21            |
| <input type="checkbox"/> TED198 | Natural Fit LT Super Grip Handrims.....    | 370.00 21            |

## MOUNTING HANDRIMS:

### Standard Mount Handrim ..... Standard

|                                  |                                |           |
|----------------------------------|--------------------------------|-----------|
| <input type="checkbox"/>         | Close Mount Handrim.....       | No Charge |
| <input type="checkbox"/> TER40   | Cut Tabs Aluminum handrim..... | 40.00     |
| <input type="checkbox"/> TERTi40 | Cut Tabs Titanium handrim..... | 80.00     |

### Overall Width by Camber Angle

| Wheel Type | 0°   | 3° | 6°  | 9°    | 12°   |
|------------|------|----|-----|-------|-------|
| Spoke      | 6 ¾" | 8" | 10" | 12 ¼" | 14 ¼" |

Add the amount above to the seat width to obtain overall width at its widest point

## CAMBER SYSTEM: (CHOOSE STYLE A OR B/ CAMBER) <sup>26</sup>

### ☐ A Camber Tube with single camber..... Standard

Camber: ☐ 0° ☐ 3° ☐ 6° ☐ 9°

|                                |                                                                                                                 |        |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------|--------|
| <input type="checkbox"/> Ted64 | Extra Camber Tube.....                                                                                          | 200.00 |
|                                | <input type="checkbox"/> 0° <input type="checkbox"/> 3° <input type="checkbox"/> 6° <input type="checkbox"/> 9° |        |

### ☐ B Adjustable Camber System (Choose Two)..... No Charge

Primary Angle..... No Charge

☐ 0° ☐ 3° ☐ 6° ☐ 9° ☐ 12°

Secondary Angle..... No Charge

☐ 0° ☐ 3° ☐ 6° ☐ 9° ☐ 12°

|                                |                                                                                                                                              |        |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|--------|
| <input type="checkbox"/> TED46 | Additional Camber Inserts (Choose Two).....                                                                                                  | 125.00 |
|                                | <input type="checkbox"/> 0° <input type="checkbox"/> 3° <input type="checkbox"/> 6° <input type="checkbox"/> 9° <input type="checkbox"/> 12° |        |

## CAMBER INSERT RELEASE TYPE:

### Nut & Bolt ..... Standard

|                                |                          |       |
|--------------------------------|--------------------------|-------|
| <input type="checkbox"/> TED48 | Quick Release Lever..... | 30.00 |
|--------------------------------|--------------------------|-------|

## FOOTREST/LEGREST ACCESSORIES:

### Calf Strap ..... Standard

|                                 |                                                  |        |
|---------------------------------|--------------------------------------------------|--------|
| <input type="checkbox"/> TED28  | Frame Protector Set (2 front, 1 rear).....       | 60.00  |
| <input type="checkbox"/> TED136 | Chair Armor Black Neoprene Frame Protectors..... | 65.00  |
| <input type="checkbox"/> TED137 | Chair Armor Blue Neoprene Frame Protectors.....  | 65.00  |
| <input type="checkbox"/> TED138 | Chair Armor Red Neoprene Frame Protectors.....   | 65.00  |
| <input type="checkbox"/> TED103 | Luggage Carriers.....                            | 100.00 |

## SEAT ACCESSORIES/ CUSHION:

|                                 |                                              |                     |
|---------------------------------|----------------------------------------------|---------------------|
| <input type="checkbox"/> TED29  | Seat Pouch.....                              | 35.00               |
| <input type="checkbox"/> TED144 | Heavy Duty Seat Pouch.....                   | 100.00              |
| <input type="checkbox"/> TED72  | Positioning Strap – Hook & Loop (K0978)..... | 40.00               |
| <input type="checkbox"/> TED74  | Positioning Strap – Auto Style (K0978).....  | 50.00               |
| <input type="checkbox"/> TED14  | 2" Soft Cushion.....                         | 75.00 <sup>27</sup> |

## SEAT RIGIDIZER/SOLID SEAT:

|                                 |                                                      |                     |
|---------------------------------|------------------------------------------------------|---------------------|
| <input type="checkbox"/> TED31  | Welded aluminum rigidizer bars w/ABS solid seat..... | 175.00              |
| <input type="checkbox"/> TED112 | Screw-On Seat Terminator Everyday.....               | 75.00 <sup>28</sup> |

## FOOTNOTES

<sup>25</sup> Not available with 22" wheels.

<sup>26</sup> Single camber (choice A) is recommended for those who do not change camber. Adjustable camber (choice B) is not available w/heavy duty option.

<sup>27</sup> Cushions are filled with 2" soft foam which offer minimum protection. If a therapeutic cushion is needed please order from Invacare Customer Service 1-800-333-6900.

<sup>28</sup> Screw-on seat upholstery not available with tapered seat width.

## ANODIZING ACCENTS: STD. FORKS AND BACK ANGLE BRACKET

### ANOBK Bold Black..... Standard

|                                 |                     |           |
|---------------------------------|---------------------|-----------|
| <input type="checkbox"/> ANOSLV | Classic Silver..... | No Charge |
| <input type="checkbox"/> TED140 | Shadow Blue.....    | 50.00     |
| <input type="checkbox"/> TED141 | Victory Red.....    | 50.00     |

### Arm Heights

(Top of Seat Rail to Top of Armrest)

|           |               |                               |
|-----------|---------------|-------------------------------|
| Swingaway | 10.5" – 12.5" | (Adjustable in 1" Increments) |
| T-Arm     | 9.5"-11.5"    | (Adjustable in 1" Increments) |

## ARM TYPE

|                                |                                                                   |                      |
|--------------------------------|-------------------------------------------------------------------|----------------------|
| <input type="checkbox"/> TED78 | Swingaway Arms (E0973): <b>not available with HD option</b> ..... | 175.00 <sup>29</sup> |
| <input type="checkbox"/> TED17 | T-Arm Adjustable Height. Desk (E0973).....                        | 175.00 <sup>30</sup> |

## WHEEL LOCKS (CHOOSE ONE)

### Hideaway Undermount Wheel Locks ..... Standard

|                                |                                       |                     |
|--------------------------------|---------------------------------------|---------------------|
| <input type="checkbox"/>       | Push or Pull To Lock Wheel Locks..... | No Charge           |
| <input type="checkbox"/> TED81 | Wheel Lock Extensions.....            | 39.00 <sup>31</sup> |
| <input type="checkbox"/> TED22 | Hill Holders (E0974).....             | 150.00              |
| <input type="checkbox"/>       | Omit Wheel Locks.....                 | No Charge           |

## MISCELLANEOUS ACCESSORIES

|                                 |                                                                                        |                     |
|---------------------------------|----------------------------------------------------------------------------------------|---------------------|
| <input type="checkbox"/> TED21  | Spoke Guard <input type="checkbox"/> Clear <input type="checkbox"/> Black (K0065)..... | 75.00               |
| <input type="checkbox"/> TED23  | Rear Anti Tippers (E0971).....                                                         | 100.00              |
| <input type="checkbox"/> TED26  | Crutch Holder.....                                                                     | 125.00              |
| <input type="checkbox"/> TED52  | Anti Tipper w/3" Swivel Wheel.....                                                     | 150.00              |
| <input type="checkbox"/> TED94  | Alignment Gauge.....                                                                   | 75.00               |
| <input type="checkbox"/> TED145 | Water bottle and cage.....                                                             | 50.00               |
| <input type="checkbox"/> RACE22 | Top End T-Shirt <input type="checkbox"/> L <input type="checkbox"/> XL.....            | 10.00 <sup>32</sup> |

## FRAME COLOR AND UPHOLSTERY

### UPHOLSTERY: Black..... Standard

|                               |                                |                      |
|-------------------------------|--------------------------------|----------------------|
| <input type="checkbox"/> UPH1 | Woodland Camouflage.....       | 100.00 <sup>33</sup> |
| <input type="checkbox"/> UPH2 | Digital Desert Camouflage..... | 100.00 33            |

## FRAME COLOR:

## CONTEMPORARY COLORS

|                                          |                                        |                                     |
|------------------------------------------|----------------------------------------|-------------------------------------|
| <input type="checkbox"/> Wet Black       | <input type="checkbox"/> Black Prism   | <input type="checkbox"/> White Snow |
| <input type="checkbox"/> Silver Metallic | <input type="checkbox"/> Stealth Black |                                     |

## WAY COOL COLORS

|                                             |                                        |                                            |
|---------------------------------------------|----------------------------------------|--------------------------------------------|
| <input type="checkbox"/> Electric Blue      | <input type="checkbox"/> Midnight Blue | <input type="checkbox"/> Grasshopper Green |
| <input type="checkbox"/> Electric Blueberry |                                        |                                            |

## TOO HOT COLORS

|                                        |                                           |                                         |
|----------------------------------------|-------------------------------------------|-----------------------------------------|
| <input type="checkbox"/> Tangerine Red | <input type="checkbox"/> Electric Red     | <input type="checkbox"/> Grape Madness  |
| <input type="checkbox"/> Safety Yellow | <input type="checkbox"/> Sunkissed Orange | <input type="checkbox"/> Bubblegum Pink |

**NOTE: CAD drawings** are computer aided design that translate the measurements to a 3D preview. Our experienced draftsman completes each CAD request and will send it back for final approval. The charge for RACE101 is \$125 with standard order form selections. A separate purchase order is required to begin a CAD and another purchase order is needed to order the product.

## FOOTNOTES

<sup>29</sup> Swingaway tubular padded arms interfere with adj. height push handles. Fixed height, bolt-on or integral push handles are compatible with Swingaway Arms.

<sup>30</sup> Rigid side guards are included with the T-Arms.

<sup>31</sup> Push or Pull to Lock Wheel Locks only.

<sup>32</sup> Dealer discount does not apply.

<sup>33</sup> Camouflage upholstery not available for screw-on seat, mesh or fold over back upholstery. Camouflage upholstery is ballistic nylon adjustable tension with anchors and straps for the back upholstery. If cushion is ordered, it will also be camouflage.

# Top End® Terminator™ Everyday Ultra Lightweight Active Performance Rigid Wheelchair

| *Overall Chair Weight: Width & Depth Matrix |       |       |       |       |       |       |       |       |       |       |
|---------------------------------------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| Width (across)<br>Depth (down)              | 12"   | 12"   | 13"   | 14"   | 15"   | 16"   | 17"   | 18"   | 19"   | 20"   |
| 12"                                         | 18.43 | 18.53 | 18.63 | 18.73 | 18.83 | 18.93 | 19.03 | 19.13 | 19.23 | 19.33 |
| 13"                                         | 18.53 | 18.63 | 18.73 | 18.83 | 18.93 | 19.03 | 19.13 | 19.23 | 19.33 | 19.43 |
| 14"                                         | 18.63 | 18.73 | 18.83 | 18.93 | 19.03 | 19.13 | 19.23 | 19.33 | 19.43 | 19.53 |
| 15"                                         | 18.73 | 18.83 | 18.93 | 19.03 | 19.13 | 19.23 | 19.33 | 19.43 | 19.53 | 19.63 |
| 16"                                         | 18.83 | 18.93 | 19.03 | 19.13 | 19.23 | 19.33 | 19.43 | 19.53 | 19.63 | 19.73 |
| 17"                                         | 18.93 | 19.03 | 19.13 | 19.23 | 19.33 | 19.43 | 19.53 | 19.63 | 19.73 | 19.83 |
| 18"                                         | 19.03 | 19.13 | 19.23 | 19.33 | 19.43 | 19.53 | 19.63 | 19.73 | 19.83 | 19.93 |
| 19"                                         | 19.13 | 19.23 | 19.33 | 19.43 | 19.53 | 19.63 | 19.73 | 19.83 | 19.93 | 20.03 |
| 20"                                         | 19.23 | 19.33 | 19.43 | 19.53 | 19.63 | 19.73 | 19.83 | 19.93 | 20.03 |       |

## Section 1

| Standard Chair Includes:           |
|------------------------------------|
| Aluminum Frame                     |
| Tubular footrest w/o plate         |
| 4" Semi-pneumatic casters/forks    |
| 3° Camber                          |
| 24" High performance wheels        |
| 24" High pressure tires            |
| Medium height folding back         |
| Adjustable tension seating sling   |
| Adjustable tension back upholstery |

## Top End Terminator Component Weights – Section 2

|                                         |       |
|-----------------------------------------|-------|
| <b>Seating</b>                          |       |
| Omit adjustable tension seat upholstery | -.85  |
| No cushion                              | 0     |
| 2" cushion                              | +1.10 |
| Screw-on seat upholstery                | +.05  |
| <b>Back Type</b>                        |       |
| Folding                                 | 0     |
| Fixed                                   | -.80  |
| <b>Push Handles</b>                     |       |
| Adjustable height bolt-on               | +1.60 |
| Fixed height bolt-on                    | +1.70 |
| Integral: welded                        | +.75  |
| Integral: bolt-through                  | +.75  |
| <b>Back Upholstery</b>                  |       |
| Adjustable tension w/straps             | 0     |
| Fold over                               | +.09  |
| Omit back upholstery                    | -.97  |
| <b>Armrests</b>                         |       |
| Swing away padded tubular               | +2.64 |
| T-Arm, Adj. Height Desk                 | +6.25 |
| <b>Footrest</b>                         |       |
| Tubular                                 | 0     |
| Adjustable angle one piece footplate    | +.60  |
| Plastic footrest cover                  | +.35  |
| Carbon fiber footrest cover             | -.30  |
| <b>Front Casters (pair)</b>             |       |
| 3" Micro                                | -.13  |
| 3" Softroll                             | +.19  |
| 4" Semi pneumatic                       | 0     |
| 5" Semi pneumatic                       | +.31  |
| 4" Aluminum softroll                    | +.44  |
| 5" Aluminum softroll                    | +.92  |
| One-sided forks                         | -.05  |
| <b>Wheel Locks</b>                      |       |
| Push or pull-to-lock                    | +1.00 |
| Undermount scissor                      | +1.25 |
| Wheel lock extensions                   | +.40  |
| Hill holders/grade aids                 | +.51  |
| <b>Rear Suspension</b>                  | +.40  |

|                                    |       |
|------------------------------------|-------|
| <b>Sideguards</b>                  |       |
| Fold-in Style carbon fiber         | +.95  |
| Fold-in Style black aluminum       | +1.45 |
| Folding                            | +1.45 |
| Removable                          | +1.45 |
| Removable side guards with fender  | +1.94 |
| Fixed aluminum                     | +1.45 |
| Fixed carbon fiber                 | +1.00 |
| <b>Rear Tires (pair)</b>           |       |
| High pressure, Primo or Kenda      | 0     |
| Schwalbe                           | +.75  |
| Pneumatic 1 3/8"                   | +1.40 |
| Pneumatic Flat Free                | +1.45 |
| Knobby                             | +3.10 |
| Solid 1" flat free                 | +2.15 |
| <b>Rear Wheels (pair)</b>          |       |
| 22" spoke                          | -.60  |
| 24" Sun high performance spoke     | 0     |
| 25" Sun high performance spoke     | +.30  |
| 26" Sun high performance spoke     | +.65  |
| 24" Spinergy Spox                  | -.80  |
| 24" Spinergy LX                    | -1.20 |
| 25" Spinergy Spox                  | -.60  |
| 25" Spinergy LX                    | -1.00 |
| 26" Spinergy Spox                  | -.25  |
| <b>Handrim</b>                     |       |
| Aluminum                           | 0     |
| Titanium                           | +.20  |
| Plastic coated                     | +1.40 |
| Natural fit standard or super grip | +1.70 |
| Natural fit LT                     | +1.25 |
| <b>Axles</b>                       |       |
| Standard quick release             | 0     |
| Titanium quick release             | -.17  |
| Quad quick release                 | +.22  |
| <b>Accessories</b>                 |       |
| Frame protectors/Chair armor       | +.10  |
| Luggage carriers                   | +.75  |
| Calf strap                         | +.10  |
| Seat pouch                         | +.35  |
| Positioning strap – hoop and loop  | +.10  |
| Positioning strap – auto style     | +.50  |
| Spoke guards – black               | +1.40 |
| Spoke guards – clear               | +1.60 |
| Anti-tip bars                      | +1.63 |
| Crutch holder                      | +1.45 |

\*Weights are approximate.

Section 1 \_\_\_\_\_ + Section 2 \_\_\_\_\_ = Total Weight \_\_\_\_\_



### Invacare Corporation

**TOP END**  
Sports and Recreation Products  
4501 63rd Circle North  
Pinellas Park FL USA  
33781-5914

www.invacare.com

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Sun® is a registered trademark of Sun Components Inc. Natural Fit is a registered trademark of Three Rivers Holding, LLC. Spinergy™ is a registered trademark of Spinergy Corp. Kenda is a registered trademark of Kenda Rubber Industrial Co. Ltd. Schwalbe is a registered trademark of RALF BOHLE GmbH in Reichshof Germany. Specifications and prices are subject to change without notification.  
©2015 Invacare Corporation