

# Invacare eForms

For additional information, visit [www.invacare.com/eForms](http://www.invacare.com/eForms)

[www.invacare.com](http://www.invacare.com)

Customer Service: 800-333-6900

Fax: 800-678-4682

To order Literature: 800-828-6282

## Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

**Please consider the environment before printing this page with your order form.**  
*For more information on features and benefits of electronic order forms, please visit [www.invacare.com/eForms](http://www.invacare.com/eForms)*



**TRANSACTION TYPE:**

- Order       Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_  
 Company Name \_\_\_\_\_ Phone # \_\_\_\_\_  
 Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_  
 Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

**Special Note**

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

**FRAME SELECTIONS**

**FRONT FRAME**

- MVPS MVP Swingaway Front Frame ..... \$2,043.00 K0005  
 MVPF60 MVP 60° Fixed Front Frame ..... \$2,060.00 K0005  
 MVPF80 MVP 80° Fixed Front Frame ..... \$2,060.00 K0005

**REAR FRAME**

- CURVRF Curved Rear Frame ..... No Charge  
 SQURRF Square Rear Frame ..... No Charge  
 EXTRF Extended Rear Frame (with Angle Adj. Back) ..... No Charge  
 EXTRF offers angle adjustable back canes standard, but a back cane type must be selected. CURVRF and SQURRF do not have standard angle adjustable back cane hardware, so if angle adjustable back canes are required, one must select either AAB1216 or AAB1720.

**CROSSBRACE TYPE**       Clear

- CONV Conventional Crossbrace ..... Standard  
 LSR Less Seat Rails ..... No Charge  
 Sling Seat not available with LSR. LSR option for depths greater than 18" is a special.

**FRAME MODIFICATIONS**       Clear

- HDPG Heavy-Duty Package ..... \$873.00 K0108  
 HDPG not available with EXTRF or LSR. HDPG available on 18" - 22" widths and 16" - 22" depths.  
 1SHFRAME 80° 1" Shorter Than Seat Depth ..... No Charge  
 Available with MVPF80 only.

**Heavy Duty Includes:**

- Weight Limitation 350 lbs.
- Reinforced seat and back upholstery
- Reinforced frame
- Triple-construction steel crossbraces

**SEAT WIDTH**

- W15 15" Wide ..... No Charge  
 W16 16" Wide ..... No Charge  
 W17 17" Wide ..... No Charge  
 W18 18" Wide ..... No Charge  
 W19 19" Wide ..... No Charge  
 W20 20" Wide ..... \$300.00 E2201  
 W22 22" Wide (Requires HDPG) ..... \$300.00 E2201

**SEAT DEPTH RANGE (FOR EXTRFW/LSR)**       Clear

- DR0915 9" - 15" Depth Range ..... No Charge  
 DR1016 10" - 16" Depth Range ..... No Charge  
 DR1117 11" - 17" Depth Range ..... No Charge  
 DR1218 12" - 18" Depth Range ..... No Charge  
 DR1319 13" - 19" Depth Range ..... No Charge

Width and Depth Matrix										
Depth Width	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
15"	X	X	X	X	X	X	S	S	S	S
16"	X	X	X	X	X	X	S	S	S	S
17"	X	X	X	X	X	X	S	S	S	S
18"	X	X	X	X*	X*	X*	X*	X*	S	S
19"	S	S	S	X*	X*	X*	X*	X*	S	S
20"	S	S	S	X*	X*	X*	X*	X*	X*	X*
22"	S	S	S	S	S	S	S	S	S	S
X - Sizes available without HDPG						S - Sizes available through				
* - Sizes available with HDPG						Specials				

**NOTE: All specifications and dimensions are approximate.**  
 For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

**FRAME AND SEAT**

**SEAT DEPTH (FOR CURVE, SCURRF AND EXTR WITH CONV)  Clear**

- D13 13" Seat Depth ..... No Charge
- D14 14" Seat Depth ..... No Charge
- D15 15" Seat Depth ..... No Charge
- D16 16" Seat Depth ..... No Charge
- D17 17" Seat Depth ..... No Charge
- D18 18" Seat Depth ..... No Charge
- D19 19" Seat Depth ..... No Charge
- D20 20" Seat Depth ..... \$300.00 E2203
- D21 21" Seat Depth ..... \$300.00 E2203
- D22 22" Seat Depth ..... \$300.00 E2204

**SEAT DEPTH SETTING (FOR EXTRFW/LSR)  Clear**

- SDS09 9" Seat Depth Setting ..... No Charge
- SDS10 10" Seat Depth Setting ..... No Charge
- SDS11 11" Seat Depth Setting ..... No Charge
- SDS12 12" Seat Depth Setting ..... No Charge
- SDS13 13" Seat Depth Setting ..... No Charge
- SDS14 14" Seat Depth Setting ..... No Charge
- SDS15 15" Seat Depth Setting ..... No Charge
- SDS16 16" Seat Depth Setting ..... No Charge
- SDS17 17" Seat Depth Setting ..... No Charge
- SDS18 18" Seat Depth Setting ..... No Charge

**FRONT SEAT-TO-FLOOR**

**LSR option reduces seat-to-floor heights shown by 1 1/2". Refer to Seat-to-Floor Height charts for correct seat-to-floor with Swing Away or Fixed Front Frames.**

- FSTF1475 14 3/4" Front STF Height ..... \$100.00 K0056
- FSTF1525 15 1/4" Front STF Height ..... \$100.00 K0056
- FSTF155 15 1/2" Front STF Height ..... \$100.00 K0056
- FSTF1575 15 3/4" Front STF Height ..... \$100.00 K0056
- FSTF16 16.00" Front STF Height ..... No Charge
- FSTF165 16 1/2" Front STF Height ..... \$100.00 K0056
- FSTF1675 16 3/4" Front STF Height ..... \$100.00 K0056
- FSTF17 17.00" Front STF Height ..... No Charge
- FSTF1725 17 1/4" Front STF Height ..... No Charge
- FSTF18 18.00" Front STF Height ..... No Charge
- FSTF1825 18 1/4" Front STF Height ..... No Charge
- FSTF185 18 1/2" Front STF Height ..... No Charge
- FSTF1875 18 3/4" Front STF Height ..... No Charge
- FSTF19 19" Front STF Height ..... No Charge
- FSTF1925 19 1/4" Front STF Height ..... No Charge
- FSTF195 19 1/2" Front STF Height ..... No Charge
- FSTF1975 19 3/4" Front STF Height ..... No Charge
- FSTF20 20" Front STF Height ..... No Charge
- FSTF2025 20 1/4" Front STF Height ..... No Charge
- FSTF205 20 1/2" Front STF Height ..... No Charge
- FSTF21 21" Front STF Height ..... No Charge
- FSTF2125 21 1/4" Front STF Height ..... \$100.00 K0056

Front Seat -to-Floor Heights							
Fixed Front Frame (Non - Swingaway)							
Fork	Spacer	AT900 3"	51DC 5"	61DC 615DC 6"	1259 6"	81DC 8"	815DC 8"
3"		14 3/4 15 3/4	15 1/2 16 1/2				
3"	1"	15 1/2 16 1/2	16 1/2 17 1/4				
3"	2"	16 3/4 17 3/4	17 1/2 18 1/2				
5"		16 1/4 17 1/4	17 18	17 1/4 18 3/4			
5"	1"	17 18	18 19	18 1/2 19 1/2			
6"		16 3/4 17 3/4	16 16 3/4 17 17 1/2 17 3/4 18 1/2	16 1/2 17 1/4 17 1/2 18 1/4 19 1/4			
8"					17 1/2 18 1/2 18 3/4 19 3/4	18 1/4 19 1/4 20 1/4	18 1/2 19 1/2 20 1/2

Front Seat - to - Floor Height							
Swingaway Front Frame							
Fork	Spacer	AT900 3"	51DC 5"	61DC 615DC 6"	1259 6"	81DC 8"	815DC 8"
3"		15 1/4 16 1/2	16 17 1/4				
3"	1"	16 17 1/4	16 3/4 18				
3"	2"	17 1/4 18 1/2	18 19 1/4				
5"		16 3/4 18	17 1/2 18 3/4	18 19 1/4			
5"	1"	17 1/2 18 3/4	18 1/2 19 3/4	19 20 1/4			
6"		17 1/4 18 1/2	16 1/2 17 1/4 17 3/4 18.00 18 1/2 19 1/4	17 17 3/4 18 1/4 18 3/4 19 20			
8"					18 19 19 1/4 20 1/4	18 3/4 19 3/4 20 21	19 20 20 1/4 21 1/4

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

**ARM AND BACK SELECTIONS**

**ARM TYPE (PAIR)**

- AT99 Swingaway Padded Arms ..... No Charge  
**Not available with EXTRF.**
- AT99N Cantilever Non Locking Bolt On Arms ..... No Charge  
**Not available with Back Type AAB1216 or AAB1720.**  
**Not available with HDPG.**
- AT55D Cantilever Locking Desk Lgth. Bolt On Arms ..... \$195.00  
**Not available with Back Type AAB1216 or AAB1720.**  
**Not available with HDPG.**  
**Only available with FBH16, FBH18, and FBH20.**
- AT55F Cantilever Locking Full Lgth. Bolt On Arms ..... \$195.00  
**Not available with Back Type AAB1216 or AAB1720.**  
**Not available with HDPG.**  
**Only available with FBH16, FBH18, and FBH20.**
- DESKADJCV Conventional Adjustable Ht Desk Lgth Arms ... \$195.00 E0973  
**Not available with LSR (Less Seat Rail).**
- DESKADJSS Spacesaver Adjustable Ht Desk Lgth Arms ..... \$195.00 E0973  
**Not available with EXTRF.**  
**Not available with Back Type AAB1216, or AAB1720.**  
**Not available with LSR (Less Seat Rail).**
- FULLADJCV Conventional Adjustable Ht Full Lgth Arms .... \$195.00 E0973  
**Not available with LSR (Less Seat Rail).**
- FULLADJSS Spacesaver Adjustable Ht Full Lgth Arms ..... \$195.00 E0973  
**Not available with EXTRF.**  
**Not available with Back Type AAB1216, or AAB1720.**  
**Not available with LSR (Less Seat Rail).**
- TAD T - Arm Adjustable Height, Desk Length ..... \$225.00 E0973  
**Bolt-On Arm Attachment.**
- TAF T - Arm Adjustable Height, Full Length ..... \$225.00 E0973  
**Bolt-On Arm Attachment.**
- LESSARMS Less Arms.....Less \$24.00

**T-ARM HEIGHT**

- L6-10 Low Range 6" - 10" (T-Arm Height) ..... No Charge
- H8-14 High Range 8" - 14" (T-Arm Height) ..... No Charge

**BACK TYPE**

- A1418 10° Bendback with Push Handles (14" - 18") . No Charge  
**Not available with AT99N, AT55D, or AT55F arms.**
- 1014 Sportster Back without Push Handles (10" - 14") ... No Charge  
**Not available with AT99N, AT55D, or AT55F arms.**
- A1216 Straight Back with Push Handles (12" - 16") . No Charge  
**Not available with AT99N, AT55D, or AT55F arms.**
- S1418 Straight Back with Push Handles (14" - 18") . No Charge  
**Not available with AT99N, AT55D, or AT55F arms.**
- FBH12 Fixed Back Height 12" ..... No Charge  
**Not available with HDPG.**  
**Not available with AT99N, AT55D, or AT55F arms.**  
**Not available with CURVRF, SQURRF, or with 22" widths.**
- FBH14 Fixed Back Height 14" ..... No Charge  
**Not available with HDPG.**  
**Not available with AT99N, AT55D, or AT55F arms.**  
**Not available with CURVRF, SQURRF, or with 22" widths.**
- FBH16 Fixed Back Height 16" ..... No Charge  
**Not available with AT99N, AT55D, or AT55F arms.**  
**Not available with CURVRF, SQURRF, or with 22" widths.**
- FBH18 Fixed Back Height 18" ..... No Charge  
**Not available with AT99N, AT55D, or AT55F arms.**  
**Not available with CURVRF, SQURRF, or with 22" widths.**
- FBH20 Fixed Back Height 20" ..... No Charge  
**Not available with HDPG.**  
**Not available with AT99N, AT55D, or AT55F arms.**  
**Not available with CURVRF, SQURRF, or with 22" widths.**  
**Not available with 12" or 13" wide ADJTENS upholstery.**
- AAB1216 Adjustable Angle Back (12" - 16") ..... \$165.00 K0108  
**Not available with HDPG.**  
**Available with SQURRF only.**
- AAB1720 Adjustable Angle Back (17" - 20") ..... \$165.00 K0108  
**Not available with HDPG.**  
**Available with SQURRF only.**

Rear Seat - to - Floor Heights		
Rear Wheel	Curved Rear Frame	Square Rear Frame Extended Rear Frame
16COM	11 1/2 - 14 1/4	12 1/4 - 15 1/4
18COM	12 1/2 - 15 1/4	13 1/4 - 16 1/4
20COM	13 3/4 - 16 1/4	14 1/2 - 17
22COM	14 3/4 - 17 1/4	15 1/2 - 18
24COM7	15 3/4 - 18 1/4	16 1/2 - 19
20SP	13 3/4 - 16 1/4	14 1/2 - 17
22SP	14 3/4 - 17 1/4	15.55 - 18
24SP	15 3/4 - 18.25	16 1/2 - 19
26SP	16 3/4 - 19 1/4	17 1/2 - 20
SPIN24	15 3/4 - 18 1/4	16 1/2 - 19
SPIN25	15 3/4 - 19 1/4	16 1/2 - 20
SPIN26	16 3/4 - 19 3/4	17 1/2 - 20
SPLX24	15 3/4 - 18 1/4	16 1/2 - 19
SPLX25	15 3/4 - 19 1/4	16 1/2 - 20
FN1624	15 3/4 - 18 1/4	16 1/2 - 19
FN1625	15 3/4 - 19 1/4	16 1/2 - 20
FN1626	16 3/4 - 19 1/4	17 1/2 - 20

Arm Heights	
AT99 AT99N Range adjusts in 1" increments	8" - 10" 5" - 10"
AT55D AT55F Range adjusts in 1" increments	5" - 13" 5" - 13"
DESKADJCV, DESKADJSS FULLADJCV, FULLADJSS Range adjusts in 1" increments	10" - 14" 10" - 14"
TAD Range adjusts in 1/2" increments TAF Range adjusts in 1" increments	Low Range 6" - 10" High Range 8" - 14"
All measurements from top of seat rail to top of armrest. On AT99, AT55D and AT55F, add 1" for upholstered arm pads. For less seat rail chairs, add 1-1/2" to above measurements (top of upper frame rail to top of armrest).	

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

**ADJUSTABLE BACK HEIGHT**

- ABH10 10" High Adjustable Back ..... No Charge
- ABH11 11" High Adjustable Back ..... No Charge
- ABH12 12" High Adjustable Back ..... No Charge
- ABH13 13" High Adjustable Back ..... No Charge
- ABH14 14" High Adjustable Back ..... No Charge
- ABH15 15" High Adjustable Back ..... No Charge
- ABH16 16" High Adjustable Back ..... No Charge
- ABH17 17" High Adjustable Back ..... No Charge
- ABH18 18" High Adjustable Back ..... No Charge
- ABH19 19" High Adjustable Back ..... No Charge
- ABH20 20" High Adjustable Back ..... No Charge

**BACK UPHOLSTERY TYPE**

- Clear
- SLIPBACK Slip On Back Upholstery ..... Standard  
**Not available with AT99N, AT55D, or AT55F arms.**
- ADJTENS Adjustable Tension Back Upholstery ..... \$125.00 K0108  
**Not available with AT99N, AT55D, or AT55F arms.**
- L8601 Omit Back Upholstery ..... Less \$10.00

**CASTER SELECTION**

**CASTER OPTIONS (PAIR)**

- AT900 3" Roller Blade Casters ..... \$38.00  
**Not available with HDPG**
- 5ALUM 5" Aluminum Urethane Casters ..... \$92.00
- 51DC 5" Urethane Casters ..... No Charge  
**Not available with HDPG**
- 61DC 6" Urethane Casters ..... No Charge  
**Not available with HDPG**
- 615DC 6" x 1.5" Urethane Casters ..... No Charge  
**Not available with HDPG**
- 1259 6" x 2" Semi Pneumatic Casters ..... \$60.00
- 81DC 8" x 1" Urethane Casters ..... No Charge
- 815DC 8" x 1.5" Urethane Casters ..... No Charge

**FORK TYPE (PAIR)**

- Clear
- A1103 3" Fork (Single Hole Fork) ..... No Charge
- A1105 5" Fork (3 Hole Fork) ..... No Charge
- A1106 6" Fork (3 Hole Fork) ..... No Charge
- A1108 8" Fork (2 Hole Fork) ..... No Charge

**SPACER OPTION**

- ASP00 No Spacer ..... No Charge
- ASP01 1" Spacer ..... No Charge
- ASP02 2" Spacer ..... No Charge

**REAR WHEELS**

**REAR WHEEL TYPE (PAIR)**

- 16COM 16" Composite Wheels ..... No Charge
- 18COM 18" Composite wheels ..... No Charge
- 20COM 20" Composite Wheels ..... No Charge
- 22COM 22" Composite Wheels ..... No Charge
- 24COM7 24" Composite 7 Spoke Wheels ..... No Charge
- 20SP 20" Spoke Wheels ..... No Charge  
**Not available wit HDPG or 19" Seat-to-floor**
- 22SP 22" Spoke Wheels ..... No Charge  
**Not available wit HDPG or 19" Seat-to-floor**
- 24SP 24" Spoke Wheels ..... No Charge  
**Not available wit HDPG or 19" Seat-to-floor**
- 26SP 26" Spoke Wheels ..... \$48.00  
**Not available wit HDPG or 19" Seat-to-floor**
- FN1624 24" Sun® Fusion 16 Spoke Wheel ..... \$325.00
- FN1625 25" Sun® Fusion 16 Spoke Wheel ..... \$325.00
- FN1626 26" Sun® Fusion 16 Spoke Wheel ..... \$325.00
- SPLX24 24" Spinegy® LX Wheels ..... \$895.00
- SPLX25 25" Spinegy® LX Wheels ..... \$895.00
- SPIN24 24" Spinegy® Spoke Wheels ..... \$695.00  
**Not available wit HDPG or 19" Seat-to-floor. Only available with ALUMSP**
- SPIN25 25" Spinegy® Spoke Wheels ..... \$695.00  
**Only available with ALUMSP.**
- SPIN26 26" Spinegy® Spoke Wheels ..... \$695.00  
**Only available with ALUMSP.**
- 24HPS 24" High Performance Spoke Wheels ..... \$200.00
- 26HPS 26" High Performance Spoke Wheels ..... \$137.00
- 3SCF 559 Three Spoke Composite Wheel ..... \$895.00
- L8701 Omit Rear Wheels ..... Less \$50.00

**SPINERGY SPOKE COLOR**

- Clear
- SPCBK Black ..... Standard  
**If color is not chosen, spokes will be black**
- SPCBW White ..... No Charge
- SPCY Yellow ..... No Charge
- SPCR Red ..... No Charge
- SPCBL Blue ..... No Charge

**TIRE OPTIONS (PAIR)**

- Clear
- PNEU Pneumatic Tires ..... \$60.00
- PNEUFF Flat Free Inserts ..... \$65.00 E2213  
**Not available with 26SP, SPIN24, SPIN25, and SPIN26.**
- URETH Urethane (Treaded) Tires ..... No Charge  
**Not available with 20SP, 22SP, 24SP, 26SP, SPIN24, SPIN25, and SPIN26.**
- PRIMO Primo Tires ..... \$60.00  
**Only available with 26SP, SPIN24, SPIN25, and SPIN26.**

**HANDRIM OPTIONS**

- Clear
- ALUMPC Aluminum Anodized Handrims ..... No Charge  
**Not available with 16COM, 18COM, SPIN24, SPIN25, SPIN26, or 3SCF.**
- ALUMSP Aluminum Anodized Spinegy Handrims ..... No Charge
- A1512 Black Plastic Coated Handrims ..... \$72.00
- L8702 Omit Handrims ..... Less \$15.00
- NFSDLT Natural-Fit® LT Handrims Standard Grip ..... \$325.00  
**Only available with 24COM7, 24SP, 26SP, SPIN24, SPIN25, SPIN26, 24HPS, 26HPS, SPLX24, SPLX25, FN1624, FN1615, FN1626, 3SCF.**
- NFSPLT Natural-Fit® LT Handrims Super Grip ..... \$370.00  
**Only available with 24COM7, 24SP, 26SP, SPIN24, SPIN25, SPIN26, 24HPS, 26HPS, SPLX24, SPLX25, FN1624, FN1615, FN1626, 3SCF.**
- NFITSD Natural-Fit® Handrims Standard Grip ..... \$295.00  
**Only available with 24COM7, 24SP, 26SP, SPIN24, SPIN25, SPIN26, 24HPS, 26HPS, SPLX24, SPLX25, FN1624, FN1615, FN1626.**
- NFITSP Natural-Fit® Handrims Super Grip ..... \$340.00  
**Only available with 24COM7, 24SP, 26SP, SPIN24, SPIN25, SPIN26, 24HPS, 26HPS, SPLX24, SPLX25, FN1624, FN1615, FN1626.**

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

**CAMBER ANGLE**

- C00 0° Camber ..... No Charge
- C03 3° Camber ..... No Charge
- C06 6° Camber ..... No Charge
- C09 9° Camber ..... No Charge

**MOUNTING HANDRIMS (SPOKE WHEEL ONLY)**  Clear

- IN Mounting Inside Handrim ..... No Charge
- OUT Mounting Outside Handrim ..... No Charge

**WHEEL LOCKS (PAIR)**

- TSAVER Thumb Saver Wheel Locks ..... No Charge
- A1598 Hill Holder Clamp On Wheel Locks ..... \$170.00 E0974
- A1567 Pull To Lock Clamp On Wheel Locks ..... No Charge
- 9691 Push To Lock Wheel Locks ..... No Charge

**AXLE OPTIONS**

**AXLE PLATE (PAIR)**

- STDPT Standard Axle Plates ..... No Charge
- EXTPT Extended/Amputee Axle Plates ..... \$72.00 E0959  
**Not available with FULLADJCV or DESKADJCV.**

**AXLES**  Clear

- A1577 Quick Release Axles ..... Standard  
**Not available with HDPG.**
- PAXLE Permanent Axles ..... No Charge
- A1573 Quad Release Axles ..... \$89.00  
**Not available with HDPG.**  
**Not available with 16COM, 18COM, SPIN24, SPIN25, and SPIN26.**

**POSITIONING STRAPS**

**SEAT POSITIONING STRAPS (EACH)**  Clear

- A1515 Push Button Seat Positioning Strap ..... \$50.00 E0978
- A1321 Hook and Loop Seat Positioning Strap ..... \$40.00 E0978

**CHEST POSITIONING STRAPS (EACH)**  Clear

- 9515 Auto Style Chest Positioning Strap ..... \$50.00 E0960
- 7321 Hook and Loop Chest Positioning Strap ..... \$40.00 E0960

**UPHOLSTERY**

**BACK UPHOLSTERY COLOR**  Clear

**Monogramming and emblems available. Please check with Customer Service**

- U240 Black Nylon Upholstery ..... No Charge
- 24P Wet Black ..... No Charge
- 29P White Snow ..... No Charge
- 60P Silver Metallic ..... No Charge
- 115P Black Prism ..... No Charge
- 105P Midnight Blue ..... No Charge
- 116P Sunkissed Orange ..... No Charge
- 30P Sunny Yellow ..... No Charge
- 61P Electric Red ..... No Charge
- 119P Bubblegum Pink ..... No Charge
- 120P Grape Madness ..... No Charge
- 117P Tangerine Red ..... No Charge
- 62P Electric Blue ..... No Charge
- 127P Grasshopper Green ..... No Charge
- Touch Up Paint ..... \$20.00

Touch Up Color

**Touch up paint is not available in all colors.**

Discontinued frame colors are available through our specials department.  
C-070109030544

**ACCESSORIES**

**BACK ACCESSORIES**

- BP02 Back Pack ..... \$50.00
- 1483 Non Adjustable Rigidizing Handle Bar ..... \$85.00
- AASH Adj. Angle Stroller Handles (Outside Mount) ..... \$235.00  
**Not available with AAB1216 or AAB1720. AASH release buttons are mounted on the outside and AASIN release buttons are mounted on the inside.**
- AASHIN Adj. Angle Stroller Handles (Inside Mount) ..... \$235.00  
**Not available with AAB1216 or AAB1720. AASH release buttons are mounted on the outside and AASIN release buttons are mounted on the inside.**

**SEAT ACCESSORIES**

- SP-01 Seat Pouch ..... \$35.00
- 9662L Drop Hook Solid Folding Seat ..... \$126.00 K0108
- L8651 Omit Seat Upholstery ..... Less \$10.00

**ARM ACCESSORIES (PAIR)**

**Only available with AT99 or LESS ARMS.**

- A1571 Fabric Clothing Protectors ..... \$46.00 K0108
- 1671-SM Rigid Side Guard, Small ..... \$122.00 K0108
- 1671-LG Rigid Side Guard, Large ..... \$122.00 K0108
- 1671-SMH Quad Rigid Side Guard, Small ..... \$122.00 K0108
- 1671-LGH Quad Rigid Side Guard, Large ..... \$122.00 K0108

**MISCELLANEOUS ACCESSORIES**

- A1209 Tire Pump ..... \$35.00
- A1307 Spoke Protectors (For 20", 22", 24", 26") (pair) ..... \$90.00 K0065
- A1541 Caster Swivel Locks ..... \$80.00
- A1576 Wheel Lock Extension (6") (pair) ..... \$39.00 E0961
- A1580 Rear Anti Tippers (pair) ..... \$60.00 E0971
- A1569 Crutch and Cane Carrier ..... \$46.00 E2207

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

**FRONT RIGGING OPTIONS**

**FRONT RIGGING (PAIR)**

Clear

For MVPS front frame only.

- SWING1319 60° Swingaway Footrests ..... No Charge
- SWING914 60° Swingaway Footrests MFXs ..... No Charge
- 70NTAPER 70° Swingaway Footrests ..... No Charge
- 70MFX 70° Swingaway Footrests MFXs ..... No Charge  
**Standard with AT5543 adjustable angle footplate.**
- 70STAPER 70° Tapered Footrests ..... No Charge  
**Must have starting seat-to-floor height of 17" or above.**
- 90SWING 90° Swingaway Footrests ..... No Charge  
**Standard with AT5543 adjustable angle footplate.**
- AT5044 Elevating Legrests ..... \$255.00 E0990
- ATL4A Center Pivot Articulating Legrests ..... \$320.00 K0053

**FOOTPLATE TYPE (MVPS) (PAIR)**

Clear

- AT5543 Adjustable Angle Flip Up Footplates ..... \$165.00 K0040
- 1651 Flip Up Composite Footplates ..... No Charge
- A1562 Aluminum Footplates ..... \$59.00
- TUBULAR Individual Tubular Footplates ..... \$59.00
- A1564 Extra Large Aluminum Footplates ..... \$59.00 K0041

**FOOTPLATE TYPE (MVPF60 AND MVPF80)**

Clear

- AT5540 1 Piece Tubular Flip Up (each) ..... No Charge
- AT5545 1 Piece Adjustable Angle Flip Up ..... \$165.00 K0040
- AT5548 1 Piece Plastic Coat Tubular Flip Up ..... \$40.00

**MAX KNEE-TO-HEEL LENGTH**

Clear

- KTH14 14" Knee-To - Heel ..... No Charge
- KTH1516 15" - 16" Knee-To - Heel ..... No Charge
- KTH1618 16" - 18" Knee-To - Heel ..... No Charge

**FOOTREST/LEGREST ACCESSORIES (PAIR)**

- A1337 Calf Strap ..... \$24.00 K0038
- FRMPRBK ChairArmor Black Neoprene Frame Protectors . \$65.00  
**Not available with SWING914, 70MFX, 90SWING, AT5044 or ATL4A.**
- FRMPRB ChairArmor Blue Neoprene Frame Protectors ... \$65.00  
**Not available with SWING914, 70MFX, 90SWING, AT5044 or ATL4A.**
- FRMPRR ChairArmor Red Neoprene Frame Protectors .... \$65.00  
**Not available with SWING914, 70MFX, 90SWING, AT5044 or ATL4A.**
- FRMPRS ChairArmor Silver Neoprene Frame Protectors . \$65.00  
**Not available with SWING914, 70MFX, 90SWING, AT5044 or ATL4A.**
- ALPT Longer Pivot And Slide Tubes ..... No Charge
- A1600BK Heel Loops With Ankle Straps ..... \$51.00 E0951
- A9342 "H" Strap ..... \$51.00 K0039
- HLOOP Heel Loops Without Ankle Straps ..... \$51.00 E0951

Knee-to-Heel Range		
Front Rigging Type	Minimum	Maximum
70° Tapered	15 3/4"	19 3/4"
70°	13 1/2"	18 1/2"
70° MFX	8 3/4"	14 3/4"
Elevating	15 1/4"	21 1/4"
90°	5 1/4"	10 1/4"
60°	14 3/4"	18 1/2"
60° MFX	11.00"	15.00"
Range adjusts in 1" increments.		
Measurement taken from top of seat rail to back of footrests.		

**SEATING**

**SEATING CUSHION SELECTION**

Clear

For additional information (size availability, additional seating options and pricing) see form # 07-027 Rehab Seating and Positioning Products price list.

- EC Absolute Cushion ..... \$104.00 E2601  
**Available in even sizes only. See form 07-027 for more information.**
- ITSG InTouch™ Stabilite™ Gentle Contour Cushion .... \$340.00 E2607
- ITSM InTouch Stabilite Max Contour Cushion ..... \$340.00 E2605
- ITFG InTouch Flovair™ Gentle Contour Cushion ..... \$360.00 E2607
- ITFM InTouch Flovair Max Contour Cushion ..... \$360.00 E2607

**SEATING BACK SELECTION**

Clear

For additional information (size availability, additional seating options and pricing) see form # 07-027 Rehab Seating and Positioning Products price list.

- PBR Personal Back 10 Reg. Height ..... \$565.00 E2615
- PBT Personal Back 10 Tall Height ..... \$565.00 E2615
- PPR Personal Back 10 Plus Reg. Height ..... \$645.00 E2615
- PPT Personal Back 10 Plus Tall Height ..... \$645.00 E2615
- ITPR InTouch Propel™ Back ..... \$438.00 E2611
- ITPCSR InTouch PCS Regular ..... \$635.00 E2620 E2621
- ITPCST InTouch PCS Tall ..... \$635.00 E2620 E2621

**INTOUCH PCS™ ACCESSORIES**

- PCSPS Pelvic Stabilizers, 3" (Pair).....\$230.00
- PCSTL Thoracic Lateral Supports, 4" (Pair).....\$230.00

For additional seating needs please refer to the following:  
Invacare Rehab Seating and Positioning Products Price List Form  
#07-027.

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**

**Features**

Basic Chair Includes:

- Swingaway or Fixed Front Folding Frame
- Curved, Square or Extended Rear Frame
- Crossbrace Selection
- Seat Widths (15" - 19" Wide)
- Seat Depth (13" -19" Deep)
- Arm Type (AT99 or AT99N)
- A1418, 1014, A1216, S1418, FBH12, FBH14, FBH16, FBH18 or FBH20
- Frame Finish
- Slip on Back Upholstery
- Caster Selection (51DC, 61DC, 615DC, 81DC, 815DC)
- Fork Type
- Wheel Options 16COM,18COM,20COM,22COM.,24COM7, 20SP,22SP,24SP
- URETH Urethane (treaded) Tires
- ALUMPC Aluminum Anodized Handrims
- STDPT Standard Axle Plates
- Permanent Axles or Quick Release Axles
- Camber Selection (0 ,3, 6 or 9 Degree)
- Front Rigging (SWING1319, SWING914, 70NTAPER, 70MFX, 70STAPER, 90SWING)
- Footplate Type (1651 or AT5540)
- Knee to Heel Length Selection
- Seat Upholstery
- Wheel Locks (TSAVER, A1567 or 9691)

**Weight Limitation 250 Pounds**

Invacare recommends all Invacare wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

**Available Through Our Specials Department**

- One Arm Drive SM - 010502123533
- Oxygen Holder SM - 010227014237
- 80° Swingaway Front Riggings SM - 200411121756
- Attendant Operated Wheel Locks SM - 040302051903
- Hide away undermount wheel locks. C - 030402014153
- Longer Pivot Tube for 90SWING. SM- 201229021218

All specials are subject to final approval. Some special quotes are not compatible for all frame configurations. Please contact Customer Service for details and prices.

**Specifications**

Frame Style	Folding
Seat Widths	15" - 22"
Seat Depths	9" - 22"
Minimum Overall Width	Seat width plus 10"
Adjustable Angle Back	-15°, 10°, -5°, 0°
Adjustable Height Back	10" - 18"
Camber Angles	0°, 3°, 6°, 9°
Tubing O.D.	1" Aluminum
Weight Limit	250 lbs.
Weight Of Chair	17 lbs. (16" x 16" without wheels and accessories)
Heavy Duty Package	350 lbs.
Warranty	Lifetime on frame and crossbraces

Comments

**NOTE: All specifications and dimensions are approximate.**  
 For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

**Invacare Corporation**

[www.invacare.com](http://www.invacare.com)

**USA**

One Invacare Way  
 Elyria, OH USA  
 44035-4190  
 800-333-6900

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Spinergy is a registered trademark of Spinergy Inc. FloGel is a registered trademark of S.N.F. Specifications and prices are subject to change without notification. ©2009 Invacare Corporation

