

Top End® Pro™ BB Wheelchair

Price List and Order Form

Effective January 15, 2010 Revised January 1, 2015

For accuracy of ordering, please fax order to Top End at 727-522-1007 or e-mail to TopEndOrders@invacare.com

For questions call 800-532-8677 or 727-522-8677.

To order Literature: 800-828-6282

www.topendwheelchair.com



Provider Name _____ Date _____ Top End Job Number _____

Account # _____ P.O. # _____ Phone # _____ Fax # _____

Address _____ City _____ State _____ Zip Code _____

Provider Primary Contact: _____ Email: _____ Phone No _____

Athlete Name _____ Email _____ Phone No _____

Address _____ City _____ State _____ Zip Code _____

Indicate: ☐ Order or ☐ Quote Ship order to: ☐ Provider ☐ Customer Order will ship to Provider unless otherwise noted.

☐ SPBB: Top End Pro Basketball Aluminum Wheelchair 2,395.00


1. SEAT WIDTH: ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20"

2. SEAT DEPTH: ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18"

3. ADJUSTABLE SEAT HEIGHT: choose tall or short

TALL: Front Seat Height 19-21": set at 19" Standard
Rear Seat Height 15-21" : set at 17"

☐ SHORT: Front Seat Height 17-19": set at 18" No Charge
Rear Seat Height 13-17" : set at 16"



Wheel Size Guide			
Rear Seat Height	16.0" or less	16.5" - 19.0"	19.0" - 21.0"
Recommended Wheel Size	24"(540)	25"(559)	26"(590)

wheel size guide

REAR WHEEL SIZE: ☐ 24"(540) ☐ 25"(559) ☐ 26"(590)

FOOTREST: Platform Footrest (w/ Velcro strap) Standard¹

☐ TER124 Footrest Click strap upgrade 75.00

5. AXLE POSITION: *CENTER OF GRAVITY*²

Adjustable Axle Position: 3.5" to 5.5" Set at _____ (5.0 Std)³

6. WHEEL BASE: ☐ 15" (Short) ☐ 17" (Med) ☐ 19" (Long)

(choose one: if not chosen, manufacturer will select based on seat depth)

Wheel Base Guide			
Seat Depth	14"	15", 16"	17", 18"
Recommended Wheel Base	15" (Short)	17" (Med)	19" (Long)

wheel base guide

ANTI-TIP/5TH WHEEL: Fixed Single Anti-tip Standard

☐ TER115 Double Anti-tip (swivel), Aluminum 120.00

8. CAMBER: (choose one: if not chosen standard will be sent)


☐ 15° ☐ 18° (Standard) ☐ 20°

FOOTNOTES

¹ Platform Footrest includes Adjustable Angle Plate with Toe Stop, and Velcro Strap. Footrest adjusts Up/Down, Fore/Aft, Angle..

² Measure the user's existing chair from the Front of the Backpost, horizontally (as shown in side view pg 1), to the center of the rear axle.

³ Adjustable Center of Gravity: Non-folding, angle adjustable backrest. If side guards are desired, Clamp-on style is the only compatible option.



WHEELS:

Sun® High Performance Wheels Standard

☐ TER131 Spinergy® SLX Upgrade (4mm/24 spokes) 895.00

Spoke color: (choose one-if color is not chosen, spokes will be black)

☐ Black ☐ Yellow ☐ White ☐ Red ☐ Blue

☐ TER143 Sun® Equalizer Sport Wheel Upgrade

(48 spoke, cross lacing pattern with 12 tab pushrims for 25" & 26")

(36 spoke, cross lacing pattern with 6 tab pushrim for 24") 250.00

TIRES:

Primo high pressure clincher Standard

☐ TER101 Kenda® Tires: ☐ Red or ☐ Yellow 75.00

AXLES:

Quick release Standard

14. PUSHRIMS, ANODIZED (Aluminum)

Standard mount Standard

☐ Close mount NC

☐ TER40 Cut Tabs Aluminum pushrim 40.00

☐ TER129 12-Tab upgrade*

*Available on 25" or 26" Spinergy wheels only, no 24" 150.00

11. BACK HEIGHT: (Select one)

☐ Short 9" - 11" ☐ Low 10" - 14" ☐ Med 14" - 18

BACKREST/SIDEGUARDS:

NON-FOLDING BACKREST/ADJUSTABLE ANGLE:

Fixed BackrestStandard

SIDEGUARD: Compatible with Standard Fixed Backrest

☐ TER95 Clamp-On Fixed aluminum side guards 125.00

FOLDING BACKREST OPTION:

☐ TER3 Folding back 250.00

SIDEGUARD: Compatible with Folding Backrest

☐ TED 121 Fold-in aluminum sideguards 125.00

SAFETY AND COMFORT OPTIONS:

Calf Strap Standard

☐ TER35 Frame protector – rear 20.00
(required for tournament play)

☐ TER21 Spokeguards: ☐ Clear or ☐ Black 75.00
(If color not chosen, black will be sent)

☐ TER14 Cushion, soft ☐ 2" ☐ 3" 75.00

☐ TER33 Cushion, hard ☐ 2" ☐ 3" ☐ 4" 125.00

☐ TER148 Basketball carrying bag \$75.00

☐ TER126 Click strap large (pad length ~12") 86.00

☐ TER125 Click strap medium (pad length ~10") 80.00

☐ TER124 Click strap small (pad length ~8") 75.00

☐ TED72 Positioning Strap – Hook & Loop 40.00

☐ TER128 Metal O-Ring positioning belt 60.00

☐ RACE22*Top End T-Shirt, ☐ L ☐ XL 10.00

*No discount item

UPHOLSTERY:

Screw-on Seat Upholstery Standard

Adjustable Tension Back UpholsteryStandard

☐ TER149 Back upholstery upgrade: 3"leather flap aluminum plugs
..... \$100.00

UPHOLSTERY COLOR:

BlackStandard

☐ UPH1 Woodland Camouflage 100.00

☐ UPH2 Digital Desert Camouflage 100.00

FRAME COLORS:

CONTEMPORARY COLORS

☐ Wet Black ☐ White Snow ☐ Silver Metallic

☐ Stealth Black

WAY COOL COLORS

☐ Royal Blue ☐ Midnight Blue ☐ Grasshopper Green

☐ Electric Blueberry

TOO HOT COLORS

☐ Tangerine Red ☐ Raspberry Red ☐ Grape Madness

☐ Safety Yellow ☐ Orange ☐ Bubblegum Pink



Invacare Corporation
TOP END
Sports and Recreation Products
4501 63rd Circle North
Pinellas Park, FL USA
33781-5914

www.invacare.com

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Spinergy is a registered trademark of Spinergy, Inc. Sun is a registered trademark of Sun Components Inc.
Kenda is a registered trademark of Kenda Rubber Industrial Co. Ltd.
Specifications and prices are subject to change without notification.
©2015 Invacare Corporation