

SPECIFICATIONS

Item #	Handle to Floor (highest setting)	Handle to Floor (lowest setting)	Patient Height (highest)	Patient Height (lowest)	Product Weight (each)	Weight Limit
10410 - 10414	39"	30"	6'2"	5'0"	1.5 lbs.	300 lbs.
10403/G	37"	28"	6'2"	5'0"	2.5 lbs.	300 lbs.
10404/G	30.5"	21.5"	5'2"	4'2"	2.3 lbs.	300 lbs.
10405/G	42"	33"	6'6"	5'10"	2.6 lbs.	300 lbs.
10407	21"	14"	4'6"	3'2"	.93 lbs.	250 lbs.
10403 HD	37"	28"	6'2"	5'0"	5.7 lbs.	500 lbs.
10404 HD	30.5"	21.5"	5'2"	4'2"	4.6 lbs.	500 lbs.
10405 HD	42"	33"	6'7"	5'10"	6.8 lbs.	500 lbs.