



INVACARE® TOP END®
FORCE™ -3
HANDCYCLE



The **Invacare® Top End® Force™ -3 Handcycle** is a great handcycle at a great price. This bike is ready to ride with 27 speeds, performance wheels and Shimano® components. Perfect whether you are just looking to get started in handcycling or ready to force the pace at your next competition.

Features

15" or 17" seat width

Stealth Black Powder Coat Paint

60" Wheelbase

Adjustable back angle, seat and footrests

27 speed drive train

Shimano® Deore Rapid Fire™ shifter/brake

Shimano® 170 mm cranks

Avid® disk brakes

Shimano® cable driven upper derailleur

Adjustable crank height: diagonal and vertical

Alex DM18 26-559 wheels

0° camber

High pressure cruiser clincher tires

Seat and back cushions

Chest and footrest safety straps

Drafting bumper

Threaded axles

Cable springs

Frictionless ovalized ergonomic hand pedals

Safety flag


Chain guard

Reflectors

Tektro® parking brake

250 lb. weight capacity

Model #: FRC3



1. Wrap around frame with adjustable seat upholstery
2. Cable driven shifting for upper chainrings
3. Front disk brake and welded footrest clamps
4. Backrest adjusts forward and backward
5. Adjustable back angle allows 35-65° recline
6. Shimano components for hands-on shifting and braking



Invacare Corporation

www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2014

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Shimano, Rapid Fire and Tektro are registered trademarks of Shimano Inc. Avid is a registered trademark of SRAM Corporation. BionX is a registered trademark of BionX International Corporation. Specifications are subject to change without notification.
Form No. 13-692
131225

Options

Heart rate monitor

Alignment gauge

Backpack hydration system

Water bottle and cage

Cordless computer

Helmet

Safety lights, mirror

V-Crank upgrade

Trunk power position backrest

BionX® Power Boost

Additional safety flag

Click straps

Bike rack

Indoor training roller

Extra cable kit

Tool and tire repair kit

Wheel bag

Top End Bike Jersey

Factory Assembly



Yes, you can.®