

lifetime limited warranty

Your Drive brand product is warranted to be free of defects in materials and workmanship for the lifetime of the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

In the event of a defect covered by this warranty, we will, at our option, repair or replace the device.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.



© 2006 Medical Depot, Inc. All rights reserved.
Drive is a trademark of Medical Depot, Inc.
Port Washington N.Y. 11050 USA Made in China

exam lamp



EU Authorized Representative
Drive International, LLC
Globe Lane - Stockport
DUKINFIELD SK16 4RE, GB

www.drivemedical.com

drive
MEDICAL DESIGN & MANUFACTURING

assembly instructions

WARNING:

BULB IS NOT TO EXCEED A MAXIMUM OF 60 WATTS.

Set-Up Standard pole

1. Thread pole (A) onto adapter on base (C).
2. Thread Flex arm (B) onto pole set.

WARNING:

The lamp is equipped with a 3-prong grounding plug for protection against shock hazards. DO NOT under any circumstances cut or remove the grounding prong from the plug. The plug must be used directly with the correct 3-prong receptacle. DO NOT use an extension cord or adapter. To make alterations call a professional electrician.



REF. 13405 Shown

