

SPECIFICATIONS:

Pump:

- Power 120 VAC
- Frequency: 60 HZ
- Pressure Range 30 mmHg to 80 mmHg (inner mattress pressure)
- Air volume 8-LPM (Litres per Minute)
- Cycle Time: 10 Minutes
- Functions: Static and Dynamic Features
- Power Cord Length: 14.76 Ft.
- Pump Dimensions: 11" (W) x 8" (H) x 3.5" (D)
- UL: 60601-1, 1st Edition-2006
- Maximum Current: 0.5 Amps
- Warranty: 12 Months Limited

Mattress:

- Topper Material: zippered all around and quilted urethane coated nylon fabric (2-way stretch) is; water resistant, vapor permeable, machine washable and low shear
- Base Material: zippered PVC mesh
- 4-Air cells at head section remain static to enhance pillow function
- Dimensions: 36" (W) x 80" (L) x 8" (H)
- Maximum Support Capacity: 350 lbs.
- Warranty: 12-Months Limited