

Competitor II Manual Height Adjustable Bed

HCPCS:



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E



FIGURE F

- Large glow in the dark hand pendant has easy-to-use controls for adjusting head and foot deck sections to comfortable positions built in battery backup and power indicator light. NEW Lock out feature can be used to prevent accidental actuation of head or foot adjustment
- In the event of a power failure, a 9 volt battery (not included) located in the hand pendant can lower the head and foot sections to a horizontal position. Two additional Batteries can be installed into the motor housing to increase adjustment speed and extend battery life. No crank needed
- Slat Deck provides strength, no-sag mattress support, no links to connect, easy to clean
- Patient weight capacity is 350 lbs. Total weight capacity 450 lbs. evenly distributed. Includes occupant, mattress, side rails, trapeze, accessories, personal affects, etc.
- Ultra light head and foot frame weigh only 41 lbs. each making one person transport and assembly easy
- #15561 Competitor II Manual Hgt Adj. Bed frame only, casters optional
- Full and Half Rails meet IEC Entrapment guidelines when used on the #15561/15571 beds and are backward compatible for use with the #15560 & 15570 Competitor beds
- Product Safety Standards: Meets IEC 60601-2-52 60601-1 60601-1-2 60601-1-11

PRODUCTS	Description
15561	Compet II Man Hgt Adj Bed, 1/cs
15561FR	Compet II Man Hgt Adj Bed w/FRail, 1/cs
15561FR-PKG	Compet II ManHgtAdjBedw/FRail,Matt, 1/c
15561FR-PKG-FM	Compet II ManHgtAdjBedw/FRail,FMatt,1/c
15561HR	Compet II Man Hgt Adj Bed w/HRail, 1/cs
15561HR-PKG	Compet II ManHgtAdjBedw/HRail,Matt, 1/c
15561HR-PKG-FM	Compet II ManHgtAdjBedw/HRail,FMatt,1/c
15571F	Compet II Hi/Lo Foot bed end, 1/EA
15571H	Compet II Hi/Lo Head bed end, 1/EA
SPECIFICATIONS	Description
Weight Capacity	350 lbs (158.8 kg) Patient Wgt 450 lbs (204 kg) Total Wgt Includes occupant
Carton Shipping Weight	155 lbs. (70.3 kg)
Weight	138 lbs. (62.6 kg)
ACCESSORIES	Description
15561EXT	Compet II 4" Extension 1/EA, EA
15561TRANSC	Comp II Transporter w/4 casters 1/CV, CS
15571TRANS	Comp II Transporter Only 1/EA, EA