

Spirit Plus™

CAR SEAT



ACCESSORIES



Incontinence Covers



3-inch Extensor Thrust Wedge



Buckle Guard



Retainer Clip Guard



The Spirit Plus™ Car Seat provides positioning support for users with mild to complex positioning needs. The Spirit Plus is the only Car Seat on the market today offering swing-away padded trunk supports, hip supports and abductor.

AVAILABLE COLORS



Convertible Pink (CP)



Go Kart Green (GG)



Rally Blue (RB)



Roadster Red (RR)



Speedway Gray (SG)

BENEFITS

User Need Accommodation

Optimal Support

User Safety & Security

Accessory & Expansion Options

FEATURES

- Two sizes with seat depth extender (12–16") allow for built-in growth
- Height & weight capacities support users up to 66" in height and 130 lb. in weight
- Recline bar adds 25° of extra tilt for added comfort
- Trunk & hip supports and abductor help maintain proper positioning and postural support
- Durable, machine-washable padded cover
- Memory foam back/seat padding reduces sores and fatigue
- Open shell design with low-profile sides and swing-away supports for easy transfers
- Optional incontinence covers, anti-thrust wedge and buckle guard available

For more information, call us today!

877-224-0946 | www.inspiredbydrive.com