



Belt-positioning booster seat

Weight range: 44–175 lbs

Height range: 44–72 in.

Maximum height: 62 in.

15in. or 18 in. seat depths

Unlimited shoulder width

Unlimited thigh width

Super low profile seat
simplifies transfer of
occupant into vehicle

Scan for
more details!




The Churchill™ with POSITIONING HARNESS

AVAILABLE & CUSTOM OPTIONS:



BUCKLE
GUARD



CHEST CLIP
GUARD



VEST



EZ-UP
HEADREST




EZ-UP CAP



STAY-PUT
POMMEL



HIP
INFLECTION
WEDGE



QUICK-CHANGE
INCONTINENT
COVER

...MORE OPTIONS AVAILABLE!



www.merrittcarseat.com



(317) 409-0148



info@merrittcarseat.com