

Elevate Upright Walker

HCPCS:



The Drive Elevate Upright Walker brings users increased balance and stability while promoting better posture and providing a more natural line of sight than traditional mobility aids. Height-adjustable upholstered forearm pads keep arms supported and level for comfortable use, while large, ten-inch wheels make it perfect for indoor and outdoor use.

- Elevate Upright Walker is the stand-up mobility aid that's easy to take on the go.
- Padded armrests and upright handle adjust to 44"
- Durable and comfortable padded seat
- Simple one-handed, side-to-side folding
- Includes removable zipper storage bag with strap
- Dual front fork for optimal stability and
- Loop lock brakes for added safety and control
- Adjustable length dual handbrake system
- Large 10 multi-terrain front wheels for
- Height adjustable backrest
- *Sit to stand assist handles currently on blue model only*

PRODUCTS	Description
UPWKRBL10	UPRIGHT WALKER, ROLLATOR, 1EA, BLUE
UPWKRGY10	UPRIGHT WALKER, ROLLATOR, 1EA, GREY
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	24.8 lb
Material	Aluminum
Seat Dimensions	18" (W) x 10" (D)
Handle Height	39" - 44"
Weight Capacity	300 lb
Seat To Floor Height	21"
Casters	Front: 10 Rear: 8
CA Prop 65	no
Handle Extension Range	4" - 6"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]