



Designed for versatility.

The Golvo® patient lift is unique in the fact that it combines the functional advantages of an overhead lift with the flexibility of a mobile lift. The lift strap and great range of accessories make the Golvo the perfect lift in a variety of situations.

Key features at a glance

- Flexible lift strap design that provides several functional benefits for patients and caregivers
- Parallel base-width adjustment makes it easier to lift the patient to/from chairs, beds and the floor. Optional LowBase version provides even more space under low furniture
- Excellent lift interval – ideal for ambulation and lifting from the floor
- Armrests are standard and facilitate patient mobilization
- Designed in aluminum – strong, durable and environmentally friendly
- Light-weight, environmentally-friendly NiMH-battery that is easy to charge and easy to handle
- Can be equipped with quick-release hook system to quickly and easily change the sling bar or install a scale
- Dual control options to operate the lift – by hand control or control panel
- Intelligent diagnostic system informs caregivers that the lift is safe to use and when service is required

The perfect lift in a variety of care situations.

The unique lift strap design and flexible lifting interval makes the Golvo® mobile lift ideal for many different care environments, such as rehabilitation clinics, ICU wards, childcare centers and the operating room.

Together with Hill-Rom's wide range of functional Liko® slings and accessories, the Golvo lift is ideal for those with high demands and wide-ranging needs.

The Golvo patient lift is equipped with an intelligent diagnostic system that helps caregivers know when the lift is safe to use and when it's time for service.

The lightweight, environmentally-friendly NiMH-battery is easy to charge and easy to handle.



	<u>Golvo® 8000 Lift</u>	<u>Golvo® 8008 Lift</u>	<u>Golvo® 8008 LowBase Lift</u>
Product number	2000034	2000035	2000039
Maximum load	440 lbs	440 lbs	440 lbs
Weight	92 lbs	97 lbs	101 lbs
Lifting speed	1.61"/s	1.61"/s	1.6"/s
Alternative lifting speed	1.3"/s	1.3"/s	1.3"/s
Maximum lifting height	68.3"	73.2"	73.2"
Base height	4.1"	4.1"	2.4/4.1"
Width	25.7"	28.9"	29.3"
Length	43.9"	46.6"	46.6"
Material	Anodized aluminum	Anodized aluminum	Anodized aluminum
Emergency lowering	Mechanical and electrical	Mechanical and electrical	Mechanical and electrical
Base-width adjustment	Electrical	Electrical	Electrical

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2013 Liko R&D AB. ALL RIGHTS RESERVED.

ORDER NUMBER 179881 rev 1 E 05/09/2013 ENG – US

Zinc Number US-AC-1212-0077

www.hill-rom.com

USA **800-445-3730**
CANADA **800-267-2337**