



## Easy lifting in small spaces.

The Viking® S and XS mobile lifts are small and narrow mobile lifts designed for care environments where lifting space is particularly restricted. Like all lifts in the popular Viking series, these lifts are made of lightweight, environmentally friendly and easy-to-clean aluminum, providing excellent weight-to-capacity ratio.

### Key features at a glance

- Designed in aluminum – strong, durable and environmentally friendly
- Suitable for most common lifting needs in many types of care environments
- Ideal when space is limited – for example in home-care situations or extended care facilities
- Three different height settings (Viking S only) to adapt the lift to different lifting situations
- Can be used together with Hill-Rom's wide range of slings and lifting accessories

# The right lift for many daily situations.

The Viking® mobile lifts are designed for most common lifting situations and can be used for lifts to and from the toilet, between bed and wheelchair, and from the floor.

The Viking S and Viking XS lifts can be used with most Liko® slings and other accessories, providing excellent patient comfort and security in every situation.



	<b>Viking® S Lift</b>	<b>Viking® XS Lift</b>
<b>Product number</b>	2040006	2040007
<b>Maximum load</b>	350 lbs	350 lbs
<b>Weight</b>	59 lbs	56 lbs
<b>Lifting speed</b>	1.18"/s	1.18"/s
<b>Maximum lifting height</b>	68.9"	62.0"
<b>Base height</b>	4.3"	4.3"
<b>Width</b>	25.4"	20.9"
<b>Length</b>	49.6"	43.5"
<b>Material</b>	Aluminum	Aluminum
<b>Emergency lowering</b>	Mechanical and electrical	Mechanical and electrical
<b>Base-width adjustment</b>	Manual	Manual

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2013 Liko R&D AB. ALL RIGHTS RESERVED.

ORDER NUMBER 179882 rev 1 E 05/09/2013 ENG – US

Zinc Number US-AC-1212-0079c

[www.hill-rom.com](http://www.hill-rom.com)

USA 800-445-3730

CANADA 800-267-2337