



INVACARE® TOP END®
PRELIMINATOR™
RACING WHEELCHAIR



Designed for the up and coming wheelchair racer, the **Invacare® Top End® Preliminator™ Racing Wheelchair** offers exceptional stability and freedom to compete. The Preliminator Racing Wheelchair comes as a set configuration for schools, sports clubs, and weekend racers or a more customized configuration for the more competitive user.

The Invacare® Top End® Preliminator™ Racing Wheelchair is for active youth who want to participate in track racing with their family and friends. This sports chair is designed to help a child or young adult to begin their journey into an active lifestyle.

- Adjustable seat depth, foot plate, and back angle allow for growth
- Telescoping, removable anti-tipper provides safety for novice racers
- Easy to order based on one measurement which allows for quick delivery
- Accessories include bike rack, click strap, heart rate monitor, helmet, racing suit, safety light, training roller, water bottle, wheels, and a wheel bag

**Invacare® Top End® Preliminator™
Mini Racing Wheelchair**

**10" Seat Width
Gecko Green
Model #: 10028**



**Invacare® Top End® Preliminator™
Prep Racing Wheelchair**

**12" Seat Width
Safety Yellow
Model #: 10029**



**Invacare® Top End® Preliminator™
Junior Varsity Racing Wheelchair**

**14" Seat Width
White Snow
Model #: 10030**

**Invacare® Top End® Preliminator™
Varsity Racing Wheelchair**

**16" Seat Width
Simple Red
Model #: 10031**



**Invacare® Top End® Preliminator™
Adult Racing Wheelchair**

**18" Seat Width
Simple Blue
Model #: 10032**

Racing Wheelchair Model	10028 Mini Preliminator	10029 Prep Preliminator	10030 Junior Varsity Preliminator	10031 Varsity Preliminator	10032 Adult Preliminator
Seat Dimensions	10" W x 10-12" D	12" W x 10-12" D	14" W x 13-15" D	16" W x 13-15" D	18" W x 13-17" D
Front seat-to-floor height	Adjustable from 17-21" in ½" increments	Adjustable from 17-21" in ½" increments	Adjustable from 17-21" in ½" increments	Adjustable from 17-21" in ½" increments	Adjustable from 17-21" in ½" increments
Rear seat-to-floor height	Adjustable from 13.5-15" in ½" increments	Adjustable from 13.5-15" in ½" increments	Adjustable from 13.5-15" in ½" increments	Adjustable from 13.5-15" in ½" increments	Adjustable from 13.5-15" in ½" increments
Back Height	Adjustable from 10-14" in 1" increments	Adjustable from 10-14" in 1" increments	Adjustable from 10-14" in 1" increments	Adjustable from 10-14" in 1" increments	Adjustable from 10-14" in 1" increments
Weight Capacity	250 lb	250 lb	250 lb	250 lb	250 lb
Color	Gecko Green	Safety Yellow	White Snow	Simple Red	Simple Blue
Limited Warranty	Frame – 3 years Parts – 1 year	Frame – 3 years Parts – 1 year	Frame – 3 years Parts – 1 year	Frame – 3 years Parts – 1 year	Frame – 3 years Parts – 1 year

Racing Wheelchair Model	10033 Custom Preliminator
Seat Dimensions	10-18" W x 10-17" D (width and depth measured in 1" increments)
Front seat-to-floor height	Adjustable from 17-21" in ½" increments
Rear seat-to-floor height	Adjustable from 13.5-15" in ½" increments
Back Height	Adjustable from 9-11" 10-14" 12-16" and 14-18"
Steering	Horizontal or vertical steering handle with triangle track control bottom or top mounted
Weight Capacity	250 lb
Color	Tangerine Red, Raspberry Red, Simple Orange, Safety Yellow, Gecko Green, Grasshopper Green, Midnight Blue, Royal Blue, Electric Blueberry, Grape Madness, Bubble-Gum Pink, White Snow, Stealth Black, Wet Black, Black Prism, Silver Metallic, Stealth Grey
Limited Warranty	Frame – 3 years Parts – 1 year

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End
www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2017
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and ©.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted.
Form No. 16-279
161235

