


# SWING-AWAY AMP

**HCPCS CODING:** E1020  
**RESIDUAL LIMB SUPPORT**



## PRODUCT FEATURES

- + Designed to install on all popular wheelchairs within minutes.
- + Provides patients a fully adjustable amputee support that will coincide with any wheelchair cushion.
- + Allows patients to be transferred in and out of a wheelchair within seconds.
- + Installs onto either left or right side of wheelchair.
- + Fits tubing sizes from 7/8" to 1".
- + Foam construction has soft entry layer for comfort and firmer base layer for support covered with incontinent proof cover.
- + Amputee pad has adduction to help center residual limb.



## VERTICAL ADJUSTABILITY



- Toggle lever allows Amp pad to be positioned at multiple angles.
- Up to four inches of height adjustment to accommodate varying seat cushion thicknesses.
- Pad can be adjusted three inches in depth to accommodate wheelchair depth conversion kits.

## QUICK RELEASE SWING-AWAY HARDWARE

PATENT PENDING

- Quick release unit allows for fast and simple swing away.
- Swing away feature allows patients to be easily transferred in and out of a wheelchair.
- Pad can be removed and replaced from wheelchair.



## STANDARD COVER



**COMFORT-TEK FABRIC:** the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multi-dimensional stretch and helps alleviate pressure on bony prominences.

### SWING-AWAY AMP 10 INCH EXTENSION

model number	width x length	thickness
<b>AMPSA810</b>	8 x 10	3"
<b>AMPSA910</b>	9 x 10	3"
<b>AMPSA1010</b>	10 x 10	3"
<b>AMPSA1110</b>	11 x 10	3"
<b>AMPSA1210</b>	12 x 10	3"

### SWING-AWAY AMP 14 INCH EXTENSION

model number	width x length	thickness
<b>AMPSA814</b>	8 x 14	3"
<b>AMPSA914</b>	9 x 14	3"
<b>AMPSA1014</b>	10 x 14	3"
<b>AMPSA1114</b>	11 x 14	3"
<b>AMPSA1214</b>	12 x 14	3"

WARNING + SWING AWAY AMP SHOULD NOT BE USED FOR PATIENT TRANSFERS.