

Invacare® Top End® Pro™ 2 Wheelchair



Product Description:

The Pro™-2 Basketball Wheelchair is in stock and ready to ship. Designed for aspiring athletes who want durability and quick, easy adjustments, all at a great price. Adjustments that can take up to four hours in other chairs, take only 15 minutes in the Pro 2 Basketball wheelchair. Plus, it is so easy to order that you only need to measure seat width. Now that is easy!

About New Smart Forms:

Top End is excited to introduce a new way to order Top End products. Our new Smart Forms provide an easy and electronic way to fill out order forms, save copies, and submit to them via e-mail. Please see below for the features and useful information.

How to use:

- Use any Adobe Acrobat application to electronically fill out the order form
- Use the default selections to order a product with the most common configuration
- Use the **"Reset"** button located in the forms footer to reset selections to their defaults
- Use the **"Clear All"** button located in the forms footer to clear selections, including defaults. This feature allows for a blank order form to be printed and completed by hand
- Use the **"E-mail to Top End order box"** button located on the last page of the form to e-mail the completed form directly to the topendorders@invacare.com e-mail inbox
- Forms can be completed electronically, printed, and e-mailed or faxed to topendorders@invacare.com or 727-522-1007 if access to Acrobat Reader for Windows is unavailable

(Acrobat Reader for Windows with java-script enabled is required to save and attach to e-mail)

Features:

- Tabulated options for an organized and standardized work-flow
- New more specific option codes which will appear on confirmations to improve customer / manufacturing communication
- Common special requests have been given option codes
- Built in exclusions to help eliminate errors
- Ability to save the form after it is filled out
- Ability to e-mail the form inside Acrobat Reader rather than saving a copy and attaching to an e-mail
- Ability to clear the entire form so it can be printed and completed by hand in the field when a computer is not available
- Ability to reset form to default selections

Invacare® Top End® Pro™ 2 Wheelchair



PROVIDER / CUSTOMER INFORMATION:			
Provider Information:		End User Information:	
Date:	Name:	Height:	
IVC Acct #:	Ship to Address:	Weight (lbs.):	
Contact:	Address Line 2:	Leg Length (in.):	
Phone #:	City:	Hip Width (in.):	
Fax/Email:	State:	Zip:	Phone #:
☐ Ship to Provider	☐ Ship to End User		Email:

ORDER REQUEST TYPE:		
☐ Quote	☐ Order	PO #

BASE MODEL:		
☒ 10014	Top End Pro 2 - All Sport	\$2,195

Assembly: <i>(see definitions)</i>		
☒ 10048	Assembled	\$0
Shipping:		
Shipping Fee <i>(Continental US Only - See definitions)</i>		\$150

SIZING:		
SEAT WIDTH:		
☐ 15005	10" Seat Width - Neon Green	\$0
☐ 15009	12" Seat Width - Safety Yellow	\$0
☐ 15013	14" Seat Width - White Snow	\$0
☐ 15017	16" Seat Width - Simple Red	\$0
☐ 15021	18" Seat Width - Royal Blue	\$0
SEAT DEPTH: <i>(See Definitions)</i>		
☒ 15357	Adjustable Seat Depth	\$0
SEAT HEIGHT:		
☒ 15096	Adjustable Front Seat Height: 18" - 20" Adjustable Rear Seat Height 15" - 20"	\$0
BACK HEIGHT:		
☒ 15136	Adjustable Back Height 11" to 14"	\$0

BACK TYPE:		
☒ 15159	Angle Adjustable, Non-Folding Backrest	\$0

FOOTREST:		
☒ 60003	Clamp-On Tubular, Height Adjustable Footrest	\$0

CAMBER ANGLE:		
☒ 45008	15 Degree Camber	\$0

WHEELBASE: <i>(Front to Rear Axle)</i>		
☒ 75068	17" Wheelbase	\$0

CENTER OF GRAVITY: <i>(Front of Backrest to Center of Axle)</i>		
☒ 45073	Adjustable 2.5" to 5" Center of Gravity	\$0



WHEEL SETS:		
WHEEL SIZE:		
☒ 50191	25" (559)	\$0
WHEEL TYPE:		
☒ 50196	High Performance	\$0
HANDRIM TYPE:		
☒ 50297	6 Tabs, Aluminum, Anodized	\$0
HANDRIM MOUNTING:		
☒ 50308	Standard Mount	\$0
TIRES:		
☒ 50216	Primo Grey	\$0

SIDE GUARDS:		
☒ 20014	Screw-On, Aluminum Side Guards	\$0

ANTI-TIP:		
☒ 20030	Single Telescoping, Removable, Height Adjustable Anti-Tip	\$0

Invacare® Top End® Pro™ 2 Wheelchair



UPHOLSTERY:		
MATERIAL & COLOR:		
☒ 30001	Nylon, Black Upholstery	\$0
SEAT UPHOLSTERY:		
☒ 30019	Screw-On Seat Upholstery	\$0
BACK UPHOLSTERY:		
☒ 30012	Adjustable Tension back upholstery	\$0
CUSHIONS:		
☒ 30005	2" Soft Seat Cushion	\$0

ACCESSORIES:		
STANDARD:		
☒ 25066	Calf Strap	\$0
☒ 25065	Rear Frame Protector	\$0
☒ 25060	Hook and Loop Positioning Strap	\$0
UPGRADABLE:		
POSITIONING:		
☐ 60030	Click Strap, Single Pad for Footrest	\$75
☐ 25001	Click Strap, 8" Small	\$75
☐ 25002	Click Strap, 10" Medium	\$80
☐ 25003	Click Strap, 12" Large	\$90
WINGS:		
☐ 25074	Clamp-On Basketball Wing Set	\$600
APPAREL:		
☐ 25084	Lemon Yellow Top End T-Shirt	\$20
☐ 25085	Grey Top End T-Shirt	\$20
☐ 25086	Grey Top End Racer Back Tank Top	\$20
<i>Please select a size for your shirt selection:</i>		
☐ 25087	Extra Small	
☐ 25088	Small	
☐ 25089	Medium	
☐ 25090	Large	
☐ 25091	Extra Large	
☐ 25092	XX Large	
☐ 25093	XXX Large	

DEFINITIONS:	
Assembly and Shipping:	Arrives factory assembled, rear wheels and some accessories will need to be installed. Shipping rate is \$150 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
Adjustable Seat Depth:	All Invacare Top End Pro-2 All Sport Wheelchairs have adjustable seat depths. For seat widths of 10" and 12" the seat depth will be adjustable to depths of 10" and 12". For seat widths of 14", 16", and 18" the seat depth will be adjustable to depths of 13", 15" and 17".
Clamp-on Basketball Wing:	The Clamp-on Basketball Wing set components extend the protective bumper around the front of the chair. This is typically desired for use in Wheelchair Basketball. This option is not assembled on the chair and requires minimal assembly upon arrival.
Accessories:	All Additional Accessories may be shipped in a separate box than the Pro-2. There will be no additional shipping charges when accessories are purchased with a product.
Upgradable Accessories:	All Add-on accessories may be boxed separately and require self installation.

Warranty Info:

Invacare warrants the frames when purchased new and unused to be free from defects in materials and workmanship for a period of two (2) years from the date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Invacare warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Invacare or a dealer, with a copy of the seller's invoice required for coverage under this warranty

E-mail to Top End Order Box

Invacare Corporation Top End Division

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.

©2015 Invacare Corporation

1181428 REV G KPIPARO 11-4-2019

Clear All

Reset Form

Page of