

# Forearm Crutches

FAX TO 310-305-1718 OR SUBMIT ONLINE









Dealer: \_\_\_\_\_ Marked For: \_\_\_\_\_ Account #: \_\_\_\_\_

Date: \_\_\_\_\_ PO #: \_\_\_\_\_ Ordered By: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_ Email: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims. Referred HCPCS codes are also not guarantee of payment.

BASE UNIT	FC100-2G	FC200-2G	FC300-2G
	\$128.44	\$128.44	\$140.22
<b>Product Size</b>	Small	Medium	Large
<b>Weight Capacity</b>	120 lb.	160 lb.	185 lb.
<b>User Height</b>	2'6"-3'5"	3'2"-4'5"	4'4"-5'5"
<b>Floor-to-Grip Height</b>	15-22"	19-26"	24-34"
<b>Grip-to-Top of Cuff</b>	5.25-6.25"	6.25-7.25"	7.5-8.5"
<b>Cuff Diameter</b>	2.6"	2.6"	3.5"
<b>Cuff Style</b>	Full	Full	Full
<b>Product Weight</b>	15 oz.	16 oz.	18 oz.
<b>HCPCS</b>	E0110	E0110	E0110
<b>Color</b>	 Knight Blue (B)  Wizard Purple (P)	 Knight Blue (B)  Wizard Purple (P)	 Knight Blue (B)  Wizard Purple (P)