


Softform[®]

Premier Active 2

with RX control unit

Our latest hybrid support surface Softform Premier Active 2 provides alternating pressure support that self-adjusts to the needs of the patient.



A new hybrid therapeutic support surface

For personalized pressure injury management



Yes, you can.®

A Reactive and Active Support Surface



Pressure redistribution

- › Modular construction with a high density foam base, alternating air cells for active therapy and high quality, teardrop castellated foam top layer.
- › Easy 'step-up and step-down' therapy with the addition of the RX control unit for further support in pressure injury management.



Clinically proven

- › Excellent pressure redistribution and shear management as found in independent test with Berlin Cert.*

* Berlin Cert is an expert independent test facility with specific experience in support surfaces.



Quality assured

- › Built upon the renowned Softform legacy worldwide, Softform Premier Active 2 is a tried and tested solution.
- › Softform is one of our most popular therapeutic support surfaces and a trusted solution in pressure injury management.

Technical data ▼

Softform Active 2 RX



Width ▼

36, 39, 42, 48"



Length ▼

80, 84"



Height ▼

6"



Total Product Weight ▼

35 lbs¹



Max. User Weight ▼

500, 650 lbs



Max. Cleaning Temperature ▼

95° C

Control Unit



Length ▼

10.8"



Width ▼

6.1"



Height ▼

4.1"



Total Product Weight ▼

3.9 lbs

1. Product weight may vary depending on size.

Tel: (800) 333-6900
www.invacare.com

Invacare Corporation
Worldwide Headquarters
1 Invacare Way
Elyria, OH 44035-4190



Yes, you can.®

©2023 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Please refer to the User Manual for instructions for use. Form No. 23-022 Job 230163 US-10/2023.