

TiLite 2GX Order Form

Effective November 15, 2019 (New Items in Red)

FOLDING FRAME

Page 1 of 7

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:	User Weight:	Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	State: Zip:
Ship Date:	Order #:	Phone:	Fax:

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Please Email or Fax to the Permobil Manual Customer Service Team at tilitesales@permobil.com or 866.586.2413 / 509.586.2413



FRAME STYLE

G2FS1	<input type="checkbox"/> 2GX Titanium Frame (Full Adjustability. 265 Pound Weight Limit. Light Weight, High Strength and Vibration Dampening. Not Available with Frame Transit Tie-Down.)	\$4,265	K0005
XFS1	<input type="checkbox"/> Swing Away	STD	
XFS3	<input type="checkbox"/> Hemi Height Frame (Available with a Front Seat Height of 13.5" - 20" (34-51cm). Not Available with Depth Adjustable Back/Frame, 12" or 13" Seat Depth or Heavy Duty Frame.)	NCO	
XFS2	<input type="checkbox"/> Fixed Front	NCO	K0005
XHDR1	<input type="checkbox"/> Heavy Duty Frame (350 Pound Weight Limit. Includes Tension Adjustable by Straps Back Upholstery. Available with 18" to 22" Seat Widths and 16" to 20" Seat Depths. Not Available with 90° front end)	\$435	
XPWR1	<input type="checkbox"/> Power Adaptable Frame (No charge when 20" Rear Seat Width A and 20" Seat Depth B are selected. Requires LiteSpeed Billet Aluminum Front Wheels and 0 degree Camber. Not available with Single-Sided Forks.)	\$340	
XRFR1	<input type="checkbox"/> Reinforced Frame (No charge when 20" Rear Seat Width and 20" Seat Depth are selected.)	\$340	

SMARTDRIVE POWER ASSIST SYSTEM



XMAXSD1	<input type="checkbox"/> Permobil Smart Drive MX2+ (SwitchControl comes standard. Not available with rear seat heights less than 13.5").	\$7,131	E0986
XMAXSD6	<input type="checkbox"/> NEW PushTracker E2 OR XMAXSD7 <input type="checkbox"/> PushTracker	NCO	
XMAXSD8	<input type="checkbox"/> SwitchControl Button with Mono Jack Option	\$150	
XMAXSD4	<input type="checkbox"/> Carrying Bag	\$78	

TRANSIT OPTIONS

XWTB	<input type="checkbox"/> Transport Bracket Option (Consists of four Tie-Down Brackets. Unoccupied Use Only.) (Not Available when Swing Away frame and Heavy Duty Frame are both selected. Not available with Swing Away Tubular Armrest or Swing Away -Flip Back Tubular Armrests. Requires a minimum seat depth of 14" or 15" seat Depth with Hemi Height Frame.)	\$350	
------	--	-------	--

FRAME FINISHES

Natural Finishes

XFF1	<input type="checkbox"/> Satin (Available with Titanium Frame only.)	STD	
XSCBT1	<input type="checkbox"/> Satin Finish Buffing Pads	\$20	
XFF3	<input type="checkbox"/> Tattoos (Tattoos on front frame tubes. Titanium frame finished in Satin. Aluminum frame tattooed in blasted aluminum with clear coat. Tattoos will be trimmed and adjusted to fit your frame.)	\$335	
XFF45	<input type="checkbox"/> Need For Speed		
	XFF44 <input type="checkbox"/> Flower Garden		
XFF43	<input type="checkbox"/> Mosaic		



High Gloss Paint	<input type="checkbox"/> TFF8 Super Wet Black	<input type="checkbox"/> TFF10 TiLite Red	<input type="checkbox"/> TFF29 Neon Pink
Metallic Paint	<input type="checkbox"/> TFF14 Ocean Blue Metallic	<input type="checkbox"/> TFF15 Tangerine Metallic	<input type="checkbox"/> TFF30 Acid Green Metallic
Pearlescent Paint	<input type="checkbox"/> TFF17 Midnight Blue Pearl	<input type="checkbox"/> TFF18 White Pearl	
Matte Paint	<input type="checkbox"/> TFF19 Black Lava	<input type="checkbox"/> TFF27 Titanium Grey	
Candy Paint	<input type="checkbox"/> TFF22 Ruby Red	<input type="checkbox"/> TFF23 Azure Blue	<input type="checkbox"/> TFF25 Electric Plum
Matte Textured Paint	<input type="checkbox"/> TFF40 Hyper Blue	<input type="checkbox"/> TFF41 Flame Orange	<input type="checkbox"/> TFF42 OD Green STD

SEAT WIDTH A Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

XSW1	<input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19"	STD	
XSW2	<input type="checkbox"/> 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" (21" and 22" Seat Widths A are only available with Heavy Duty Frame.)	\$495	E2201



Permobil.com
2701 West Court Street | Pasco, Washington 99301
509.586.6117 Tel | 800.545.2266 Toll-Free
509.586.2413 Fax | 866.586.2413 Toll-Free Fax



TiLite 2GX Order Form

Effective November 15, 2019 (New Items in Red)

FOLDING FRAME

Page 2 of 7

FRAME MEASUREMENTS



SEAT DEPTH ③ Measure from 1/2" in front of back post to the front edge of seat sling. (12"-14" Seat Depths ③ not available with 25" (559) and 26" (560) Rear Wheels.)

XSD1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19"..... **STD**

XSD2 ☐ 20"..... **\$525**

XSD3 ☐ Custom Frame Depth: _____" Frame Depth Longer Than Seat Depth (Use 1" increments with a maximum of 2". A

Custom Frame Depth with a maximum of 20" is \$475.)

FRONT SEAT HEIGHT ④ Measure from floor to top of seat tube at front edge of seat sling. Adjustable in 1/2" increments. (13.5" - 14"

Front Seat Height ④ requires Hemi Height Frame Style. 20.5" - 21.5" Front Seat Height ④ not available with Hemi Height Frame Style.)

FRONT SEAT HEIGHT RANGES..... **STD**

Front Wheel Size	Front Seat Height ④	Hemi Front Seat Height ④
3" Front Wheel	14.5" - 18.5"	13.5" - 17.5"
3" Front Wheel w/ Glide™	16.5" - 17.5"	15.0" - 17.0"
4" Front Wheel	15.5" - 19.0"	14.5" - 18.0"
4" Front Wheel w/ Glide™	17.0" - 18.0"	15.5" - 17.0"
5" Front Wheel	16.5" - 19.5"	16.0" - 18.5"
5" Front Wheel w/ Glide™	17.5" - 18.5"	16.5" - 18.0"
6" Front Wheel	17.5" - 20.0"	16.5" - 19.0"
6" Front Wheel w/ Glide™	19.0"	17.5" - 18.0"
8" Front Wheel	19.5" - 21.0"	18.5" - 19.5"

XFSH1 ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17.0" ☐ 16.5" ☐ 16.5" ☐ 17.0"
☐ 17.5" ☐ 18.0" ☐ 18.5" ☐ 19.0" ☐ 19.5" ☐ 20.0" ☐ 20.5" ☐ 21.0"..... **NCO**

REAR SEAT HEIGHT ⑤ Measure from floor to top of seat tube at back post. Adjustable in 3/8" increments.

Rear Wheel Size	Rear Seat Height
22"	13.0" - 18.0"
24"	14.0" - 19.0"
25"	14.0" - 19.0"
26"	15.0" - 20.0"

XRSH1 ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5"
☐ 19" ☐ 19.5" ☐ 20"..... **NCO**

BACK SUPPORTS

SSH3 Back Support Installation ☐ Install Back OR ☐ Do Not Install Back..... **NCO**

ROHO® AGILITY® BACK SUPPORT MIN CONTOUR (Comes standard with quick release hardware.)

Shell Height ☐ 10" ☐ 13" ☐ 16" ☐ 18"

BRS5 ☐ Foam..... **\$510** **E2611**

BRS6 ☐ Foam + Lumbar..... **\$615** **E2613**

BRS7 ☐ Foam w/ Air Insert..... **\$720** **E2620**

BRS8 ☐ Foam w/ Air Insert + Lumbar..... **\$720** **E2615**

ROHO® AGILITY® BACK SUPPORT MID CONTOUR (Comes standard with quick release hardware.)

Shell Height ☐ 10" ☐ 13" ☐ 16" ☐ 18"

BRS9 ☐ Foam w/ Air Insert..... **\$730** **E2620**

BRS10 ☐ Foam w/ Air Insert + Lumbar..... **\$730** **E2615**

TiLite 2GX Order Form

Effective November 15, 2019 (New Items in Red)

FOLDING FRAME

Page 3 of 7

BACK SUPPORTS

COMFORT ACTA-EMBRACE™

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16"

SBR516 ☐ Acta-Embrace Back with 2-PT Quick Release Compass Hardware.....\$560 E2613

CPSB1 ☐ Comfort Privacy Shield.....\$40

COMFORT ACTA-RELIEF

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16"

SBR513 ☐ Acta-Relief Back with 2-PT Quick Release Compass Hardware.....\$730 E2615

CPSB1 ☐ Comfort Privacy Shield.....\$40

SEAT BACK HEIGHT ③/SEAT BACK/PUSH HANDLES Measure from top of back post to top of seat tube at rear of frame.

SEAT BACK TYPE

XSB13 ☐ Fixed Angle, Adjustable-Height Seat Back.....STD

XLB1 ☐ Folding Stabilizer Bar (Required with Transit Tie-Down Option unless Depth Adjustable Backrest is selected. Includes Tension Adjustable by Strap Back Upholstery.).....\$390

XSB14 ☐ Fixed Angle, Adjustable Height Seat Back w/8-Degree Bend. (Not available with Heavy Duty Frame.).....NCO

XSB18 ☐ Folding Adjustable-Angle and Height Seat Back (Not available with Heavy Duty Frame or Ride Designs® Corbac™.).....\$230

XSB16 ☐ Depth Adjustable Back/Frame (Not Available with Hemi Height Frame, Heavy Duty Frame, 19" or 20" Seat Depths ③, Swing Away Tubular or Swing Away Flip-Back Armrest. Seat Depths ③ of 12"-17" allow 3" of rearward adjustment. Seat Depth ③ of 18" allow 2" of rearward adjustment. Requires Tension Adjustable Seat Upholstery.).....NCO

PUSH HANDLES

XSB6 ☐ Integrated Push Handles.....NCO

XSB7 ☐ Integrated Fold-Down Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with Heavy Duty Frame.).....\$115

XSB9 ☐ Bolt-On Height Adjustable Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with Heavy Duty Frame, Folding Adjustable-Angle and Height Seat Back, Swing Away Tubular Armrests or Swing Away - Flip Back Tubular Armrests.).....\$335

XSB20/XSB14 **FIXED ANGLE ADJUSTABLE SEAT BACK HEIGHT ③** (11"-15" Fixed Seat Back Height not available with Heavy Duty Frame, Swing Away Tubular Arms 12"-15", Bolt-On Height Adjustable Push Handles or Folding Stabilizer Bar.)

Push Handles (Selected on Pg. 2)		Short	Tall	Set at: (Use 1" Increments)
None	<input type="checkbox"/>	11.0" - 15.0"	<input type="checkbox"/>	16" - 20"
Integrated	<input type="checkbox"/>	11.0" - 15.0"	<input type="checkbox"/>	16" - 20"
Integrated Fold-Down	<input type="checkbox"/>	11.0" - 15.0"	<input type="checkbox"/>	16" - 20"
Bolt-On Height Adj.	<input type="checkbox"/>	-----	<input type="checkbox"/>	16" - 20"

XSB18 **FOLDING ADJUSTABLE- ANGLE AND HEIGHT SEAT BACK ③** (Not available with a Folding Stabilizer Bar or Heavy Duty Frame)

Push Handles (Selected on Pg. 2)		Short	Medium	Tall	Set at: (Use 1" Increments)
None	<input type="checkbox"/>	11.0" - 13.0"	<input type="checkbox"/>	13.5" - 16.5"	<input type="checkbox"/>
Integrated	<input type="checkbox"/>	12.5" - 14.5"	<input type="checkbox"/>	13.5" - 16.5"	<input type="checkbox"/>
Integrated Fold-Down	<input type="checkbox"/>	15.0" - 17.0"	<input type="checkbox"/>	16.0" - 19.0"	<input type="checkbox"/>
				18.5" - 20.5"	<input type="checkbox"/>

FRAME MEASUREMENTS

CAMBER ⑥

XCBR8 ☐ 0° ☐ 2° ☐ 4° ☐ 6°.....STD

FRONT ANGLE ① Front Angle ① is measured from inside of front frame tube to top of bottom frame tube.

Swing Away Front End

XFA3 ☐ 70° ☐ 80° ☐ 85° (Not available with 6" Front Wheels.) ☐ 90° (Requires aluminum footplate. Not available with Heavy Duty Frame, 5" or 6" Front Wheels.).....STD

XFA0 ☐ Omit Hangers.....NCO

Fixed Front End

XFA6 ☐ 70° ☐ 80° ☐ 85° (85° Not available with 6" Front Wheels.).....STD

SEAT TO FOOTREST ①/FOOTRESTS Measure from front edge of seat sling to top rear of footrest.

If your Front Angle ① is 70°, 80°, OR 85° select one Seat to Footrest ① length and one Footrest type:
Standard Footrests Measurement must be at least 2" shorter than Front Seat Height ③. 70° Fixed Front has a minimum seat to footrest of 14.5"

☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"

XFTR5 ☐ Flip-up Footrest.....STD

XFTR9 ☐ Depth Adjustable Flip-up Footrest.....\$110 K0045

XFTR10 ☐ Depth & Angle Adjustable Flip-up Footrest (Available for 15"-20" Seat to Footrest ① only.).....\$270 K0040, K0045

XFTR11 ☐ Angle Adjustable Flip-up Footrest (Available for 15"-20" Seat to Footrest ① only.).....\$170 K0040, K0045

TiLite 2GX Order Form

Effective November 15, 2019 (New Items in Red)

FOLDING FRAME

Page 4 of 7

FRAME MEASUREMENTS

High Mount Footrests (Not Available with Luggage Carriers, or Impact Guards.)

	<input type="checkbox"/> 7" <input type="checkbox"/> 7.5" <input type="checkbox"/> 8" <input type="checkbox"/> 8.5" <input type="checkbox"/> 9" <input type="checkbox"/> 9.5" <input type="checkbox"/> 10" <input type="checkbox"/> 10.5" <input type="checkbox"/> 11" <input type="checkbox"/> 11.5" <input type="checkbox"/> 12" <input type="checkbox"/> 12.5" <input type="checkbox"/> 13" <input type="checkbox"/> 13.5" <input type="checkbox"/> 14" <input type="checkbox"/> 14.5"		
XFTR20	<input type="checkbox"/> High Mount Flip-up.....	\$190	K0037
XFTR15	<input type="checkbox"/> Angle Adjustable High Mount Flip-up Footrest.....	\$325	K0040, K0037
XFTR16	<input type="checkbox"/> Depth & Angle Adjustable High Mount Flip-up Footrest.....	\$375	K0040, K0037

If your Front Angle ① is 90°, select one Seat to Footrest ① length and one Footrest type:

	<input type="checkbox"/> 7" <input type="checkbox"/> 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13"		
XFTR5	<input type="checkbox"/> Flip-up Footrest.....	STD	
XFTR9	<input type="checkbox"/> Depth Adjustable Flip-up Footrest.....	\$110	K0037
XFTR10	<input type="checkbox"/> Depth & Angle Adjustable Flip-up Footrest.....	\$270	K0040, K0037
XHLP1	<input type="checkbox"/> Heel Loops (Not available with 13" or narrower Seat Width ① or 1-Piece Footrests.).....	\$50	E0951

CUSTOM CONFIGURATION OPTIONS

FRONT WHEELS (See Front Seat Height Ranges Chart on Page 2 for Front Seat Height ② limitations.)

XFW24	<input type="checkbox"/> 3"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub.....	\$190	
XFW14	<input type="checkbox"/> 4"x.75" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame or Heavy Duty Frame).....	NCO	
XFW3	<input type="checkbox"/> 4"x.75" Light-Up Micro-Caster (Not available with Power Adaptable Frame or Heavy Duty Frame.).....	\$80	
XFW26	<input type="checkbox"/> 4"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Assist, Power Add-on System or Heavy Duty Frame.).....	\$95	
XFW27	<input type="checkbox"/> 4"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub.....	\$200	
XFW6	<input type="checkbox"/> 5"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or with 90° Front Angle ①.).....	NCO	
XFW28	<input type="checkbox"/> 5"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or 90° Front Angle ①.).....	\$95	
XFW30	<input type="checkbox"/> 5"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub (Not available with 90° Front Angle ①.).....	\$215	
XFW9	<input type="checkbox"/> 6"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or 85° or 90° Front Angle ①.).....	NCO	
XFW32	<input type="checkbox"/> 6"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame, 85° or 90° Front Angle ①.).....	\$100	
XFW33	<input type="checkbox"/> 6"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub (Not available with 85° or 90° Front Angle ①.).....	\$225	
XFW16	<input type="checkbox"/> 8"x1" Plastic Wheel w/Poly Tire (Not available with Power Assist, Power Add-on, Heavy Duty Frame or 12"-14" Seat Depths. Available with 70° Front Angle ① only.).....	NCO	

FRONT FORKS

XFK1	<input type="checkbox"/> TiLite Standard Fork.....	STD	
XGL1	<input type="checkbox"/> Out-Front Glide™ Suspension Fork (Available in Black anodize only. Not available with 8" Front Wheels or Heavy Duty Frame). Rider's Weight: <input type="checkbox"/> 125 lbs. <input type="checkbox"/> 125-180 lbs. <input type="checkbox"/> 180-265 lbs.	\$485	E1015

ONE ARM DRIVE

XOAD1	<input type="checkbox"/> One Arm Drive (Available with 24" Mag Wheels only. Rear Seat Height, Camber or axle restrictions may apply.)..... Drive Handrim Side (Available with standard aluminum handrim only.) <input type="checkbox"/> Install on Left <input type="checkbox"/> Install on Right	\$1,200	E0958
-------	--	---------	-------

STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height ③ limitations.)

XRW14	TiLite SHADOW™.....	STD	
	<input type="checkbox"/> 22"(501)		
	<input type="checkbox"/> 24"(540)		
	<input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth ③.)		
	<input type="checkbox"/> 26"(590) (Requires 15" minimum Seat Depth ③.)		

PERFORMANCE REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height ③ limitations.)

XRW7	Spinergy® SPOX <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White.....	\$845	
	<input type="checkbox"/> 22"(501) - Std. Hub		
	<input type="checkbox"/> 24"(540) - Std. Hub		
	<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth ③.)		
	<input type="checkbox"/> 26"(590) - Sport Hub (Requires 15" minimum Seat Depth ③. Not available with Quad Release Axles.)		
XRW8	Spinergy® LX (Available with Heavy Duty Frame only when 0 or 2 degree camber is selected.) <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Bright Orange <input type="checkbox"/> Flamingo Pink <input type="checkbox"/> Purple.....	\$1,060	
	<input type="checkbox"/> 22"(501) - Std. Hub		
	<input type="checkbox"/> 24"(540) - Std. Hub		
	<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth ③.)		
	<input type="checkbox"/> 26"(590) - Sport Hub (Requires 15" minimum Seat Depth ③. Not available with Quad Release Axles.)		
XRW16	Spinergy® Carbon Blade (Not available with Heavy Duty Frame.).....	\$1,480	
	<input type="checkbox"/> 24"(540)		
	<input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth ③.)		
XRW17	X-Core Mag (Not available with Heavy Duty Frame, Short Tabs, Quad Release Axles, Spoke Guards or Omit Handrims Option. Requires Aluminum Silver Anodized Handrims.).....	NCO	
	<input type="checkbox"/> 22"(501)		
	<input type="checkbox"/> 24"(540)		

TiLite 2GX Order Form

Effective November 15, 2019 (New Items in Red)

FOLDING FRAME

Page 5 of 7

CUSTOM CONFIGURATION OPTIONS

AXLES

XAXL1	<input type="checkbox"/> Stainless Quick Release	STD
XAXL2	<input type="checkbox"/> Titanium Quick Release	\$255
XAXL3	<input type="checkbox"/> Quad Release (Not available with Mag or 26"(590) Spinergy Spox or Spinergy LX Rear Wheels.)	\$155

PNEUMATIC REAR TIRES

XRTR1	<input type="checkbox"/> Treaded (Not available with 25"(559) Rear Wheels.)	STD
XRTR2	<input type="checkbox"/> Airless Inserts	\$90
XRTR3	<input type="checkbox"/> 25"(559) Primo Treaded (Standard with 25"(559) Rear Wheels.)	STD
XRTR4	<input type="checkbox"/> Primo Slick	NCO
XRTR30	<input type="checkbox"/> 24"(540), 25"(559) and 26"(590) Primo Sentinel	\$125
XRTR8	<input type="checkbox"/> 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution	\$140
XRTR9	<input type="checkbox"/> 24"(540) and 25"(559) Schwalbe® Slick	\$145
XRTR24	<input type="checkbox"/> 24"(540) Kenda Kobra Knobby	\$95

E2213

SOLID REAR TIRES

XRTR22	<input type="checkbox"/> 22"(501), 24"(540) and 25"(559) Primo Express (Not available with Spinergy Carbon Blade Rear Wheels.)	\$110
XRTR23	<input type="checkbox"/> 24"(540) and 25"(559) Primo Active	\$135

HANDRIMS

XRIM1	<input type="checkbox"/> Aluminum – Silver Anodized	STD
XRIM21	<input type="checkbox"/> Aluminum – Black Anodized	NCO
XRIM12	<input type="checkbox"/> Titanium – Satin Finish (Not available with 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.)	\$470

Handrims by OUT-FRONT®

XRIM10	<input type="checkbox"/> The Surge (All-In-One Oval with Gription™ Strip, Larger Profile)	\$400
XRIM11	<input type="checkbox"/> The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile)	\$430
XRIM20	<input type="checkbox"/> Q-Grip (Exclusive High-Friction, Anti-Slip Coating)	\$330

The Natural-Fit® (Original, Larger Profile) (Not available with 20"(451) Rear Wheels.)

TRIM14	<input type="checkbox"/> Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)	\$395
TRIM15	<input type="checkbox"/> Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)	\$430
TRIM16	<input type="checkbox"/> Without Thumb Grip	\$360

The Natural-Fit® LT (Reduced Weight, Smaller Profile) (Not available with 20"(451) Rear Wheels.)

TRIM17	<input type="checkbox"/> Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)	\$415
TRIM18	<input type="checkbox"/> Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)	\$475
TRIM19	<input type="checkbox"/> Without Thumb Grip	\$390

Additional Handrim Options

XRIM5	<input type="checkbox"/> Plastic Coated – Black	\$160
XRIM22	<input type="checkbox"/> Omit Handrims	NCO

Tab Length

XRIM7	<input type="checkbox"/> Long Tabs	STD
XRIM8	<input type="checkbox"/> Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.)	\$25

WHEEL LOCKS

Wheel Locks by OUT-FRONT® (Composite wheel locks not always compatible with 12"-13" Seat Depth. Uni-locks will be substituted at no charge.)

XWLK7	<input type="checkbox"/> Compact Composite Push to Lock	NCO
XWLK8	<input type="checkbox"/> Compact Composite Pull to Lock	NCO
XWLK10	<input type="checkbox"/> Extension Handles for Compact Composite Push/Pull to Lock	\$85
XWLK9	<input type="checkbox"/> Composite Scissor Lock (Not available on Aero X2 Swing Away.)	\$145

E0961

Additional Wheel Locks

XWLK3	Uni-Lock Set At: <input type="checkbox"/> Push OR <input type="checkbox"/> Pull	\$165
XWLK4	<input type="checkbox"/> Extension Handles for Uni-Lock	\$85
XWLK15	<input type="checkbox"/> Grade Aids Push to Lock (Not available with solid rear tires)	\$175

E0961

E0974

BACK UPHOLSTERY (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

XBUP1	<input type="checkbox"/> Padded Nylon (Black Only. Not available with Heavy Duty Frame, Folding Back, Folding Stabilizer Bar, Depth Adjustable Frame, Fold-Down Push Handles or Ultrasuede Backrest Stripe.)	STD
Tension Adjustable by Straps (Included with Heavy Duty Frame. Required with Aluminum Fendered Side Guards and Fold-Down Push Handles.)		
XBUP3	<input type="checkbox"/> Padded Nylon	\$340
XBUP5	<input type="checkbox"/> AIR Breathable Upholstery (Not available with Ultrasuede Backrest Stripe.)	\$340
XUPC3	Ultrasuede Backrest Stripe <input type="checkbox"/> Black (Not available with Padded Nylon or AIR Breathable Back Upholstery.)	\$105
XBUP0	<input type="checkbox"/> Omit Back Upholstery	NCO

E2611

E2611

SEAT UPHOLSTERY

XSUP4	<input type="checkbox"/> Tension Adjustable (Black only)	STD
XSUP0	<input type="checkbox"/> Omit Seat Upholstery	NCO

TiLite 2GX Order Form

Effective November 15, 2019 (New Items in Red)

FOLDING FRAME

Page 6 of 7

CUSTOM CONFIGURATION OPTIONS

SEAT CUSHION (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

XCSH2	<input type="checkbox"/> Comfort Curve Cushion (Stretch-Air cover)	\$120	E2601
XCSH23	<input type="checkbox"/> Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion	\$360	E2607
XCSH14	<input type="checkbox"/> ROHO Hybrid Elite SR Single Compartment (Fits seat widths 14" to 19")	\$560	E2622
XCSH14	<input type="checkbox"/> ROHO Quadtro Select (Fits seat widths 14" to 19") Profile: <input type="checkbox"/> Low <input type="checkbox"/> Mid <input type="checkbox"/> High	\$525	E2624

SIDE GUARDS

Aluminum

XSDG2	<input type="checkbox"/> Rigid Removable <input type="checkbox"/> Small OR <input type="checkbox"/> Std	\$200	
XSDG3	<input type="checkbox"/> Quad Removable w/Hole <input type="checkbox"/> Small OR <input type="checkbox"/> Std	\$215	
XSDG6	<input type="checkbox"/> Fendered Adjustable (Requires Tension Adjustable by Straps Back Upholstery. Not available with Angle-Adjustable Fold-Down Seat Backs.)	\$335	

Carbon Fiber

XSDG7	<input type="checkbox"/> Rigid Removable <input type="checkbox"/> Small OR <input type="checkbox"/> Std	\$310	
XSDG8	<input type="checkbox"/> Quad Removable w/Hole <input type="checkbox"/> Small OR <input type="checkbox"/> Std	\$320	
XSDG11	<input type="checkbox"/> Fendered Adjustable (Requires Tension Adjustable by Straps Back Upholstery. Not available with Angle-Adjustable Fold-Down Seat Backs.)	\$565	

ARMRESTS (Desk Arms come with built-in Side Guards.)

XARM1	Swing Away Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar with Short or Medium Seat Back Height ranges or Depth Adjustable Back.)	\$225	E0973
	<input type="checkbox"/> 9" - 12" OR <input type="checkbox"/> 12" - 15"		
XARM4	<input type="checkbox"/> Swing Away - Flip Back Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar or Depth Adjustable Back.)	\$315	E0973
XARM2	Removable Height Adjustable DeskArm with Rigid Side Guard (Requires 14" minimum Seat Depth Ⓐ .)	\$335	E0973
	<input type="checkbox"/> 9" to 12" OR <input type="checkbox"/> 12" to 15"		
	<input type="checkbox"/> Standard Length Pad (10") OR <input type="checkbox"/> Full Length Pad (14")		
XARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 16" minimum Seat Depth Ⓐ .)	\$375	E0973
	<input type="checkbox"/> 9" to 12" OR <input type="checkbox"/> 12" to 15"		
	<input type="checkbox"/> Standard Length Pad (10") OR <input type="checkbox"/> Full Length Pad (14")		
XARM8	Standard 2-Point Flip-Back Height Adjustable Desk Arm With Rigid Side Guard	\$295	E0973
	<input type="checkbox"/> 9" to 12" OR <input type="checkbox"/> 12" to 15"		
	<input type="checkbox"/> Standard Length Pad (10") OR <input type="checkbox"/> Full Length Pad (14")		
XARM9	Standard 2-Point Flip-Back Height Adjustable Desk Arm With Rigid Side Guard and Transfer Loop	\$335	E0973
	<input type="checkbox"/> 9" to 12" OR <input type="checkbox"/> 12" to 15"		
	<input type="checkbox"/> Standard Length Pad (10") OR <input type="checkbox"/> Full Length Pad (14")		

Available Armrest Option

XARM5	<input type="checkbox"/> Covered in Black Ultrasuede (Not available with Desk Arms)	\$110	
-------	---	-------	--

SEAT BELTS

XBLT1	<input type="checkbox"/> Velcro Adjustable, Black Nylon Webbing	\$55	E0978
XBLT2	<input type="checkbox"/> Bodypoint® 2" Push-Button Buckle (Auto-Style.)	\$95	E0978
XBLT4	<input type="checkbox"/> Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths Ⓐ .)	\$180	E0978
XBLT5	<input type="checkbox"/> Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths Ⓐ .)	\$170	E0978
XBL10	<input type="checkbox"/> Bodypoint® Optional Band Clamp Mounting Hardware	\$20	K0108
XBLT16	<input type="checkbox"/> Bodypoint® EvoFlex® Pelvic Stabilizer Push Button	\$185	E0978

CALF STRAP

XCLF2	Velcro Adjustable	\$40	K0038
	<input type="checkbox"/> Small (Recommended for 9"- 12" Footrest Widths)		
	<input type="checkbox"/> Medium (Recommended for 11"- 15" Footrest Widths)		
	<input type="checkbox"/> Large (Recommended for 14"- 17" Footrest Widths)		
	<input type="checkbox"/> Xtra Large (Recommended for 16"- 19" Footrest Widths)		
XCLF3	Bodypoint® Padded Velcro Adjustable	\$100	K0038
	<input type="checkbox"/> Medium (Recommended for 12"- 15" Footrest Widths)		
	<input type="checkbox"/> Large (Recommended for 16"- 20" Footrest Widths)		
XCLF1	<input type="checkbox"/> Buckle Adjustable, Black Nylon Webbing	\$30	K0038

RESIDUAL LIMB SUPPORT

XRSL1	<input type="checkbox"/> Residual Limb Support (8" width x 10" length)	\$395	E1020
XRSL2	<input type="checkbox"/> Residual Limb Support Pair (8" width x 10" length)	\$790	E1020

BACK PACK

XBCK3	<input type="checkbox"/> Black Permobil Back Pack	\$95	
-------	---	------	--

SEAT POUCH

XPCH3	<input type="checkbox"/> Black Permobil Seat Pouch	\$60	
XPCH4	<input type="checkbox"/> Bodypoint® Mobility Bag	\$100	

TiLite 2GX Order Form

Effective November 15, 2019 (New Items in Red)

FOLDING FRAME

Page 7 of 7

CUSTOM CONFIGURATION OPTIONS

SPOKE GUARDS

XGRD1 ☐ Clear with TiMotion Logo (Not available with MagRear Wheels.) \$135 K0065

ANTI-TIPS

XTIP1 ☐ Rear Aluminum - Black \$160 E0971

XTIP2 ☐ User Friendly Flip-up \$335 E0971

XTIP9 ☐ Tip Assist \$105 E0971

ACCESSORIES

XUEST1 Upper Extremity Support Tray ☐ 12" to 14" OR ☐ 13" to 15" Seat Widths **A** (Requires Armrests.) \$290 E0950

XLUG1 ☐ Luggage Carriers (Not available with High Mount Swing Away Footrests or 1-Piece Flip-up Footrest.) \$145

XMPCT2 Neoprene Impact Guards with Ultrasuede Stripe (Not available with High Mount Footrests, 90° Hangers or Elevating Leg Rest.)
☐ Black \$105

XCH1 ☐ Crutch Holder \$225

XAFP1 4" Logo Patch ☐ Army ☐ Navy ☐ Air Force ☐ Marines \$105

GENERAL

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobilus.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Straint® is a registered trademark of Q'Straint Inc.

QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR